UC15Implant System

urisimplants.com





OMNI Fixture

OMNI Straight Fixture

O4 OMNI Tapered Fixture

05 Cover Screw

06 Healing Abutment



Surgical Instruments

10 Surgical KIT

17 Surgical Drilling Protocol

Surgical Guide Instruments

20 Pylon Kit

22 Pylon Plus Kit

29 Pylon Kit Guide Drilling Protocol

31 Pylon Crestal Sinus Kit

37 Sinus Guide Drilling Protocol

38 Pylon Anchor Kit

Surgery Instruments

41 Instrument

42 Prosthetic Instrument

PROSTHETICS

Fixture Level Impression

46 Basic Abutment

52 Impression Post [Pick Up/Transfer Type]

Fixture Analog

55 Angled Abutment

58 Milling Abutment

61 Temporary Abutment

62 Abutment Screw

Abutment Level Impression

64 Direct Abutment

67 Direct Abutment Components

68 Protective Cap

69 Impression Cap

Burn Out Coping

70 Direct Abutment Analog

73 Multi-Unit Abutment

74 Straight Abutment Angled Abutment

75 Multi-Unit Abutment Screw

77 Multi-Unit Abutment Components

78 Multi-Unit Healing Cap

Multi-Unit Ti Cylinder

Multi-Unit Burn Out Cylinder

Multi-Unit Temporary Cylinder

Multi-Unit Cylinder Screw

Multi-Unit Impression Post [Pick Up/Transfer Type]

Multi-Unit Abutment Analog

Overdenture

80 Ball Abutment

Ball Abutment Components

81 Retainer Cap

Retainer

O-ring

Ball Abutment Analog

PROSTHETICS

83 T:Loc Abutment

T:Loc Components

84 T:LOC Component Set

85 T:LOC Titanium Cap

Processing Replacement Male

Retention Replacement male

Block out spacer

T:LOC Impression Coping

T:LOC Analog

Digital Components

87 Fixture Level Digital Components

88 Uris Base

89 Scan Body

Fixture Digital Analog

90 Multi-Unit Abutment Digital Components

91 Multi-Unit Base

Multi-Unit Scan Body

Multi-Unit Digital Analog

Multi-Unit Conversion Scan Body

92 Prosthetic Driver Kit

Packages & Instructions

94 Fixture Packaging

95 Instruction for Use

97 Abutment Packaging

98 Prosthetic Manual

100 General Principles of Surgical Tool Management

OMNI Implant System

Versatile Prosthetic Options

powered by TruAbutment







Custom Healing



Custom Abutment



T:LO Custo Remova



ASC Angulated Screw Channe

11° Morse Taper

- Creates hermetic seal
- Elimination of microgar
- Reduction of micro-movemen and screw loosening.

Platform Switching

- Reduction of peri-implant bone loss
- Greater magnitude with increased diameter
- Aids in maintaining biological width without
disruption during clinical functions
- Beveled feature facilitates bone growth
above the shoulder

Cutting Edge

- Triple cutting edge
- sen-tapping.

SEM x500

Preserves more of the peri-implant bone, stabilizes more of the soft tissues reduces the microgap size found in the abutment-implant connection and proper geometry for narrower mesio-distal edentulous spaces.





Blood pockets

- 6.5 ~12 % Reduction: Reduction of stress / pressure points
- Creation of "blood pockets" promotes bone growth during osseointegration
- Gradual distribution of force along the cortical plate




Thread Design

- The bite of the thread features two-stage ram
- Allows for gradual bone condensing for wide range of bone conditions
- Even insertion while protecting bone structure.



Apex design

- Aggressive triple cutting edge allows for adjustable implant orientation during manua insertion for optimal final placement.
- Round apex to protect sinus floor, nerve canal or other anatomical structures



- Hybrid design provides best features of both straight and tapered implants.
- Coronal taper with unique design aids in initial stability, load dissipation and reduces stress/strain on the crestal bone.
- The body and apical design assist in guiding the implant to the desired position while providing bone density distribution.

OMNI Fixture

Length	8.5	10.0	11.5	13.0	14.5
D Ø3.1					
	OMF 300 85	OMF 30100	OMF 30115	OMF 30130	OMF 30145
Length	8.5	10.0	11.5	13.0	14.5
D Ø3.5					
	OMF 35085	OMF 35100	OMF 35115	OMF 35130	OMF 35145



OMF 45070

OMF 45085

Length	7.0	8.5	10.0	11.5	13.0	14.5
	OMF 400 70	OMF 40085	OMF 40 100	OMF 40115	OMF 40130	OMF 40145
		3	***************************************			
D Ø4.0	<u></u>					
Length	7.0	8.5	10.0	11.5	13.0	14.5

OMF 45100

OMF 45115

OMF 45130

OMF 45145

Length	7.0	8.5	10.0	11.5	13.0	14.5
D Ø5.0						
	OMF 500 70	OMF 500 85	OMF 50100	OMF 50115	OMF 50130	OMF 50145
Length	7.0	8.5	10.0	11.5	13.0	14.5
D Ø5.5						
	OMF 550 70	OMF 550 85	OMF 55100	OMF 55115	OMF 55130	OMF 55145



OMF 65100

OMF 65070

OMF 65085





8.5

Hex 2.5



Cover

Hex driver: 1.27 | Torque: 5~10 Ncm






UCS 16



UCS 20

URIS Implant System

14.5







Cuff/Height	1.0 / 2.0	2.0 / 3.0	3.0 / 4.0	4.0 / 5.0	5.0 / 7.0	2.0 / 2.0
D Ø4.0						
	UHAN 40 12	UHAN 40 23	UHAN 40 34	UHAN 40 45	UHAN 40 57	UHAN 4022
Cuff/Height	1.0 / 2.0	2.0 / 3.0	3.0 / 4.0	4.0 / 5.0	5.0 / 7.0	2.0 / 2.0
D Ø4.5						T
	UHAN 4512	UHAN 4523	UHAN 45 34	UHAN 45 45	UHAN 45 57	UHAN 4522



1.0 / 2.0	2.0 / 3.0	3.0 / 4.0	4.0 / 5.0	5.0 / 7.0	2.0 / 2.0
					V
UHA 4012	UHA 4023	UHA 40 34	UHA 40 45	UHA 40 57	UHA 40 22
1.0 / 2.0	2.0 / 3.0	3.0 / 4.0	4.0 / 5.0	5.0 / 7.0	2.0 / 2.0
	UHA 4012	UHA 4012 UHA 4023	UHA 4012 UHA 4023 UHA 4034	UHA 4012 UHA 4023 UHA 4034 UHA 4045	UHA 4012 UHA 4023 UHA 4034 UHA 4045 UHA 4057





UHA 4523










Cuff/Height	1.0 / 2.0	2.0 / 3.0	3.0 / 4.0	4.0 / 5.0	5.0 / 7.0	2.0 / 2.0
D Ø5.5						Y
Cuff/Height	UHA 5512	UHA 5523	UHA 5534	UHA 5545	UHA 55 57	UHAN 5522 2.0 / 2.0
D Ø6.5	1.0 / 2.0	2.0 / 5.0	3.0 / 4.0	4.0 / 3.0	3.0/1.0	2.0 / 2.0
	T					
	UHA 65 12	UHA 65 23	UHA 65 34	UHA 65 45	UHA 65 57	UHA 6522
Cuff/Height	1.0 / 2.0	2.0 / 3.0	3.0 / 4.0			2.0 / 2.0
D Ø7.5	T					T
	UHA 7512	UHA 7523	UHA 75 34			UHA 75 22

SurgeryInstruments

Scale 1:1.5







Point Drills through the cortical bone to create an ideal path for the next drills at the selected site. 800~1200 rpm (must use with saline water)

Round Levels out the uneven bone and to remove the remaining gingiva residue. 800~1200 rpm






RD 1726

Drill

Side Cut Drills with the side blade to modify the drilling path.

800~1200 rpm

(Unit: mm)





Pilot Expands cortical bone to let the final drill enter the path easily. 800~1200 rpm




PD 26070


Initial Expands the prepared site and used according to the chosen fixture's length.

Drill	800~1200 rpm		Ü	Ü		
(Unit: mr	n)					
Length	7.0	8.5	10.0	11.5	13.0	14.5
D Ø2.0	0.5		T ozce i		14.5 13.5	↑ 115 ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑ ↑
	ISD 200 70	ISD 20085	ISD 20100	ISD 20115	ISD 20130	ISD 20145

Final The last drill used to place the Uris fixture and is available for each fixture diameter/length and bone density. 800~1200 rpm









Fixture

Driver

*Caution: Must torque after fully engaging into the fixture securely to avoid any breakages.





Parallel

Confirms and indicates the bone preparation location.



Confirms the path, the gingiva depth, the hex direction of the prosthetic connection after the placement.







Surgical Drilling Protocol

Required O Use at your discretion Requires a drill that is 1.5 mm shorter than the fixture

Implant Ø3.0 x 11.5mm




Implant Ø3.5 x 11.5mm




Implant Ø4.0 x 11.5mm

O Use at your discretion

Requires a drill that is 1.5 mm shorter than the fixture



Implant Ø4.5 x 11.5mm




Implant Ø5.0 x 11.5mm

O Use at your discretion


Required

Requires a drill that is 1.5 mm shorter than the fixture




Required












PBP 5027

PBP 3526



(Unit: mm)

Removes soft tissue during a flap-less surgery. Path Increases the accuracy of the drilling path.




Tissue Punch / Path Drill / Initial Drill

PTP 3227

Low speed: 50~100 rpm; within 5 seconds with high torque | High speed: 800~1200 rpm (use irrigation)

Initial



Final Used as the last drill to place an Uris fixture and is selected according to the fixture diameter and bone density. Low speed: 50~100 rpm; within 5 seconds with high torque


High speed: 800~1200 rpm (use irrigation)

	igii speed. 600° 12	oo ipiii (use iiii)	gation					
(Unit: mm)								
Length 1	7.0	8.5	10.0	11.5	13.0	14.5	16.0	17.5
Length 2	17.1	18.6	20.1	21.6	23.1	24.6	26.1	27.6
D Ø3.0	6 L2	(D)	0)	\$11×0+0	114000	0	0	0 -
	NSD 300 70	NSD 300 85	NSD 30100	NSD 30115	NSD 30130	NSD 30145	NSD 30160	NSD 30175

Length 1	7.0	8.5	10.0	11.5	13.0	14.5	16.0	17.5
Length 2	17.1	18.6	20.1	21.6	23.1	24.6	26.1	27.6
D Ø3.4		- 10 × × × × × × × × × × × × × × × × × ×	(D) a story (D)	00 - 00	0)	90 × × H+	-10 + 10 O	0
	NSD 340 70	NSD 340 85	NSD 34100	NSD 34115	NSD 34130	NSD 34145	NSD 34160	NSD 34175

Length 1	7.0	8.5	10.0	11.5	13.0	14.5	16.0	17.5
Length 2	17.1	18.6	20.1	21.6	23.1	24.6	26.1	27.6
D Ø3.9		0)	(D) (C) (C)	0 (0.08 x 115	(D) (D) (D)	- 10 x 5 10 10 10 10 10 10 10 10 10 10 10 10 10	-014 610) O	O CONTROL
	NSD 390 70	NSD 390 85	NSD 39100	NSD 39115	NSD 39130	NSD 39145	NSD 39160	NSD 39175







SocketCuts the ridge of the extracted site by the fixture's diameter size, preventing drill slippage caused by remaining bone residue. 800~1200 rpm


(Unit: mm)

Diameter Ø3.5 Ø4.5 Ø5.0 Ø4.0 ESD 35 ESD 40 ESD 45 ESD 50


Contours the bone around the coronal aspect to facilitate full seating of the abutment after removing the surgical guide. 800~1200 rpm



Guide Drilling Protocol



Guide Sleeve Offset 11.0



*Caution: Not recommended on hard bone. Do NOT torque over 25 Ncm and should not drill deeper than 6.5 mm. Needs to be torqued with a ratchet driver as the final step.





Ratchet

Recommended Torque Value: 30 Ncm ~ 45 Ncm









Recommended Torque Value: 30 Ncm ~ 45 Ncm

Stopper



GFDR 20095












Sleeve




NS 53GR

Guide Sleeve Offset 12.5







Drill

Initial Enlarges the osteotomy site. Do NOT drill deeper than 0.5~1.0 mm directly below the maxillary sinus floor. Low speed: 50~100 rpm; within 5 seconds with high torque High speed: 800~1200 rpm (use irrigation)

(Unit: mm) Length 12.0 13.0 14.0 15.0 16.0 D Ø2.0 PSID 20120 PSID 20130 PSID 20140 PSID 20150 PSID 20160

Drill

Step Used after the D2.0 initial drill to enlarge the osteotomy site. Do NOT drill deeper than 0.5~1.0 mm directly below the maxillary sinus floor. Low speed: 50~100 rpm; within 5 seconds with high torque High speed: 800~1200 rpm (use irrigation)

(Unit: mm)					
Length	12.0	13.0	14.0	15.0	16.0
D Ø2.7	9)	-27 X 13	Q)	Q)	Q) 27846
	PSMD 27120	PSMD 27130	PSMD 27140	PSMD 27150	PSMD 27 16 0





Drill

Sinus Used to approach the sinus membrane and should drill 1~2 mm deeper than the step drill. Low speed: 50~100 rpm; within 5 seconds with high torque (*Drill while applying pressure) High speed: PROHIBITED

(Unit: mm)



Membrane Elevator

Carrier

Delivers the membrane elevator. Needs to be engaged onto a handle.

Membrane **Elevator**

Lifts sinus membrane with hydraulic pressure.

Needs to be engaged onto a carrier.



PHME 6537

Final

Used with a surgical guide after bone grafting. Low speed: 50~100 rpm; within 5 seconds with high torque



Membrane Elevator

Tube

18.0

PSFD 43180

A transparent silicon tube. Outer D4.0/Inner D2.0/Length 300 mm | Autoclave before use | Single use only





Bone

Condenser

Inserts bone graft materials into inferior wall once the wall is completely opened.

Used after removing the surgical guide.



Measures depth with measurement markings at the tip.





PSDG 28195

Bone Condenser


Handle




PSBC 9095

Pylon Crestal Sinus Kit Guide Drilling Protocol


Implant Ø5.0




Guide Sleeve Offset 9.5



Guide Sleeve Offset 11.0



Guide Sleeve Offset 12.5













AD 1413



Anchor Holds the surgical guide onto the bone. Mostly used in an edentulous case.

AP 15

Fixation

Holds the surgical guide after the 2.0 drill before the fixture delivery.





PAI 2007



Drill **Extension**



DEH 4314

Torque

Wrench



TW 40

Fixation



Screw Holds the surgical guide after the fixture delivery.

Offset 9.5 11.0 12.5 PAF 095N PAF 110N PAF 125N









Anchor

Driver Tightens down the anchor pin.





PADH 23

Depth

Gauge Measures the depth of a drilling hole.

Titanium Bowl





TB 4762

URIS Implant System



1.27 Hex

Driver

*Caution: Must be properly engaged and used at the recommended torque value.





Driver the ball abutment.

Multi-Unit Compatible with Straight

Multi-Unit abutment.

Straight Driver T:LOC abutments

Compatible



BADR 24





MUDR 33



TLC-TLSD13

Angulated Screw Channel

Must be properly engaged and torqued at the recommended torque value of 20~25 Ncm within 25°.



Removal

Driver

Used to remove the abutment off the fixture. Needs to be inserted upright and tightened clockwise.



Extension

Adapter

Extends the instrument's length.








Ratchet

Fixture Level Impression Flowchart




Prosthetics

Basic Abutment

Нех Туре

Hex driver : 1.27 Torque: Narrow[20 Ncm]/Regular[30 Ncm] Packing unit : Basic Abutment_Hex Type + Abutment Screw Order Code : ex) UDA 5534H











Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø4.0		1074	4024	4034	404	4054	4064
Height	4.0 5.5 7.0	UDAN 4014H UDAN 4015H UDAN 4017H	UDAN 4024H UDAN 4025H UDAN 4027H	UDAN 4034H UDAN 4035H UDAN 4037H	UDAN 4044H UDAN 4045H UDAN 4047H	UDAN 4054H UDAN 4055H UDAN 4057H	UDAN 4064H UDAN 4065H UDAN 4067H
Cuff D Ø4.5		1.0	2.0	3.0	4.0	5.0	6.0
		4514	4524	4534	4544	1,455	4564
Height	4.0 5.5 7.0	UDAN 4514H UDAN 4515H UDAN 4517H	UDAN 4524H UDAN 4525H UDAN 4527H	UDAN 4534H UDAN 4535H UDAN 4537H	UDAN 4544H UDAN 4545H UDAN 4547H	UDAN 4554H UDAN 4555H UDAN 4557H	UDAN 4564H UDAN 4565H UDAN 4567H



Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø4.5						(A)	4564
		1524	4524	1534	154	PSSP	
Height	4.0 5.5	UDA 4514H UDA 4515H	UDA 4524H UDA 4525H	UDA 4534H UDA 4535H	UDA 4544H UDA 4545H	UDA 4554H UDA 4555H	UDA 4564H UDA 4565H
incigiit	7.0	UDA 4517H	UDA 4527H	UDA 4537H	UDA 4547H	UDA 4557H	UDA 4567H

Hex driver : 1.27 Torque: Narrow[20 Ncm]/Regular[30 Ncm] Packing unit : Basic Abutment_Non Hex Type + Abutment Screw Order Code : ex) UDA 5534N	Cuff Fixture Level
6564 (6	4014 ()
2534	

Basic Abutment

Non Hex Type

Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø5.5		7755	18254	5534	8544	P588	18264
Height	4.0 5.5 7.0	UDA 5514H UDA 5515H UDA 5517H	UDA 5524H UDA 5525H UDA 5527H	UDA 5534H UDA 5535H UDA 5537H	UDA 5544H UDA 5545H UDA 5547H	UDA 5554H UDA 5555H UDA 5557H	UDA 5564H UDA 5565H UDA 5567H
Cuff D Ø6.5		1.0	2.0	3.0	4.0	5.0	6.0
		6514	6524	6534	979	559	V
Height	4.0 5.5 7.0	UDA 6514H UDA 6515H UDA 6517H	UDA 6524H UDA 6525H UDA 6527H	UDA 6534H UDA 6535H UDA 6537H	UDA 6544H UDA 6545H UDA 6547H	UDA 6554H UDA 6555H UDA 6557H	UDA 6564H UDA 6565H UDA 6567H





Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø4.0		4014	4024	1034	404	4054	4064
Height	4.0 5.5	UDAN 4014N UDAN 4015N	UDAN 4024N UDAN 4025N	UDAN 4034N UDAN 4035N	UDAN 4044N UDAN 4045N	UDAN 4054N UDAN 4055N	UDAN 4064N UDAN 4065N
	7.0	UDAN 4017N	UDAN 4027N	UDAN 4037N	UDAN 4047N	UDAN 4057N	UDAN 4067N
Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø4.5			60	34	4544	4554	4564
		4514	452		V	V	W
Height	4.0 5.5	UDAN 4514N UDAN 4515N	UDAN 4524N UDAN 4525N	UDAN 4534N UDAN 4535N	UDAN 4544N UDAN 4545N	UDAN 4554N UDAN 4555N	UDAN 4564N UDAN 4565N
neight	7.0	UDAN 4515N UDAN 4517N	UDAN 4527N	UDAN 4537N	UDAN 4547N	UDAN 4557N	UDAN 4567N



Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø4.5		4534	4524	4534	454	PESS	4564
	4.0	UDA 4514N	UDA 4524N	UDA 4534N	UDA 4544N	UDA 4554N	UDA 4564N
Height	5.5	UDA 4515N	UDA 4525N	UDA 4535N	UDA 4544N	UDA 4555N	UDA 4565N
neignt							
Height	7.0	UDA 4515N UDA 4517N	UDA 4525N LIDA 4527N	UDA 4535N UDA 4537N	UDA 4545N UDA 4547N	UDA 4555N UDA 4557N	UDA 4

Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø5.5		255	\$524	\$5534	\$554	15555	A358
Height	4.0 5.5 7.0	UDA 5514N UDA 5515N UDA 5517N	UDA 5524N UDA 5525N UDA 5527N	UDA 5534N UDA 5535N UDA 5537N	UDA 5544N UDA 5545N UDA 5547N	UDA 5554N UDA 5555N UDA 5557N	UDA 5564N UDA 5565N UDA 5567N
Cuff D Ø6.5		1.0	2.0	3.0	4.0	5.0	6.0
		6514	6524	6334	999	1559	181
Height	4.0 5.5 7.0	UDA 6514N UDA 6515N UDA 6517N	UDA 6524N UDA 6525N UDA 6527N	UDA 6534N UDA 6535N UDA 6537N	UDA 6544N UDA 6545N UDA 6547N	UDA 6554N UDA 6555N UDA 6557N	UDA 6564N UDA 6565N UDA 6567N



ImpressionPost

Pick Up Type

(Unit: mm)



Height		11.0				15.0			
D Ø3.7	Hex	Non-Hex			Hex	Non-Hex			
					A	A			
	UIPPN 37 11 H	UIPPN 3711N	UIPPNS 15	UIPPNS 21	UIPPN 37 15 H	UIPPN 37 15 N	UIPPNS 19	UIPPNS 25	



Height		11.0			15.0				_
	Hex	Non-Hex			Hex	Non-Hex			
D Ø4.0 D Ø4.5 D Ø5.5	UIPP 4011H UIPP 4511H UIPP 5511H	UIPP 4011N UIPP 4511N UIPP 5511N	UIPPS 16	UIPPS 20	UIPP 4015H UIPP 4515H UIPP 5515H	UIPP 4015N UIPP 4515N UIPP 5515N	UIPPS 22	UIPPS 26	



Height	11.0			15.0			
D Ø3.7	Hex	Non-Hex		Hex	Non-Hex	A	
				\$ 17 E			
	UIPTN 37 11 H	UIPTN 3711N	UIPTNS 16	UIPTN 3714H	UIPTN 3714N	UIPTNS 19	



Height		11.0		15.0		
	Hex	Non-Hex		Hex	Non-Hex	4
		40				
D Ø4.0	UIPT 40 11 H	UIPT 4011N	UIPTS 16	UIPT 40 14 H	UIPT 40 14 N	UIPTS 19
D Ø4.5	UIPT 45 11 H	UIPT 4511N		UIPT 45 14 H	UIPT 45 14 N	
D Ø5.5	UIPT 5511H	UIPT 5511N		UIPT 55 14 H	UIPT 5514N	

Fixture

Analog





UDAG 35















Torque: Narrow[20 Ncm]/Regular[30 Ncm]
Packing unit : Angled Ab timent_Non Hex Type + Noutment Screw
Order Code : ex) U17AA 4537N











Cuff	2.0	3.0	4.0	5.0
D Ø4.0		En.	/1	//
			#	\forall
	U17AAN 4027H	U17AAN 4037H	U17AAN 4047H	U17AAN 4057H



Cuff	2.0	3.0	4.0	5.0
D Ø4.5	4			
	U17AA 4527H	U17AA 4537H	U17AA 4547H	U17AA 4557H
Cuff	2.0	3.0	4.0	5.0
D Ø5.5	#			
	U17AA 5527H	U17AA 5537H	U17AA 5547H	U17AA 5557H



Cuff	2.0	3.0	4.0	5.0
D Ø4.0				
	U17AAN 4027N	U17AAN 4037N	U17AAN 4047N	U17AAN 4057N



Cuff	2.0	3.0	4.0	5.0
D Ø4.5			h-	l a
	/1	//		//
	4	4	44	6.0
	M	W	W	W
	U17AA	U17AA	U17AA	U17AA
	4527N	4537N	4547N	4557N
Cuff	2.0	3.0	4.0	5.0
D Ø5.5				P =
	10	/1	/1	//
	//			44





U17AA 5537N





URIS Implant System


Milling Abutment Нех Туре

Milling Abutment












Cuff

D Ø5.0



1.0



3.0



























Abutment

Screw







UAS 16H

UAS 20H

Abutment Level ImpressionFlowchart








DirectAbutment



Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø4.0		4014	4024	4034	4044	4054	1900
Height	4.0 5.5 7.0	UDAN 4014 UDAN 4015 UDAN 4017	UDAN 4024 UDAN 4025 UDAN 4027	UDAN 4034 UDAN 4035 UDAN 4037	UDAN 4044 UDAN 4045 UDAN 4047	UDAN 4054 UDAN 4055 UDAN 4057	UDAN 4064 UDAN 4065 UDAN 4067
Cuff D Ø4.5		1.0	2.0	3.0	4.0	5.0	6.0
5 54.3		4514	4524	4534	4544	4554	P95F
Height	4.0 5.5	UDAN 4514 UDAN 4515	UDAN 4524 UDAN 4525	UDAN 4534 UDAN 4535	UDAN 4544 UDAN 4545	UDAN 4554 UDAN 4555	UDAN 4564 UDAN 4565
	7.0	UDAN 4517	UDAN 4527	UDAN 4537	UDAN 4547	UDAN 4557	UDAN 4567



Cuff		1.0	2.0	3.0	4.0	5.0	6.0
D Ø4.5							E-1
			500	\$534	4544	4554	4564
		4514	4524	15			- 10
		1	W	W	W	W	W
		I	T	T		I	I
	4.0	UDA 4514	UDA 4524	UDA 4534	UDA 4544	UDA 4554	UDA 4564
Height	5.5	UDA 4515	UDA 4525	UDA 4535	UDA 4545	UDA 4555	UDA 4565
	7.0	UDA 4517	UDA 4527	UDA 4537	UDA 4547	UDA 4557	UDA 4567



Height 4.0 5.5 7.0 D Ø4.0 UDAC 4004 UDAC 4005 UDAC 4007

Height	4.0	5.5	7.0
D Ø4.5			
	UDAC 450 4	UDAC 4505	UDAC 450 7

Height	4.0	5.5	7.0
D Ø5.5			
	UDAC 5504	UDAC 550 5	UDAC 550 7

Height	4.0	5.5	7.0
D Ø6.5			
	UDAC 650 4	UDAC 650 5	UDAC 650 7

ImpressionCap

Diameter	Ø4.0	Ø4.5	Ø5.5	Ø6.5
		*	*	書
	UDAIC 40	UDAIC 45	UDAIC 55	UDAIC 65

Burn Out Coping

(Unit: mm)

Single	Ø4.0	Ø4.5	Ø5.5	Ø6.5
	UDABC 40S	UDABC 45S	UDABC 55S	UDABC 65S
			-	
Bridge	Ø4.0	Ø4.5	Ø5.5	Ø6.5
	UDABC 40B	UDABC 45B	UDABC 55B	UDABC 65B

69

Direct Abutment Analog



Direct Abutment

Analog

(One min)		
4.0	5.5	7.0
(A)	500	107
40	Ü	Ü
	٥	4
UDAG 40 4	UDAG 40 5	UDAG 40 7
4.0	5.5	7.0
	63	п
454	455	457
TT .	П	
ĕ	ä	ä
UDAG 45 4	UDAG 45 5	UDAG 45 7
4.0	5.5	7.0
		B
554	555	255
T	T	T
X.	X.	X.
3000	111	107
UDAG 554	UDAG 55 5	UDAG 55 7
4.0	5.5	7.0
	63	
654	559	159
	4.0 UDAG 404 4.0 UDAG 454 4.0 UDAG 554	4.0 5.5 UDAG 404 UDAG 405 4.0 5.5 UDAG 454 UDAG 455 UDAG 554 UDAG 555



UDAG 654







UDAG 65**5**


UDAG 65**7**

URIS Implant System

Abutment Level Impression

Multi-unit Abutment & Components

Flowchart







Cuff	1.0	2.0	3.0	4.0	5.0	6.0
D Ø5.0					43	//
	//	A	4			17
	V	V	V	V	V	V
	8	#		#	#	#
	UMAN 5010	UMAN 5020	UMAN 5030	UMAN 5040	UMAN 5050	UMAN 5060



)						
Cuff	1.0	2.0	3.0	4.0	5.0	6.0
D Ø5.0						//
		43				
	()		W	W	W	- U
						¥
	UMA 5010	UMA 5020	UMA 5030	UMA 5040	UMA 5050	UMA 5060

Multi-Unit

butment	Angled Type

(Unit: mm)



2.5 / 17°	3.0 / 17°	4.0 / 17°	3.5 / 29.5°	4.0 / 29.5°	5.0 / 29.5°
/1	12	dA	di	el A	
U17MAN 50 25 H	U17MAN 50 30 H	U17MAN 50 40 H	U30MAN 50 35 H	U30MAN 50 40 H	U30MAN 50 50 H
	4	U17MAN U17MAN	U17MAN U17MAN U17MAN	U17MAN U17MAN U30MAN	U17MAN U17MAN U30MAN U30MAN U30MAN



	Cuff/Angle	2.5 / 17°	3.0 / 17°	4.0 / 17°	3.5 / 29.5°	4.0 / 29.5°	5.0 / 29.5°
	D Ø5.0						
		1	da	n/A	de	er A	
		1					
		U17MA	U17MA	U17MA	U30MA	U30MA	U30MA
		50 25 H	50 30 H	50 40 H	50 35 H	50 40 H	50 50 H
N							
W							
	Cuff/Angle	2.5 / 17°	3.0 / 17°	4.0 / 17°	3.5 / 29.5°	4.0 / 29.5°	5.0 / 29.5°
	D Ø5.0						
		1	dA	al A	ed A		
		U17MAN	U17MAN	U17MAN	U30MAN	U30MAN	U30MAN
		50 25 N	50 30 N	50 40 N	50 35 N	50 40 N	50 50 N
R							
	Cuff/Angle	2.5 / 17°	3.0 / 17°	4.0 / 17°	3.5 / 29.5°	4.0 / 29.5°	5.0 / 29.5°
	D Ø5.0						
		<i>J</i> 1	1	dA	di	de	
		1	4	W	W	W	W
		U17MA	U17MA	U17MA	U30MA	U30MA	U30MA
		50 25 N	50 30 N	50 40 N	50 35 N	50 40 N	50 50 N

Multi-Unit Abutment Screw









Cuff/Angle	3.0 / 17°	4.0 / 17°	5.0 / 17°	4.0 / 29.5°	5.0 / 29.5°	6.0 / 29.5°
D Ø5.0		/1	1		10	da
			W.			
	U17MUAN 5030H	U17MUAN 5040H	U17MUAN 5050H	U30MUAN 5040H	U30MUAN 5050H	U30MUAN 5060H
R						
Cuff/Angle	3.0 / 17°	4.0 / 17°	5.0 / 17°	4.0 / 29.5°	5.0 / 29.5°	6.0 / 29.5°
D Ø5.0			/=		<i>k</i> .	Ja
	*			W.	W.	
	U17MUAR	U17MUAR	U17MUAR	U30MUAR	U30MUAR	U30MUAR

5050H

5040H

5050H

5060H

Abutment

Screw



5030H

5040H

UAS 16H



UAS 20H

Multi-Unit Components



Scale 1:1.5

UMTIC 50N

Multi-Unit Cylinder Screw

UMCS 16

Burn Out Cylinder



Multi-Unit

Temporary Cylinder







Impression

Post










Abutment Level Impression **Overdenture - Ball Abutment**

Flowchart









Cuff	1.0	2.0	3.0	4.0	5.0	6.0
D Ø3.5		8				
	UBAN 3510	UBAN 3520	UBAN 3530	UBAN 3540	UBAN 3550	UBAN 3560





Components

Ball Abutment Analog

Retainer Cap



UBSC 35

UBSO 35

Retainer



GOR 4515K

O-ring


UBAG 35

URIS Implant Syster

Abutment Level Impression

Overdenture - T:LOC Abutment

Flowchart









T:LOC **Titanium Cap** Black Processing **Replacement Male**

TLC-TC5423







Retention **Replacement Male**

Block Out **Spacer**





TLC-RRM478B TLC-RRM47P TLC-RRM47T





TLC-PRM56K

T:LOC Component



Impression Coping

T:LOC **Analog**




UBAG 35



UBAG 35

Digital Impression & Components Flowchart







GH/Height	GH 1 / 3.5	GH 1 / 5.5	GH 2 / 3.5	GH 2 / 5.5
D Ø4.5				m
	F	H	A	H
	7	7	T	7
	UOR-R13S	LIOR-B15S	LIOR-B23S	LIOR-B25S



GH/Height	GH 1 / 3.5	GH 1 / 5.5	GH 2 / 3.5	GH 2 / 5.5
D Ø4.0			Ī	
	UON-B13B	UON-B15B	UON-B23B	UON-B25B



GH/Height	GH 1 / 3.5	GH 1 / 5.5	GH 2 / 3.5	GH 2 / 5.5
D Ø4.5			Ī	
	UOR-B13B	UOR-B15B	UOR-B23B	UOR-B25B

URIS Scan Body



FixtureDigital Analog





Multi-Unit **Base**

Multi-Unit **Scan Body**









UMSB 50

Multi-Unit **Digital Analog**










UMCSB 50







Fixture Packaging



Cover Screw

Cover Screw

intended for use in partially or fully edentulo multiple-unit restorations. The surface is SLA (Sandblasted, Large grit and Acid etched) treated and is provided sterile. It consists of two implant lines, the OMNI and the OMNI Tapered, with corresponding cover screws, healing abutments and prosthetic abutments. The OMNI Tapered implant has a tapered wall with a single thread design. The OMNI is straight walled with smaller threading at the coronal end, and bigger threading at the apical end. Both implant lines have two platform sizes, Narrow (\emptyset 3.5 mm) and Regular (\emptyset 4.0 – \emptyset 6.5 mm). Both implant lines share the

URIS OMNI System fixtures are dental implants made of Unalloyed Titanium, grade 4 (ASTM F67)

URIS OMNI System

Instruction for Use

Ø 3.5 x 8.5. 10. 11.5. 13. 14.5mm (L) Ø 4.0 x 7, 8.5, 10, 11.5, 13, 14.5mm (L) Ø 4.5 x 7, 8.5, 10, 11.5, 13, 14.5mm (L) Ø 5.0 x 7, 8.5, 10, 11.5, 13, 14.5mm (L) Ø 5.5 x 7, 8.5, 10, 11.5, 13, 14.5mm (L) Ø 6.0 x 7. 8.5. 10mm (L) Ø 6.5 x 7, 8.5, 10mm (L).

Device Description

URIS Prosthetic System is made of titanium alloy (Ti-6Al-4V ELI) intended for use as an aid in prosthetic restoration. It consists of Cover screw, Healing abutment, Direct Abutment, Basic abutment, Angled abutment, Milling abutment, Temporary abutment,

and Abutment screw. The surface of cover screw and healing abutment are anodized in yellow

Device Component	Diameters (Ø)	Lengths	Angulation
Cover Screw	2.78/3.48mm	4.875/5.375mm	-
Healing Abutment	4.0/4.5/5.5/6.5/7.5mm	Cuff Height: 1.0mm~5.0mm	-
Direct Abutment	4.0/4.5/5.5/6.5/7.5mm	Cuff Height: 1.0mm~6.0mm	-
Basic Abutment	4.0/4.5/5.5/6.5mm	Cuff Height: 1.0mm~6.0mm	-
Angled Abutment	4.0/ 4.5/5.5mm	Cuff Height: 2.0mm~5.0mm	17°
Milling Abutment	4.0/5.0/6.0/7.0mm	Hex Type: 14.1/14.85mm Non-Hex Type: 13.9/14.85mm	-
Temporary Abutment	3.7 / 4.3mm	Cuff Height: 1.0mm~3.0mm	-
Abutment Screw	1.9/2.3mm	7.2/7.7mm	-

Fixtures and cover screw are provided sterile and other prosthetics are provided non-sterile. All non-sterile products must be sterilized by end users before use

Indications for Use

URIS OMNI System is indicated for use in partially or fully edentulous mandibles and maxillae, in overdenture restorations, and final or temporary abutment support for fixed bridgework. It is

Instructions for operation and use

- 1) Before clinical use, the clinician must be well acquainted with the surgical procedure of the product, and has to inform the patient about the limitations of the implant system. The patient should also be well aware of any functional and aesthetic limitations of the implant.
- 2) Because proper selection and fixation of the implant are closely related to the life-span of the implant, the clinician must follow the indications, contraindication, cautions and recommendations.
- 3) Handling procedures must be followed in order to prevent potential damage to the implant. Damage to the implant and/or patient may occur without careful review of the patient's
- condition and establishment of proper diagnosis and restorative plans 4) The clinician must select the appropriate device based on careful review of the patient's X-ray
- 5). Check the products' expiration date and condition of the packaging for any visible damages 6) Since the product is packaged aseptically, do not use if the packaging is damaged or torn
- 7) Be sure to properly maintain the hygiene standards and preparatory state of the surgical instruments in order to prevent the use of contaminated instruments which may lead to
- 8) Inspect for any foreign-material before use.

B. Instructions and procedural sequence

During the diagnosis and planning, you must exclude any patient with local lesions or other contraindications and choose candidates who have proper bone condition to undergo implant surgery. perioral area thoroughly. After proper draping, perform local anesthesia and make an incision on the implant site and form a flap. Expose the implant site sufficiently and proceed to implant surgery. (1) Implant site preparation

To implant the fixture, various drills are used in sequence for site preparation during os To place the fixture accurately in the selected site, a hole must be made according to the size of the artificial dental prosthesis, using the respective instruments (drilling, tapping). Rotatory speed during these procedures must be adjusted taking the recipient bone condition and type of equipment used into consideration. The maximum permissible rotatory speed for the drill is generally 1,000~1,500rpm and 20~30rpm for the tap drill. The procedure should be performed using adequate normal saline to reduce the generation of heat on the bone tissue.

Pick up the fixture from the sterile vial using the Fixture Driver and Adapter and place the fixture into the osteotomy. Install the fixture at low speed (25 rpm) under profuse irrigation and the maximum torque set at 45 Ncm. Allow the implant to work its way into the osteotomy. Avoid applying

NOTE: The final recommended torque at seating should be 20~40Ncm for the URIS OMNI System. Excessively high insertion torque may cause necrosis of the peri-implant bone in the receiving site which may result in implant failure.

After the fixture has been placed, attach the cover screw using a driver below 10Ncm torque. Make sure there are no foreign bodies inside and suture the operation site.

(4) Connecting the abutment

Osseo-integration of the fixture requires 3~4 months for the mandible and 6~8 months for the maxilla. After this period, expose the implant and connect the healing abutment to

After a healing period of between 2~4 weeks, connect the impression post to obtain an impression and manufacture a dental mockup. Deliver the final prosthesis.

- 1) The operation must be performed by a well-trained, qualified dental specialist.
- 2) While performing the osteotomy, you must follow the procedure outlined in the catalog and the fixture should be adequately implanted.
- 3) Ensure that the soft tissue does not interfere with the connection between the fixture and prosthesis by verifying complete and proper seating.
 4) All instruments and tooling used during the procedure must be maintained in good condition
- and care must be taken that instrumentation does not damage implants and/or other components. Therefore inspect the condition of the instruments before every operation.
- 5) The product is provided sterile via gamma ray sterilization therefore it recommended to be opened prior to immediate use.
- 6) If the package has been damaged, discard the product since the aseptic condition has been

- A. In cases with insufficient bone tissue where severe bone resorption is predicted. Or if there is insufficient remaining bone for early-fusion in the proximal tooth extraction wound.
- B. Disorder in mastication or functional relation
- C. Pathologic condition of the alveolar bone D. Prior radiotherapy on jawbone
- F. Pathologic change of oral mucosa (vitiligo, lichen planus, stomatitis)
- G. Macroglossia
- H. If vital anatomical structures are nearby
- I. Cellulitis in surrounding soft tissues
 J. If there are not sufficient soft tissues or its condition is poor
- 2) Transient contraindications
- B. Pregnancy
- . Temporary effect of specific drugs (anticoagulant, immune-suppressant)
- D. Mental, physical fatique
- 3) Psychological contraindication
- A. Poor compliance B. Alcohol or other substance abuse
- . Neurosis, psychosis patient
- D. Troublesome patient
- A. General/nutritional condition age (obesity, cachexia, 5year survival rate)
- C. Metabolic disorder (pubertal diabetes, overt hyperglycemia (>300mg/dl))
- . Hematologic disorder (disorder of RBC, WBC, coagulation)
- F. Cardiovascular diseases (artherosclerosis, overt hypertension (>300mmHg))
- Metabolic disorder of skeletal system (osteomalacia, Paget's disease, menopausal osteoporosis) G. Connective tissue disease (dermatosclerosis, rheumatoid arthritis)
- H. Implant as potential infection focus (prosthetic valve, bacterial endocarditis)

- 1) Implant operation should be performed by skilled dental surgeon because mishandled procedure may damage the implant or recipient bone
- 2) Implant is not to be recycled and it should be used for its original purpose 3) Damaged or mishandled implant should be removed
- 4) Inappropriate implant selection and improper implantation site or unstable fixation may shorten the life-span of the implant

- 6) Handle the implant carefully to prevent any damage or deformation 7) Warning: Small diameter implants and angled abutments are not recommended for the molar region of the mouth.

(4) POTENTIAL ADVERSE EFFECTS AND COMPLICATIONS

General complications after intraoral implant surgery include local hemorrhage, edema and hematoma. Transient loss of taste sense and masticatory function may occur. Additionally, following complications may develop:

- latrogenic trauma of surrounding tissues (lower alveolar nerve injury or sensory change, injury or hemorrhage in maxillary sinus or nasal cavity)
- Wound dehiscence on sutured site
- Delayed recovery, edema due to anesthesia
- Mucositis around implant due to insufficient adhesive soft tissue
- Incomplete implant placement due to insufficient bone removal or overt compression
- General hypersensitivity reaction

MR Statement

The URIS OMNI System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of URIS OMNI System in the MR environment is unknown.

94

All dental implants (fixture) and cover screw are supplied sterile and are labeled "STERILE". All products sold sterile are for single-use before the expiration date printed on the product label. Do not use sterile products if the packaging has been damaged or previously opened. Do not re-sterilize.

End User Sterilization Information

All prosthetic abutments are provided non-sterile and must be sterilized before use. To correctly sterilize the products, use a steam sterilizer with pre-vacuum process, at a temperature of steam sterilizer at 132° C for 4 minutes, wrap and dry 20 minutes with a validated cycle according to the standard ISO 17665-1 following the autoclave manufacturer instructions.

	Pre-Vacuum Autoclave
Temperature	132° C
Exposure Time	4 minutes
Dry Time	20 minutes

Note: The validated procedures require the use of FDA-cleared sterilizers, sterilization trays, sterilization wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended. The health care facility should monitor the sterilize for the facility according to an FDA recognized sterility assurance standard such as ANSI/AAMI ST79.

 $\begin{tabular}{ll} \textbf{Storage} \\ \textbf{The product has to be stored in its original package in a dry place at room temperature.} \\ \end{tabular}$

Handling

- This product is a disposable sterilized medical instrument and should therefore not be reused.
- Packing must be opened prior to surgery in a clean area.
 Discard if wrapping has been opened even if product is unused.
- Do not use the product, if the shelf life has expired. Opened products cannot be returned to the manufacturer or distributor.

TruAbutment Korea Co., Ltd. 397, Seokcheon-ro,

Gyeonggi-do,14449,Rep.of Korea

Ojeong-gu, Bucheon-si,

Phone:+82 32 678 4688

www.urisimplants.com

Manufacturer or distributor has no responsibility for products re-sterilized by users.

LABELING SYMBOLS

Symbols may be used on some international package labeling for easy identification.

Do Not Reuse

Use By Date

LOT **Batch Code**

 \sim Date Of Manufacture

Sterilized Using Irradiation STERILE R

Catalogue Number REF

Caution, Consult Accompanying Documents

Consult Instructions For Use



i

Rx Only Prescription Only

Do Not Resterilize

Manufacturer

Keep Dry

Keep Away From Sunlight

Do Not Use If Package Is Damaged

Temperature Limit

Rev.3 / 2019.06



URIS Implant System

96

Device Description

URIS OMNI System fixtures are dental implants made of Unalloyed Titanium, grade 4 (ASTM F67) intended for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations. The surface is SLA (Sandblasted, Large grit and Acid etched) treated and is provided sterile. It consists of two implant lines, the OMNI and the OMNI Tapered, with corresponding cover screws, healing abutments and prosthetic abutments. The OMNI Tapered implant has a tapered wall with a single thread design. The OMNI is straight walled with smalle threading at the coronal end, and bigger threading at the apical end. Both implant lines have two platform sizes, Narrow (\emptyset 3.5 mm) and Regular (\emptyset 4.0 – \emptyset 6.5 mm). Both implant lines share the

Ø 3.5 x 8.5, 10, 11.5, 13, 14.5mm (L) Ø 4.0 x 7, 8.5, 10, 11.5, 13, 14.5mm (L) Ø 4.5 x 7, 8.5, 10, 11.5, 13, 14.5mm (L) Ø 5.0 x 7, 8.5, 10, 11.5, 13, 14.5mm (L) Ø 5.5 x 7, 8.5, 10, 11.5, 13, 14.5mm (L) Ø 6.5 x 7, 8.5, 10mm (L).

URIS Prosthetic System is made of titanium alloy (Ti-6Al-4V ELI) intended for use as an aid in prosthetic restoration. It consists of Cover screw, Healing abutment, Direct Abutment, Basic abutment, Angled abutment, Milling abutment, Temporary abutment, and Abutment screw. The surface of cover screw and healing abutment are anodized in yellow

Device Component	Diameters (Ø)	Lengths	Angulation
Cover Screw	2.78/3.48mm	4.875/5.375mm	-
Healing Abutment	4.0/4.5/5.5/6.5/7.5mm	Cuff Height: 1.0mm~5.0mm	-
Direct Abutment	4.0/4.5/5.5/6.5/7.5mm	Cuff Height: 1.0mm~6.0mm	-
Basic Abutment	4.0/4.5/5.5/6.5mm	Cuff Height: 1.0mm~6.0mm	-
Angled Abutment	4.0/ 4.5/5.5mm	Cuff Height: 2.0mm~5.0mm	17°
Milling Abutment	4.0/5.0/6.0/7.0mm	Hex Type: 14.1/14.85mm Non-Hex Type: 13.9/14.85mm	-
Temporary Abutment	3.7 / 4.3mm	Cuff Height: 1.0mm~3.0mm	-
Abutment Screw	1.9/2.3mm	7.2/7.7mm	-

Fixtures and cover screw are provided sterile and other prosthetics are provided non-sterile. All non-sterile products must be sterilized by end users before use.

Indications for Use

URIS OMNI System is indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations including; cemented retained, screw retained, or overdenture restorations, and final or temporary abutment support for fixed bridgework. It is intended for delayed loading.

Restorative Components:

The abutments are used to restore a dental implant, acting like the base for the prosthesis. They are available in different shapes and sizes to respond to different needs.

It should maintain at least 4mm from the abutment platform to avoid damaging the abutment screw:

There are four types of titanium abutments available:

- Basic Abutments: None of the Basic abutments are to be used as titanium base abutments or as part of a hybrid abutment (e.g. part of a two-piece abutment). Basic abutments are full size abutments and are to be used straight. No angular correction or divergence is allowed by any additional copings, or modifications.
- Angled Abutments: In general, they are angled at 17°, which can be adapted to the majority of
- Milling Abutments: Intended to be milled by hand and are intended to be used straight. No angular correction or divergence is allowed.
- Temporary Abutments: Intended for use in conjunction with the fixture in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit cement-retained restorations. Maximum duration for use of Temporary Abutment is less than six months.
- Temporary restorations should be out of occlusion.

The screws are made from Ti-6Al-4V ELI(ASTM F136), recommended for its biocompatibility, its mechanical strength and hardness. It serves to attach the abutment or prosthesis to the implant (clinical screw) or to the laboratory analogue (laboratory screw).

Recommendations for its specific use include:

The screws are for single-use only. It is not recommended to use the screws again after their removal, not even in the laboratory, due to the possible deterioration of their behavior. It is vitally important to not use clinical case screws that have been previously used in a dental laboratory. It is important to verify the compatibility of the implant model to be used. You should avoid causing any damage around the area where the implant is connected, so care must be taken if carving or machining in this area. Radiography is recommended in the height of the junction of the union with the perpendicular axis of said union, once the implant is fixed, for verification

Only the implant manufacturer's recommended torque is to be used. Ncm Abutments (Narrow Connection)

20	Basic Abutment, Direct Abutment, Angled Abutment, Milling Abutment, Temporary Ab
Ncm	Abutments (Regular Connection)
5~10	Cover Screw, Healing Abutment
20	Temporary Abutment
25~30	Basic Abutment, Direct Abutment, Angled Abutment, Milling Abutment, TruBase

1. Abutments

The instructions given are insufficient if used as the only reference for the use of the cited components. These elements should only be inserted by dentists who have been fully trained in the insertion of dental implants. The use of these products without any prior specific knowledge can lead to component failure and may require implant removal. The safety of our products is guaranteed only when they are used exclusively by trained professionals. Read the instructions carefully on the labels of the products, where you will find the basic guidelines Keep a record of the products used in the patient's personal medical booklet, stating the name of the product, the reference number, and the lot number. Please inform URIS Implants of any defects or complications related to any of its products.

All URIS OMNI System products are solely for single use. To reuse the single-use products may lead to a possible deterioration of the characteristics of the product, which in turn can lead to an elevated risk in gum or tissue infection and deterioration in the patient's health. In general, implant component's placement and prosthetic design must accommodate individual patient conditions. In case of bruxism or unfavorable jaw relationships reappraisal of the treatment option may be considered.

There is a risk of accidental inhalation and/or ingestion of the products when they are used, therefore it is necessary to carefully hold onto the products in case of intraoral applications. The patient should be made aware of any limitations in his/her treatment, and the need for maintenance, for example, the need to seek medical assistance if any symptoms or side effects arise It should be recommended to the patient to conduct regular dental check-ups for maintenance

of the URIS OMNI System products.
The products are not sterilized when sold, and therefore, it is recommended to clean and sterilize the products before their use.

* Warning: Small diameter implants and angled abutments are not recommended for the

Contraindications:

- Medically unfit for an oral surgical procedure.
 With inadequate bone volume unless an augmentation procedure can be considered.
- In whom adequate sizes, numbers or desirable position of implants are not reachable to achieve safe support of functional or eventually parafunctional loads.
- Allergic or hypersensitive to titanium alloy (grade 5).

The URIS OMNI System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment.

The safety of URIS OMNI System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

End User Sterilization Information

All prosthetic abutments are provided non-sterile and must be sterilized before use. To correctly sterilize the products, use a steam sterilizer with pre-vacuum process, at a temperature of steam sterilizer at 132° C for 4 minutes, wrap and dry 20 minutes with a validated cycle according to the standard ISO 17665–1 following the autoclave manufacturer instructions.

	Pre-Vacuum Autoclave
Temperature	132° C
Exposure Time	4 minutes
Dry Time	20 minutes

Note: The validated procedures require the use of FDA-cleared sterilizers, sterilization travs sterilization wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended. The health care facility should monitor the sterilize for the facility according to an FDA recognized sterility assurance standard such as ANSI/AAMI ST79.

The product has to be stored in its original package in a dry place at room temperature.

LABELING SYMBOLS

Symbols may be used on some international package labeling for easy identification.

Do Not Reuse Use By Date LOT **Batch Code** Date Of Manufacture Non-Sterile Catalogue Number REF **Caution, Consult Accompanying Documents** Manufacturer Consult Instructions For Use Do Not Use If Package Is Damaged Keep Dry Keep Away From Sunlight

Rev.3 / 2019.06



TruAbutment Korea Co., Ltd.

397. Seokcheon-ro. Ojeong-gu, Bucheon-si, Gyeonggi-do,14449,Rep.of Korea

Phone:+82 32 678 4688 www.urisimplants.com

General Principles

of Surgical Tool Management

- Caution Wrong cleansing and sterilizing process cause corrosion and damage to the tools and if used directly, it may be the cause of second infection.

 The recommended number of use of a drill is 20~30 times based on the bone status, and it must be replaced if the blade has been damaged or transformed
- Caution If damaged drill is used, heat necrosis may occur

 3. When managing the surgical tool, one must wear a mask and a glove to prevent infection.

- 1. To prevent contaminants such as blood, tissue cell or bone residue from attaching to the surface of the instruments, the instruments must be immersed in an antiseptic solution right after use.
- 2. When using antiseptic solution, to prevent corrosion or bronzing, one must follow directions given by the manufacturer of the concentration of the antiseptic and the duration of the instrument

 $Check \ Concentration: Completely \ liquify \ the \ concentrate \ before \ placing \ the \ instruments \ in \ the \ antiseptic \ solution.$ Immersion Duration: The instruments must not be immersed more than a day.

- 3. The instruments must be fully immersed in the antiseptic solution.
 4. To decrease in sterilizing power and to prevent corrosion, the antiseptic solution must be replaced every day.

To prevent protein from clotting in 45 degrees temperature Celsius, the instruments must be rinsed in running cold water.

Cleanse the instruments right after preliminary rinse

- 1. Must only use antiseptic solution that is FDA and CE approved, and you must follow the manufacturer's directions
- When cleansing metal instruments, corrosion free antiseptic solution and cleansing product use is recommended
- 3. For safety, one must always wear personal protection gear such as gloves, glasses, and masks
- 4. The user has an obligation to be responsible for the sterilization and management of the instrument
- With repetition of cleansing, the life expectancy of all instruments will decrease. If the instruments show corrosion, transformation or discoloring of the marking area, it means that they have exceeded
- Product with a disposable mark cannot be reused
- Tungsten carbide burs, plastic composition and NiTi instruments can be damaged with hydrogen peroxide and aluminum material instruments can be damaged by caustic soda solution.
- Do not use acid solution (pH < 6) and alkaline solution (pH > 8).

After use, if the contaminants such as residual bone or blood stain are not completely removed, it may lead to corrosion; therefore all separable instruments must all be disassembled before the cleansing process.

Cleanse / Dry

1. Contaminants must be completely removed using a soft brush.

Do not use a wire brush or stainless material brush, and do not put too much pressure.

2. Immerse the products in the antiseptic solution of their characteristics and clean with an ultrasonic cleaner. However, do not cleanse the different materials together. Also, when immersing the instruments in the ultrasonic cleaner, make sure that the instruments do not touch each other.

3. Make sure that debris is not seen with the naked eye.

- Products that are fractured or transformed must be discarded.
- One should follow the recommendations for the level of concentration or the length of time provided by the manufacturer.
- The antiseptic solution must not include aldehyde, di- or tri-ethanolamines component to control the corrosion.

 4. After cleaning, the products must be rinsed with distilled water or deionized water for at least a minute. If the antiseptic solution contains corrosion inhibitor, rinsing before placing in the sterilizer
- 5. To prevent corrosion or water stain on the instruments, completely dry with a dryer or filtered compressed air 6. To prevent corrosion, decrease in sterilizing power, and contamination, antiseptic must be supplemented even

If the instruments are not properly rinsed, residue is left behind, or is not properly dried, the sterilization process might discolor or corrode the instruments, and therefore the whole process must be gone through again

Corrosion may start if debris such as blood stain or bone residue is not completely removed. They must be cleansed right after use and the debris must be completely removed when cleaning.

Check on the instruments for faults (fracture, transformation, or corrosion). If necessary, assemble the instruments.

Contaminated instruments must be cleansed or disinfected. Transformations that may affect the safety, performance or tolerance of the instruments; in other words; bent, damaged (fractured or corroded), or faulty products (discoloration of marking area or loss) must be destroyed.

Check on the dry status of the instruments and pack in the sterilized wrapping paper.

2. On the sterilized wrapping paper, attach a direction tape to check the date of sterilization. Check on the expiration date on the sterilized wrapping paper. Wrapping paper must be able to withstand up to 141 degrees that coincides with the EN ISO 11607.

13. The product is packaged cleaned and should be sterilized before its use. To correctly sterilize the products, use a steam sterilizer with pre-vacuum process at a temperature of steam sterilizer at 132° C for 4 minutes, dry for 20 minutes with a validated cycle according to the standard ISO 17665-1 following the auto-clave manufacturer instructions.

2. Instruments and plastic components must be sterilized based on their packaging label.

- Sterilizer must coincide with the requirements of EN 13060 and EN285.
 Sterilization process must regard the ISO 11607.
- One must follow the sterilization process and maintenance process of the sterilizer provided by the manufacturer.
 Efficiency management (proper packaging, no humidity level and sterilization dashboard).

- The products must not touch the inner part of the sterilization equipment, and the sterilization degree must be lower than 150
- The products that were not properly cleansed or dried may generate corrosion. If they were not cleansed, not properly dried, or has been corroded, separate them from the rest or remove the faults. (Do not sterilize the corroded instruments with the noncorroded products together)
- For sterilization, use only salt-free water or distilled water for the solution. (Do not use tap water)
- Check if the instruments are fully dried and do not leave them in a place with high moisture

nts must be stored in a sterilized container in a dry and clean environment. If the packaging is opened or damaged, we cannot guarantee the instruments' sterilization status.