**SPECIAL PRECAUTIONS:**

1. Always exercise extreme caution when operating your Motor
2. Always inspect your Motor thoroughly before each use.
3. Always use PPE (Proper Protective Equipment) when operating and cleaning your Motor.
4. See your Occupational Health and Safety Act or other safety ordinances for specific instructions.
5. This is a precision instrument. Do not drop or apply excessive force when handling.
6. If product is dropped your dealer should be contacted, it is recommended that you do not use.
7. Contact CURION right away if any problem should occur with any of your products.
8. Sterilize before first time use and before every subsequent use thereafter.
9. Indication for use: Used by a trained dental practitioner in conjunction with other attachments for the cleaning of teeth.

HOW TO USE:

1. Attachment of E-Type Motor. Align the pins on the Motor with the holes in the Air tubing turn the collar on the tubing clockwise until snug
2. To remove the Motor, Turn the collar on the tubing counter clockwise until loose.

CYCLES:

The Motor can be subject to 250 cycles = (Performed Sterilization Process) per EN ISO 14457 requirement 5.10 before showing signs of deterioration.

FEATURES & SPECIFICATIONS:

1. The Motor conforms to ISO 14457 specifications.
2. Connection with Attachments: meets ISO-9168 E-Type Connections.

CLEANING: Preparation and Cleaning before Sterilization:

**Clean and lubricate the handpiece as soon as possible after treatment. Failure to properly maintain the handpiece may cause infection, overheating leading to burn injuries.*

▲ DO NOT IMMERSE INSTRUMENTS IN DISINFECTANT SOLUTION OR CLEAN IN ULTRASONIC DEVICES. MANUAL CLEANING ONLY

1. PREPARATION:

- a) Always wear protective gloves, a mask, and protective goggles for safety purposes and to minimize the risk of infection.
- b) Remove the bur. (if applicable)
- c) Wipe the exterior of the handpiece/accessory using a cleaning wipe.
- d) Operate the handpiece at the chair-side for at least 20 seconds to purge fluids in the handpiece. (not applicable for air motors)
- e) Remove the handpiece, carry the handpiece to the decontamination area and proceed to step 2.

2. CLEANING AND DRYING:

**Wash the handpiece and irrigation nozzle under running water following the procedures below. (Water condition: <38°C (100°F), >3.5L/min; water should be same quality as drinking water)*

- a) Clean the external surfaces of the handpiece/accessory with a soft bristled toothbrush for 15 seconds or more.
- b) Clean around the bur insertion hole for 15 seconds or more, using an interdental brush (wire thickness 0.7mm) that conforms to size 4 in ISO 16409. (if applicable)
- c) Under appropriate lighting (natural office lighting), inspect the handpiece/accessory for blood or other debris. If any visible debris remain, repeat the process until it is visually clean.
- d) Wipe the exterior of the handpiece/accessory clean using a cloth moistened with disinfectant or ethanol with a concentration of 60-90% (ie. Isopropyl)

LUBRICATION:

Recommended to use approved autoclavable handpiece lubricant.

1. Remove Motor from airline.
2. Using spray lubrication can or lubrication station, apply lubricant into Motor Handpiece.
3. Reattach Motor to airline if not using a lubrication station.
4. Run entire tool at normal speed for 10 seconds.
5. Wipe off excess oil from entire tool surface.

**AUTOCCLAVING/STERILIZATION**

IMPORTANT NOTE: This handpiece or motor do not come sterilized from the manufacturer. It is also recommended that the handpiece or motor be sterilized after each use. Approved autoclaving or sterilization packaging is recommended to preserve item in sterile condition. Please note that after sterilization, handpiece should be used at your earliest convenience to make certain that sterilization is maintained.

1. Clean before sterilization (see CLEANING above)
2. Lubricate before sterilization (see LUBRICATION Above).
3. Set parameters for Moist Heat Sterilizer. PreVac cycle 134 degrees Celsius for 8 minutes and 20 minute dry time

WARRANTY INFORMATION:

essence Slow Speed Motors are warranted to be free of manufacturing defects for a period of 12 months* from date of purchase for the standard Motor.

**Warranty is conditional on proper use and maintenance of the Motor. Warranty is limited to repair or replacement at manufacturers' discretion. Under no circumstances will the manufacturer be held responsible for incremental or consequential damages. Please contact for warranty repair.*

Disclaimer:

CURION recommends that the sterilization process outlined in the Directions for use (DFU) described in this document be followed precisely. CURION is not responsible for any products supplied in which the end user does not follow the DFU.

Bisco Dental Products (Canada) Inc. Db a Curion Dental Products

2571 Smith Street Richmond, BC., V6X 2J1 Canada

See Curion Website for other products and product information
curion.ca or Phone 1.800.667.8811

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