

Introduction:

Single-use instruments (scissors and forceps) are disposable instruments intended for a single use during surgical procedures. They are designed to minimize the risk of cross-infection between patients and healthcare professionals.

Indications:

Single-use instruments (scissors and forceps) are used in surgical procedures to cut, dissect, or hold tissues. They are intended for use in a wide range of surgical specialties.

Contraindications:

- Single-use instruments are intended for single-use only and should not be re-sterilized or reused.
- Do not use if packaging is damaged or compromised.
- Do not use after the expiration date.

Precautions:

- Single-use instruments (scissors and forceps) are sterile and must be handled with care to maintain their sterility.
- Handle the instruments gently to avoid damaging them.
- Do not use any instruments that are visibly damaged.
- Use instruments that are appropriate for the intended surgical procedure.
- Dispose of the instruments in an appropriate sharps container after use.

Use:

- The instructions for use provided are intended to guide healthcare professionals in the safe and effective use of single-use instruments, specifically scissors and forceps. These instruments are intended for use in medical and surgical procedures where precision and accuracy are required.
- Before using the instruments, it is important to carefully inspect the packaging to ensure that it is intact and not damaged. If the packaging is damaged, the instrument should not be used, and an alternative should be selected.
- To use the instruments, open the package carefully and remove the instrument. It is important to handle the instrument carefully to avoid damaging the tips or other delicate parts. The instruments should be used only for their intended purpose and according to the manufacturer's instructions.
- After use, the instruments should be properly disposed of in a designated sharps container to avoid injury to healthcare professionals and others who may come into contact with them. The instruments are designed for single use only and should not be reused or re-sterilized.
- It is important to note that the use of damaged or expired instruments can result in injury to the patient or healthcare professional. Careful attention should be paid to the expiration date of the instruments, and they should be used before this date. The instruments should also be used in accordance with the instructions provided, and any deviations from these instructions should be avoided.

Storage:

- Store the instruments in a cool, dry place, away from direct sunlight.
- Store the instruments in their original packaging until ready for use.
- Do not store the instruments in areas with excessive heat, humidity, or moisture.
- Do not expose the instruments to any chemicals or disinfectants.

WARNINGS:

- Do not use the instruments if the packaging or sterile barrier is damaged or compromised in any way.
- Never reuse or re-sterilize single-use instruments, as this may cause contamination and infection.
- Always handle the instruments with care to avoid injury, as the blades are extremely sharp.
- Do not use the instruments if they appear defective or damaged in any way.
- Only trained and qualified medical personnel should handle and use these instruments.
- Keep the instruments away from direct sunlight, moisture, and extreme temperatures to avoid compromising their sterility.
- Always dispose of used instruments as medical waste in an appropriate sharps container.

DO NOT REUSE - SINGLE PATIENT USE - DO NOT RESTERILIZE