## Agar Sangre / Agar MacConkey / Agar GC

### USO

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[ventas@frogglab.com.mx](mailto:ventas@frogglab.com.mx)

Tel: 5518010660

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1. Leavit, J.M. Naidorf and P. Shugaefsky. 1955. The undetected anaerobe in endodontic, a sensitive medium for detection of both aerobes and anaerobes. The NY.J. Dentist. 25:377-382
2. Curry, A. S., G.G., Joyce, and G.N. Mc Ewwn, Jr. 1993. CTFA Microbiology guidelines. The cosmetic, Toiletry and Fragrance . Association. Inc. Washington, DC.
3. Brown, J.H. 1919. The use of blood agar for the study of streptococci. NY Monograph No. 9 The Rockefeller Institute for Medical Research.
4. The United States Pharmacopeia (USP XXIII) and The National Formulary (NF18). 1995 Sterility test, p. 1686-1690. United States Pharmacopeia Convention Inc., Rockville, MD.
5. Vera and Power. 1980 In Lennette, Balows, Hausler and Truand (ed.) Manual of clinical microbiology, 3rd ed. ASM, Washington, D.C.
6. MacConkey, A. 1905. Lactose-fermenting bacteria in feces. J. Hyg. 5:333-379.
7. United States Pharmacopeial Convention, Inc. 1995. The United States pharmacopeia, 23 ed. The United States Pharmacopeial Convention, Rockville, M. D.
8. Food and Drug Administration. 1995 Bacteriological analytical manual, 8th ed. AOAC International. Gaithersburg, MD.
9. Gray, L.D. 1995 Escherichia, Salmonella, Shigella and Yersinia, p450-456. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
10. MacConkey J. H. 5:33. 1905. Joseph Md. State. Dept. Health. Procedures, 1960. European Pharmacopoeia. 4th Ed. 2002.
11. Song Y., Roumagna P., Weill F-X., Wain J., Dolecek C., Mazzoni C.J., Holt K.E., y Achtman M. 2010 A multiplex single nucleotide polymorphism typing assay for detecting mutations that result in decreased fluoroquinolone susceptibility in Salmonella enteric serovars Typhi and Paratyphi A. Journal of Antimicrobial Chemotherapy. Vol. 65 1631-1641.
12. Jarvik T., Smillie C., Groisman E. A. y Ochman H. 2010 Short-Term Signatures of Evolutionary Change in the Salmonella enterica Serovar Typhimurium 14028 Genome. Journal of Bacteriology. Vol. 192 No. 2 560-567.
13. Porwollik S. 2011 Salmonella: From Genome to Function. Ed. Ilustrada Horizon Scientific Press. 300 Págs.
14. Farmacopea de los Estados Unidos Mexicanos: Suplemento para Dispositivos Médicos. 3a. Ed. -- México: Secretaría de Salud, Comisión Permanente de la Farmacopea de los Estados Unidos Mexicanos, 2014
15. Johnston, J. 1945. Comparison of gonococcus cultures read at 24 and 48 hours. J. Venerea. Dis. Inform. 26:239
16. Lankford, C.E., V. Scott, M.F. Cox, and W.R. Cooke. 1943. Some aspects of nutritional variation of gonococcus. J. Bacteriol. 45:321.
17. Thayer, J.D., and J.E. Martin, Jr. 1966. Improved medium selective for cultivation of N. Gonorrhoeae and N. Meningitidis. Public Health Rep. 81:559.
18. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1 American Society for Microbiology, Washington, D. C.
19. Spray, 1930. J. Lab. Clin. Med. 16:166.
20. Fortuna, 2008 Protocolo de Atención del paciente grave: Normas, procedimientos y guías de diagnóstico y tratamientos. Ed. Médica Panamericana. Págs. 482.

**WWW.FROGGLAB.COM.MX**

**ventas@frogglab.com.mx**

**Tel: 5518010660**