



Simultaneous Multiplex Real Time PCR (SIMUL-qPCR)

## Listeria species & monocytogenes Assay

Assay Kit Insert for AOAC PTM Workflow



### INTENDED USE AND APPLICATION

Food processors and associated laboratories can use the SIMUL-qPCR System as a quick and reliable method for detecting *Listeria* species and monocytogenes in food and environmental samples. All SIMUL-qPCR System assays are designed to have the same instrument run time allowing simultaneous identification of all SIMUL-qPCR System assays. In addition, each assay utilizes the power of multiplexing several targets during the same run. The SIMUL-qPCR *Listeria* species and monocytogenes Assay incorporates a multiplex approach to identifying both *Listeria* species and monocytogenes.

The Simultaneous Multiplex Real Time PCR (SIMUL-qPCR) *Listeria* species and monocytogenes Assay has been validated by the AOAC™ Research Institute under the Performance Test Methods™ program for frankfurters, ready to eat sliced turkey, cooked eggs, fresh raw soft cheese, chicken salad, ice cream, pasteurized milk, and frozen/cooked shrimp as well as stainless steel, ceramic tile, rubber, plastic, and sealed concrete environmental samples. USDA FSIS and FDA BAM methods were used for method comparison testing. The SIMUL-qPCR *Listeria* species and monocytogenes Assay was found to be equivalent to the reference methods. The limit of detection for this SIMUL-qPCR assay is 1,000 CFU/mL of enriched sample.

### PRODUCT INFORMATION

SKU	DESCRIPTION	UOM / QUANTITY
SMRT-LSLM-096	SIMUL-qPCR <i>Listeria</i> species and monocytogenes Assay Kit	1 KIT   96 Tests/Kit

### KIT COMPONENTS

SKU	DESCRIPTION	QUANTITY
KC-SMRT-LYB-25	Lysis Buffer	96 reactions   2 bottles   1 resealable pouch
KC-SMRT-LSLM-8S	PCR Assay Tubes	96 tests   12 strips of 8 tubes/caps   1 resealable pouch



p 250 450 7995 | appliedfooddiagnostics.com

387 Hazle Street, Nuremberg, PA 18241

### PRINCIPLE

This protocol is a multifaceted approach to the detection of *Listeria* species in a variety of food products and environmental samples. A specifically formulated media is utilized for enriching samples followed by cultural (detection plate) and rapid (qPCR) detection procedures. *Listeria* Recovery and Enrichment Broth (LREB) combines nutritional components with additional ingredients that are necessary to selectively improve recovery and growth of *Listeria*. The selective agents present in LREB have been optimized to efficiently inhibit competing normal bacterial flora without affecting the growth of *Listeria* species. LREB is formulated for buffering capacity to ensure growth in a variety of matrices.

During PCR amplification, forward and reverse primers hybridize to unique sequences of *Listeria* species and monocytogenes genomic DNA. A fluorogenic probe is included in the same reaction mixture which consists of a DNA probe labeled with a 5'-dye and a 3'-quencher. During PCR amplification, the probe is cleaved and the reporter dye and quencher are separated. The resulting increase in fluorescence can be detected on the real-time PCR instrument. Two unique and specific primers and probe mixtures are present in this assay. The target for *Listeria* species is found in the FAM™ dye channel. The target for *Listeria* monocytogenes is found in the CAL Fluor® Orange 560 dye channel. An Internal Amplification Control (IAC) is included which is found in the CAL Fluor® Red 610 dye channel.

The combination of the multiple detection approaches ensures cost-effectiveness and efficiency.

### ADDITIONAL MATERIALS REQUIRED

Other necessary materials not provided include:

### ALL SAMPLES

- AFD *Listeria* Recovery and Enrichment Broth (LREB)
- MyGo Pro real-time PCR instrument and installed MyGo Pro software v3.4
- Autoclave
- Distilled / deionized water
- Sterile stomacher / blender bags or equivalent with or without filter
- Stomacher / blender or equivalent
- Incubator: at 30 ± 1°C
- Incubator: at 35 ± 1°C
- Heating blocks with inserts
- Vortexer
- Calibrated thermometer
- Capping / de-capping tools (optional)
- Mini-centrifuge (optional)
- Adjustable mechanical pipettes
- Multi-channel pipette
- 1.5 ml microcentrifuge tubes or equivalent PCR-grade plastics for cell lysis
- Microcentrifuge tube racks
- Powder-free gloves
- Routine laboratory equipment
- Positive Control
- Negative Control

### ENVIRONMENTAL SAMPLES

- AFD Quick Release Sampler Hydrated with Neutralizing and Recovery Buffer (SKU#: SC-SQRS-NRB-100)
- AFD Quick Release Sampler Hydrated with Neutralizing Buffer (SKU#: SC-SQRS-NB-100)
- AFD Quick Release Sampler Hydrated with Lethen Broth (SKU#: SC-SQRS-LB-100)
- AFD MEMP Swab Kit (SKU#: MEMP-SWB-032)
- Other commercially available swabs

### PRECAUTIONS

This product is for In Vitro diagnostic use only. Do not ingest, inhale, or allow to come into contact with skin. Observe approved biohazard precautions and aseptic techniques. Biosafety level 2 procedures should be exercised (BMBL, <http://www.cdc.gov/biosafety/publications/bmbl5/bmbl.pdf> or current site). Extreme care should be taken in handling test samples and enrichment broths. All enrichment broths, plates and other items may contain various pathogens whether or not they contain *Listeria* species. This kit is to be used only by adequately trained and qualified laboratory personnel in a laboratory setting. All laboratory specimens should be considered infectious and handled accordingly.

### PROCEDURE

#### Media Preparation - Autoclave Method

1. Use a clean bottle for each liter of medium preparation.
2. Shake container of LREB before each use.
3. Measure 37.0 g of powder into the bottle and add 1 L of distilled water.
4. Constantly stir and heat solution until powder is dissolved. The acceptable pH is 7.2 ± 0.2.
5. Sterilize the bottle(s) of prepared medium by autoclaving at 121°C for 15 min.
6. Cool bottle to room temperature. Media is stable at room temperature or can be stored at 2–8°C for up to 45 days. Keep away from light.

#### Media Preparation - Non-Autoclave Method

1. Use a sterile and clean bottle for each liter of medium preparation.
2. Shake container of LREB before each use.
3. Measure 37.0 g of powder into the bottle and add to 1 L of sterile distilled or deionized water.
4. Constantly stir and heat solution until powder is dissolved. The acceptable pH is 7.2 ± 0.2.
5. Cool the prepared medium to the appropriate temperature (35 ± 1°C) and use immediately.

#### Media Preparation – Enrichment Ready Media

1. Prepare a sample bag with pre-warmed (35 ± 1°C) laboratory grade ASTM D1193 Type 2 water with the volume respective of the pouch size.
2. Tap pouch to compact media away from the tear mark. Pinch and open the perforated tear mark and empty the contents into the sample bag.
3. Shake, stomach or hand mix until powder is dissolved.
4. Use within three to four hours while maintaining pre-warmed conditions (35 ± 1°C).

### ENVIRONMENTAL SURFACE SAMPLE PREPARATION Sponge

1. When ready to test, pre-warm the prepared LREB to 35 ± 1°C.
2. Add the pre-warmed LREB to each sponge sample. Refer to Enrichment of Samples table for acceptable media volumes.
3. Homogenize the sample for 30 seconds in a stomacher / blender or equivalent. Hand mixing is an acceptable alternative for stomaching. To hand mix, massage each sponge that is in the sealed bag for approximately one minute.
4. Incubate the sample. Refer to Enrichment of Samples table for enrichment conditions.

### Swab

1. When ready to test, pre-warm the prepared LREB to 35 ± 1°C.
2. Add 9 mL of the pre-warmed LREB to each swab sample in its vial.
3. Incubate the sample. Refer to Enrichment of Samples table for enrichment conditions.

### MEMP Assay Swab

1. Refer to Molecular Environmental Monitoring Program (MEMP) Assay for *Listeria* Detection Kit Insert.
2. Sample and process the lysate as directed in the AFD MEMP Assay.
3. Proceed to AFD MEMP assay.
4. Any samples resulting positive on the AFD MEMP Assay must be enriched and run on the SIMUL-qPCR *Listeria* Assay for confirmation. Add the retained swab and sample solution to 9 mL of LREB and incubate for 30-36 hours at 30 ± 1°C.

### FOOD SAMPLE PREPARATION

1. Aseptically sample the product and place in a sterile bag.
2. When ready to test, pre-warm the prepared LREB to 35 ± 1°C.
3. Add the pre-warmed LREB to the sample. Refer to Enrichment of Samples table for acceptable media volumes.
4. Homogenize the sample for 30 seconds in a stomacher / blender or equivalent. Hand mixing is an acceptable alternative for stomaching. To hand mix, massage each sample that is in the sealed bag for approximately one minute.
5. Incubate the sample. Refer to Enrichment of Samples table for enrichment conditions.

LREB ENRICHMENT OF SAMPLES		
Enrichment Incubation: 30 ± 1°C for 30-36 hours		
MATRIX	SAMPLE SIZE / ANALYSIS UNIT	MEDIA VOLUME
Environmental Swab*	/Swab	9 mL ± 1 mL
AFD MEMP Reserve Swab*	/Swab	9 mL ± 1 mL
Environmental Sponge*	/Sponge	90 mL ± 10 mL
Frankfurter	/125 g	1 L ± 50 mL
RTE Sliced Turkey	/125 g	1 L ± 50 mL
Soft Fresh Cheese (Raw)	/25 g	225 mL ± 15 mL
Chicken Salad	/25 g	225 mL ± 15 mL
Ice Cream	/25 g	225 mL ± 15 mL
Cooked Eggs	/25 g	225 mL ± 15 mL
Pasteurized Milk	/25 g	225 mL ± 15 mL
Frozen/Cooked Shrimp	/25 g	225 mL ± 15 mL
*Validated surfaces are stainless steel, plastic, rubber, ceramic tile, and sealed concrete		

## ASSAY PREPARATION

- At the end of the enrichment phase, proceed to the cell lysis protocol for single samples.
- The qPCR set up and data entry should be completed prior to transferring samples.
- Prepare equipment:
  - Turn on the heating blocks to 95 ± 3°C as measured by a calibrated thermometer.
  - Power on the qPCR instrument and create run file from SIMUL-qPCR template. The SIMUL-qPCR template contains the required cycle. To test for *Listeria* species, select the target for *Listeria* species. To test for *Listeria monocytogenes*, select the target for *Listeria monocytogenes*. To test for both *Listeria* species and *monocytogenes*, select both targets.
- After removing the aliquot required for lysis, return the enrichment(s) to the incubator.

## SET UP THE qPCR INSTRUMENT

- Refer to the PCR Set Up Guide and the MyGo Pro PCR Software Manual for detailed instructions.
- Select the AFD Template files to begin run configuration. The AFD template file contains all of the PCR machine settings required to perform the run. Do not change any settings under the "Experiment", "Run Profile", or "Data" tabs.
- Under the "Samples" tab, populate the sample fields according to the well placement / position.
- Include the kit lot number in the "Notes" field.
- Add the targets to the sample(s).
- After cell lysis and loading, click "Start Run".

## CELL LYSIS PROTOCOL OF ENRICHED SAMPLES

After incubation, follow the steps below for cell lysis:

- Label one 1.5 mL microcentrifuge tube or equivalent PCR-grade plastic tube per sample.
- Aseptically pipette 400 µL of lysis buffer into each tube. Return lysis buffer to storage (2-8°C).
- From an enrichment broth, pipette 5 µL into the prepared tube. Cap the tube(s) and vortex.
- Heat the closed tubes for 10 minutes at 95 ± 3°C in the heat block.
- Remove the closed tubes from the heat block.
- Allow tubes to cool for 5 minutes at room temperature.
- Proceed directly to SIMUL-qPCR *Listeria* Assay or hold lysate in a refrigerator (2-8°C) up to 48 hours before proceeding to SIMUL-qPCR *Listeria* Assay.

## SET UP THE SIMUL-qPCR ASSAY

- Arrange strips of PCR tubes according to your run file.
- Using caution, remove the caps from the strip of tubes.
- Pipette 20 µL of lysate into the sample wells of the PCR test strip, ensuring the pellet is hydrated. PCR pellets must be hydrated and resealed within 10 minutes after removing the caps from the PCR tubes.
- Place the caps onto each tube and press down to seal each lid.
- Make sure each lid is tightly secured before running on the PCR machine.
- If air bubbles are present, carefully flick reaction tubes until no air bubbles remain.
- Briefly spin down the reaction tubes in a mini-centrifuge.
- Load the qPCR instrument and start the run.

## REVIEW RESULTS OF SIMUL-qPCR ASSAY

Once the SIMUL-qPCR Assay run is complete, data is analyzed automatically by the software. The software analyzes any DNA amplification data and will display a Cq value for any sample that amplifies. Amplification in the FAM™ channel indicates *Listeria* species. Amplification in the CAL Fluor® Orange 560 channel indicates *Listeria monocytogenes*. Only a Cq value that has a typical sigmoidal curve or the beginning of the curve is considered positive for the target. When a Cq value is not obtained, the result is negative for the target provided a Cq value is present in the CAL Fluor® Red 610 channel for the IAC.

All samples that are positive should be culturally confirmed by either the FDA BAM <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm071418.htm> or USDA FSIS MLG <http://www.fsis.usda.gov/wps/wcm/connect/1710bee8-76b9-4e6c-92fc-fdc290dbfa92/MLG-8.pdf?MOD=AJPERES> method.

## PRODUCT STORAGE AND EXPIRATION

Store the sealed kit at 2 - 8°C. Once opened, protect kit components from moisture and light by keeping container(s) tightly closed after each use. Re-seal qPCR tubes in re-sealable foil pouch. The expiry date is indicated on the package.

## DISPOSAL

Dispose of all materials used and the enrichment medium by autoclaving or according to approved practices.

Ensure that all biohazard waste is disposed of according to local, municipal, provincial, state and/or federal regulations.

## TECHNICAL INFORMATION

If you have any questions or experience issues with this kit, please contact our support staff via email ([support@appliedfooddiagnostics.com](mailto:support@appliedfooddiagnostics.com)). For more information about Applied Food Diagnostics, Inc., please visit us at our website ([www.appliedfooddiagnostics.com](http://www.appliedfooddiagnostics.com)).

## QUALITY CONTROL

All products manufactured by Applied Food Diagnostics, Inc. are incorporated into a quality assurance program from the time the raw materials arrive in the factory through to marketing the end product. Each batch of end product undergoes quality control, and is only marketed if it complies with acceptance criteria. Documentation concerning production and verification of each batch is archived. A Certificate of Analysis of this quality control and Safety Data Sheets are available on the web at [www.appliedfooddiagnostics.com](http://www.appliedfooddiagnostics.com).

## TERMS AND CONDITIONS

Applied Food Diagnostics, Inc. makes no representations and warranties concerning its products other than those stated herein. All product(s) delivered hereunder by Applied Food Diagnostics, Inc., its affiliates, or any other person on its behalf shall, at the time of delivery, be manufactured to meet Applied Food Diagnostics, Inc.'s specifications and all applicable laws. All other terms, conditions and warranties, including any warranty of merchantability, quality, fitness or suitability for a particular or intended purpose, implied by common law or statute, (implied warranties) are expressly excluded.

## DISCLAIMER

This kit and its performance characteristics were developed by Applied Food Diagnostics, Inc., for laboratory use. Any deviations from this protocol are not authorized by Applied Food Diagnostics, Inc.

## LIMITED USE LABEL LICENSE

This product is covered by at least one or more claims of US patents applications, which are licensed to Applied Food Diagnostics, Inc. This product is sold strictly for the use of the buyer, and the buyer is not authorized to transfer this product, or any materials using this product, to any third party.

\*All Trademarks are property of their respective owners.



p 570 450 7995 | [appliedfooddiagnostics.com](http://appliedfooddiagnostics.com)  
387 Hazle Street, Nuremberg, PA 18241



p 570 450 7995 | [appliedfooddiagnostics.com](http://appliedfooddiagnostics.com)  
387 Hazle Street, Nuremberg, PA 18241