



**Standard Operating Procedure**

**R-CARD® KwikCount® EC**

**Rapid Test Method for *Escherichia coli* (*E. coli*)  
In Ambient Waters**

Roth Bioscience, LLC  
1303 Eisenhower Drive S.  
Goshen, IN 46526  
(574) 533-3351  
[hello@rothbioscience.com](mailto:hello@rothbioscience.com)

## 1. Scope and Application

- 1.1. This method describes a procedure with the R-CARD® KwikCount® EC (Roth Bioscience, LLC) for a rapid detection and enumeration of *Escherichia coli* (*E. coli*) within 7 to 9 hrs. Because these bacteria are natural inhabitants of the intestinal tract of warm-blooded animals, their presence in water samples are an indication of fecal pollution and the possible presence of other enteric pathogens. This test for *E. coli* can be applied to fresh water, or other waters matrices. The R-CARD® KwikCount® EC medium (Roth Bioscience LLC, Goshen, Indiana) can be used to detect and enumerate *E. coli* within 7 to 24 hrs.

## 2. Summary of Method

- 2.1. A liquid sample is pipetted on the center of the card, and covered by the top film. The liquid sample will spread laterally automatically within 1 min. The card is then incubated at 41±0.5°C for 7 -9 hr. The R-CARD® are then viewed in a dark place using a long wave UV light (366nm preferred) and any fluorescent colony is counted as *E.coli* .
- 2.2. R-CARD® KwikCount® EC contains nutrients that promote the growth of the target organisms and buffers to maintain appropriate pH. It also contains a fluorogenic/chromogenic enzyme substrate mixture that is used to cause *E.coli* to fluoresce and is detectable as early as 7 hrs incubation time (7-9 hrs incubation is the peak fluorescent intensity.) The fluorescent light will gradually decrease until about 24 hrs when it will be totally dissipated. Beginning about 14-16 hrs incubation, the *E.coli* colonies will not only fluoresce in UV light, but will also appear as blue/purple colonies in ambient light. This phenomenon of a coloring compound having two unique properties simultaneously has been trademarked as a DUOGEN®.
- 2.3. The R-CARD® KwikCount® EC method therefore incorporates and includes the advantages of utilizing a Duogen®. When the blue fluorescence has dissipated, the *E.coli* colonies will continue to be present as blue/purple colonies that may continue to grow and intensify in color. The blue/purple colonies should not be counted after 48 hrs incubation.

## 3. Method Definition

- 3.1. In the R-CARD® KwikCount® EC method, ***E. coli*** are those bacteria which produce fluorescent colonies on the medium between 7-16 hrs. incubation when using a long wave UV light source (366nm preferred).
- 3.2. R-CARD® KwikCount® EC is ready-to-use for detecting *E. coli* in liquid samples.

## 4. Interferences

- 4.1. In the case that excessive crowding is observed on the R-CARD®, samples should be recollected and rerun at a higher dilution. Very small colored colonies should not be counted unless they are isolated as individual colonies with densities below 80 CFU per CARD.

## 5. Safety

- 5.1. Analyst/technician must know and observe the normal safety procedures required in a microbiology laboratory while preparing, using, and disposing of cultures, reagents, and materials and while operating sterilization equipment.
- 5.2. Mouth-pipetting is prohibited.

## 6. Suggested Equipment and Supplies

- 6.1. Sterile pipettes (1 to 25 mL)
- 6.2. Forceps: smooth, flat, sterilizable metal forceps.
- 6.3. Microscope: 10 to 15 X magnification binocular wide-field dissecting microscope.

- 6.4. Light box
- 6.5. Bunsen burner or alcohol lamp for sterilizing forceps..

## 7. Reagents and Standards

- 7.1. Sterile deionized or distilled water
- 7.2. R-CARD® KwikCount® EC

## 8. Quality Assurance/Quality Control

- 8.1. Quality control
  - 8.1.1. Each lot of the R-CARD® KwikCount® EC medium should be evaluated by the laboratory by preparing three plates of the medium (one to serve as an un-inoculated control, one to serve as a negative growth control, and one to serve as positive control).
  - 8.1.2. *E. coli* ATCC #11775 or 25922 is used as the positive control. *Klebsiella pneumoniae* ATCC 31488 is used as the negative control, and also *Pseudomonas aeruginosa* ATCC 10145 or 27853 may be used as a negative growth control microorganism.

## 9. Procedure

- 9.1. Prepare samples as usual and make a serial dilution if necessary.
- 9.2. Wear glove and lift the top portion (film) or use sterile forceps (see photos 1-2)
- 9.3. Select dilutions of the sample to produce 20-150 *E. coli* colonies on the cards.
- 9.4. Pipette 1 mL of the sample onto the center of the card (photo 3).
- 9.5. Cover with the film, and wait about 1-2 min to allow liquid to spread automatically. There is no need to use a spreader. (photo 4). (Some samples do not automatically spread in as large an area as may be wanted, but it is a simple matter to encourage spreading by gently applying pressure on the top after it is lowered on the inoculum.
- 9.6. Incubate at 41±0.5°C for 7-24 hrs (no more than 48 hrs).




Photo 1. Open the film




Photo 2. Lift the film




Photo 3. Pipette 1 mL sample






Photo 4. Cover with the top film

## 10. Data Analysis and Calculations

- 10.1. Count the number of colonies detected by blue fluorescence present on the card between 8-16 hr incubation and record as the number of E. coli / volume of sample for that test.



UV light 8 hr  
Blue fluorescence (Duogen® effect) E.coli



Ambient light at 20 hr  
Blue-purple (Chromogenic color) E.coli

## 11. Pollution Prevention and Waste Management

- 11.1. All biohazardous waste should be sterilized at 121°C for 30 min prior to disposal. Laboratory personnel should use pollution control techniques to minimize waste generation wherever possible.