



Standard Operating Procedure

R-CARD® KwikCount® EC

Rapid Test Method for *Escherichia coli* (*E. coli*) In Ambient Waters

Roth Biosciences, LLC
220 Bloomingdale Dr.
Bristol, IN. 46507
(574) 533-3351
hello@rothbioscience.com

1. Scope and Application

- 1.1. This method describes a procedure with the R-CARD® KwikCount® EC for a rapid detection and enumeration of *Escherichia coli* (*E. coli*) within 7 to 20 hrs and does not require a confirmation step. This test for *E. coli* can be applied to fresh water, or other waters matrices. It is also useful in the analysis of other waters, such as treated-wastewater effluent testing. Because these bacteria are natural inhabitants of the intestinal tract of warm-blooded animals, their presence in water samples are an indication of fecal pollution and the possible presence of other enteric pathogens.
- 1.2. The detection limit is one colony forming unit (CFU) per sample.

2. Summary of Method

- 2.1. A liquid sample is pipetted on the center of the card, and slowly covered by the top film. The liquid sample will spread laterally automatically within 1 min. The card is then incubated at $41\pm0.5^{\circ}\text{C}$ for 7 – 24 hrs. The cards are then viewed in a dark place using a long wave UV light (366nm preferred) and any fluorescent colony is counted as *E.coli*.
- 2.2. R-CARD® KwikCount® EC contains nutrients that promote the growth of the target organisms and buffers to maintain appropriate pH. It also contains a fluorogenic/chromogenic enzyme substrate mixture that is used to cause *E.coli* to fluoresce and is detectable as early as 7 hrs incubation time (7-9 hrs incubation is the peak fluorescent intensity.) The fluorescent light will gradually decrease until about 24 hrs when it will be totally dissipated. Beginning about 14-16 hrs incubation, the *E.coli* colonies will not only fluoresce in UV light, but will also appear as blue/purple colonies in ambient light. This phenomenon of a coloring compound having two unique properties simultaneously has been trademarked as DUOGEN®.
- 2.3. The R-CARD® KwikCount® EC method therefore incorporates and includes the advantages of utilizing a Duogen®. When the blue fluorescence has dissipated, the *E.coli* colonies will continue to be present as blue/purple colonies that may continue in growth and intensify of color. The blue/purple colonies should not be counted after 24 hrs incubation.

3. Method Definition

- 3.1. In the R-CARD® KwikCount® EC method, *E. coli* are those bacteria which produce fluorescent colonies on the medium between 7-20 hrs incubation when using a long wave UV light source (366nm preferred).
- 3.2. R-CARD® E. Coli is ready-to-use for detecting *E. coli* in liquid samples.

4. Interferences

- 4.1. In the case that excessive crowding is observed on the R-CARD®, samples should be recollected and rerun at a higher dilution. Small colored colonies of this nature should not be counted unless they are isolated as individual colonies with densities below 80 CFU per CARD.

5. Safety

- 5.1. Analyst/technician must know and observe the normal safety procedures required in a microbiology laboratory while preparing, using, and disposing of cultures, reagents, and materials and while operating sterilization equipment.
- 5.2. Mouth-pipetting is prohibited.

6. Equipment and Supplies

- 6.1. Sterile pipettes (1 to 25 mL)
- 6.2. Forceps: smooth, flat, sterilizable metal forceps.
- 6.3. Microscope: A 10 to 15 X magnification binocular wide-field dissecting microscope.
- 6.4. Light box
- 6.5. Bunsen burner or alcohol lamp for sterilizing forceps if necessary.

7. Reagents and Standards

- 7.1. Sterile deionized or distilled water
- 7.2. R-CARD® KwikCount® EC

8. Quality Assurance/Quality Control

- 8.1. Quality control
 - 8.1.1. Each lot of the R-CARD® KwikCount® EC medium should be evaluated by the laboratory by preparing three plates of the medium (one to serve as an uninoculated control, one to serve as a negative growth control, and one to serve as positive control).
 - 8.1.2. *E. coli* ATCC #11775 or 25922 is used as the positive control. *Klebsiella pneumoniae* ATCC 31488 is used as the negative control, and also *Pseudomonas aeruginosa* ATCC 10145 or 27853 may be used as a negative growth control microorganism.

9. Procedure

- 9.1. Prepare samples as usual and make a serial dilution if necessary.
- 9.2. Wear glove and open the top portion (film) or use sterile forceps.
- 9.3. Select dilutions of the sample to produce 20-150 *E. coli* colonies on the cards.
- 9.4. Pipette 1 mL of the sample on the center of the card.
- 9.5. Cover with the film, and wait 1 min to allow liquid to spread automatically. There is no need to use a spreader. (Some samples do not automatically spread in as large an area as may be wanted, but it is a simple matter to encourage spreading by gently applying pressure on the top after it is lowered on the inoculum.
- 9.6. Incubate at 41±0.5°C for 7-24 hrs (no more 24 hrs).

10. Data Analysis and Calculations

- 10.1. Count the number of colonies detected by blue fluorescence present on the card between 7-20 hrs incubation and record as the number of *E. coli* / volume of sample for that test.

11. Pollution Prevention and Waste Management

- 11.1. All biohazardous waste should be sterilized at 121°C for 30 min prior to disposal. Laboratory personnel should use pollution control techniques to minimize waste generation wherever possible.