Analysis of Single Amino Acid Mutations in Intact Proteins

Robert Freeman¹, Dan DiFeo¹, Kayte Parlevliet², Hans-Jürgen Wirth² ¹SGE Incorporated, 2007 Kramer Lane, Austin Texas, USA, ² SGE Analytical Science, 7 Argent Place, Ringwood Victoria, Australia.

Abstract


With the increasing availability of high resolution mass spectrometers the analysis of intact proteins becomes a feasible alternative to the more commonly employed shot-gun techniques in proteomics. High resolution separations of proteins are hindered by the size of the analyte molecules and the resulting low diffusion rates especially inside the pore system of the stationary phase.

The stationary phase of the columns employed in these studies is a 3 µm C8 silica with 1000 Å pore size. By using a dedicated protein separation column (ProteCol™ C8 HQ1003) we demonstrate that it is possible to separate ribosomal proteins derived from a rat liver with excellent resolution and identify single amino acid mutants in a 33 kDa protein.

Background


Protein Separation and Diffusion

Diffusion is the time-limited process in the separation of macro-molecules such as proteins. Due to their size, the diffusion rate of proteins in free solution is significantly lower than the diffusion rate of small molecules (see Figure 1). In addition, the effect of restriction inside the pores of the stationary phase is far more pronounced (Figure 2 and 3). To minimize the effects of hindered pore diffusion, a 1000 Å 3 µm C8-silica column (ProteCol™ C8 HQ1003) was employed for this study.





human serum albumin).


Renkin Equation



ProteCol-C8 HQ1003







IgG Mw ~ 150kDa Transferrin Mw~80kDa


hSA Mw ~ 68kDa

Figure 3. Relative diffusion rates in free solution and inside pores as function of molecule size to pore ratio.

Ribosome Protein Identification

The rat liver 80S ribosome is made up of two subunits — the large subunit has ~49 proteins and the small subunit has ~33 proteins. Since ribosomal proteins are relatively small (Mw=6,000 to 40,000), they can be rapidly identified by accurate LCMS analysis of the intact proteins.

Dedicated Protein Columns – C8, 3 µm, 1000 Å



A glass lined HPLC column of conventional design was used with 0.5 µm pore sized PEEK™ frits within a stainless steel body to ensure an inert flow path for the sample (Figure 4). The C8-silica stationary phase with 3 µm particles with 1000 Å pore size is ideal for the separation of intact proteins.




Figure 5. Design of the ProteCol™ range of columns.

The capillary column is made of a packed 1/16" OD PEEKsil™ (polyether ethyl ketone coated fused silica) capillary, with 1/32" OD, 50 µm ID PEEKsil[™] connection capillaries permanently attached (Figure 5). The frits are made from a 100 µm thick woven mesh with 0.5 µm pore size. The column is protected by a stainless steel sleeve. Using PEEKsil™ rather than fused silica increases the robustness of the column while maintaining flexibility to easily attach the column to the LC system. The same stationary phase as above was used.

81 min 85 % B

88 min 5 % B

Methodology/Procedure

Analysis of Ribosomal Proteins

Step 1. Sample Preparation

80S ribosomal proteins were isolated from a rat liver microsomal preparation¹. One optical density unit at 260 nm of 80S ribosomal proteins was mixed with 2 volumes of 6M Guanidine HCl to denature the proteins. 1 % (v/v) formic acid was subsequently added to precipitate the nucleic acids. The mixture was centrifuged for 15 minutes at 13,000 rpm and the supernatant was collected into a sample vial ready for LCMS analysis.

Step 2. LCMS Analysis

Agilent 1100 LC with Agilent 6220 ESI-TOF LC/MS Mass detection ProteCol™ – C8 HQ1003 3 μm, 1000 Å 150 mm x 2 mm ID Column

0.25 µL/min Flow rate:

Aqueous 0.1 % (v/v) Formic Acid Acetonitrile 0.1 % (v/v) Formic Acid

80 min 45 % B 0 min 5 % B Gradient 82 min 85 % B 83 min 5 % B

Analysis of Single Amino Acid Proteins


Nine single-amino acid mutations of a bacterial 33 kDa protein toxin were analysed by RP-µLC-MS in order to confirm the correct structure of the mutant protein. Samples supplied at concentrations ranging from 0.3 to 1 mg/ml. 100 µL of samples were desalted with spin columns and resultant filtrate directly injected at a volume to equal the injection of 1 μ g i.e. 2 μ L injected of 0.5 mg/mL solution assuming no losses at filtration step.


Agilent 1200 Cap-LC with Bruker microTOF-Q MS ProteCol™ – G C8 HQ1003, 300 µm ID x 100 mm

Mobile phase B: 95 % acetonitrile, 0.1 % formic acid 33 min 70 % B 30 min 55 % B


40 min 5 % B

Data for the ribosomal proteins were acquired and reference mass corrected via a dual-spray electrospray ionization (ESI) source. Each scan or data point on the Totallon Chromatogram (TIC) is an average of 15,000 transients, producing a spectrum every second. Mass spectra were created by averaging the scans across each peak and background subtracted against the first 10 seconds of the TIC. The resulting base peak chromatogram shows very high peak capacity (Figure 6). 119 discrete protein masses were identified; 4 of which were identified as 80S ribosomal proteins. In some cases several different masses of the same protein were identified which correlated with known N- and/or C- terminal processing, highlighted in Table 1.






Based on the excellent resolution of the intact ribosomal proteins we anticipated the column would be suitable for the analysis of single amino acid mutations in recombinant proteins. The results for the single amino acid mutations intact proteins are shown in Table 2, and Figure 7.



The shift in retention time can be correlated to the change in hydrophobicity. A hydrophobicity scale published by Cowan and Whittaker was used because it is based on HPLC retention data². As shown in Figure 11, there is a good correlation between the change in hydrophobicity and the retention time. The open circles indicate mutated proteins where the measured mass differs significantly from the predicted mass.



Conclusion

The separation of proteins by reversed phase chromatography remains an uncommon technique due to the poor recovery of proteins from the commonly used 300 Å pore size stationary phases. However, with the increasing resolving power of the modern mass spectrometer and the inherent compatibility of the reversed phase solvent systems with MS analysis, we explored the potential of wide pore silicas for the accurate mass analysis of intact ribosomal proteins. We demonstrate excellent resolution of the rat liver 80S ribosome, where in some cases several different masses of the same protein were identified. Further, we were able to separate variations of a 33 kDa protein which differed by a single amino acid mutation.

The change of a single amino acid influences the hydrophobicity and this change correlates well with the shift in retention time of the intact protein.

In the separation of large molecules, especially in an adsorption mode such as reversed phase chromatography, it is important to have unrestricted access to the surfaces inside the pores even when an analyte molecule is already adsorbed at the pore entrance.

References

. Williamson et al. Eur. J. Biochem, 1997; 246: 786-793

2. R. Cowan and R.G. Whittaker. Peptide Research,1990; 3: 75-80

Acknowledgements

The preparation, chromatography and analysis of the ribosomal proteins were done by Dr Nicholas Williamson and Paul O'Donnell from Bio21 Mass Spectrometry and Peptide Synthesis, Bio21 Molecular Science and Biotechnology Institute, The University of Melbourne, Parkville, Victoria, Australia.

The single amino acid mutations in intact proteins analysis was done in collaboration with D. Steer, Monash Biomedical Proteomics Facility together with X. Yan, T. Bannam and J. Rood from the Department of Microbiology, Monash University, Victoria, Australia.

