

## Protocol IV: Spectrophotometric quantification of biofilm in colony resuspensions

This protocol describes how to quantify the biofilm markers curli and cellulose in colony re-suspensions. We have tested this procedure with *Salmonella Enteritidis* and *Salmonella Typhimurium* strains. For these strains we have not observed staining of intracellular or membrane components.

### Solutions and Reagents:

ECtracer™ is provided as 1000-fold concentrated solution. The following common reagents are required (not supplied):

- Agar plates
- Phosphate buffered saline (PBS), pH 7.4
- 96-well plate (round bottom)
- Spectrophotometer

### Assay Procedure:

- Grow bacterial colonies on an agar plates under biofilm forming conditions. Notice: no morphotyping is required for this procedure.
- Dilute ECtracer™ in PBS 1:1000
- Add 100 µl into each well of a 96-well plate.
- Pick bacterial colonies from the agar plate and resuspend thoroughly into each of the pre-filled wells.
- Place the plate in a spectrophotometer to quantify biofilm in the colony resuspensions.

### Spectrophotometer Settings:

- **ECtracer™480 (part of ECtracer™ Mix&Try Kit):** Excite at 430 nm and collect emission at 480 nm. Optional: Record an emission spectrum (450 - 700 nm) with 430 nm excitation.
- **ECtracer™520 (part of ECtracer™ Mix&Try Kit):** Excite at 470 nm and collect emission at 530 nm. Optional: Record an emission spectrum (490 - 700 nm) with 470 nm excitation.
- **ECtracer™630 (part of ECtracer™ Mix&Try Kit):** Excite at 510 nm and collect emission at 635 nm. Optional: Record an emission spectrum (530 - 800 nm) with 510 nm excitation.
- **ECtracer™680:** Excite at 540 nm and collect emission at 680 nm. Optional: Record an emission spectrum (560 - 800 nm) with 540 nm excitation.

---

#### Head Office

Ebba Biotech AB  
Medeon Science Park  
Per Albin Hanssons väg 41  
SE-205 12 Malmö  
Sweden

#### Stockholm Branch

Ebba Biotech AB  
Karolinska Institutet Science Park  
Fogdevreten 2a  
SE-171 65 Solna  
Sweden

#### Company Information

Org-nr: 559016-7093  
VAT-nr: SE 559016709301

#### Contact

Web: [ebbabiotech.com](http://ebbabiotech.com)  
Email: [info@ebbabiotech.com](mailto:info@ebbabiotech.com)  
Phone: +46 73 985 40 51