

Page | 1

### Method Z410 - Iron Fe

# **Specification**

Description: Test for determining the content of iron in fresh and marine water

Range: 0,05-10 mg/l
Resolution: 0,01 mg/l
Wavelength: 520 nm

Reagent set

Product Code Description List of components

8410 Set of reagents for method Z410, Iron Fe ✓ powder Reagent Fe

(reagents for approx. 30 tests) ✓ spatula

# Performing the measurement

Select the Z410 Iron Fe method (Methods → Select method → Z410 Iron Fe).
 How to select the method, see 8.1 Choosing method.

#### NOTE:

It is recommended to use the GUIDE system by pressing the context button **GUIDE** on the photometer. It will provide you with step-by step basic instruction how to perform measurement and a timer with beeper to count down reaction time. To enable this function press the button **GUIDE**.

2. Rinse the vial and the syringe three times with the tested water.


Take exactly 5 ml of the tested water with the syringe and pour into the vial.

#### NOTE:

Make sure no air bubbles are present in the syringe. Trapped air bubbles can affect accuracy of the measurement.



3. Insert the vial into the round vial holder and press the **ZERO** key. The display will show "-0.0-", which means the device is ready for measurement.




26 08	20	<b>(III</b>	] 12:35		
Fe	Z410 Iron Fe				
	tag 1				
Measuring					
ZERO	MEAS	GUIDE			

26 08	20	-[111]	12:35		
Fe	Z410 Iron Fe				
' -	tag 1				
-0.0- mg/l					
ZERO	MEAS	GUIDE			




Page | 2

- 4. Add 1 portion of powder Reagent Fe with the spatula 5. Before taking a measurement, wait exactly into vial, shake until the powder has dissolved.
  - 5 minutes.





6. Insert the vial into the round vial holder and press the MEAS key to take a measurement. The result – the concentration of iron - is displayed in mg/l (ppm).



26 08	20	12:40		
Го.	Z410 Ir	on Fe		
Fe	tag 1			
Measuring				
ZERO	MEAS	GUIDE		

		2/2		
26 08 20		•	12:40	
Fe	Z410 Ir	on Fe		
16	tag 1			
0.15 mg/l				
ZERO	MEAS	GUIDE	REC	

### Potential interferences

the presence of:

- above 5 ppm copper (Cu) nickel (Ni) - above 5 ppm zinc (Zn) - above 5 ppm cadmium (Cd) - above 5 ppm

may interfere with the measurement