
R_xPERFECT

Perfect Dose. Perfect Health.

4Pocket User Guide



BENEFITS

Welcome to your brand new 4Pocket prescription and vitamin reminder. 4Pocket couples an extremely simple to use pill timer with a built-in pocket sized pill box. The solution was designed for the minimalist who has been prescribed a simple dosing schedule yet needs to take their medications with them. The notification system comes in the form of a strong vibrating and/or loud auditory alarm. Perfect for the pocket or purse at an extremely thin 19mm. It's really thin and a great way to Live Longer by always taking your prescriptions, supplements and vitamins perfectly every day.

FEATURES

- Up to 4 daily vibrating and/or auditory alarms
- Snooze function
- Built-in 4 compartment pill box
- Capacity: 32 aspirin sized tablets
- Only 19mm thin
- Beautiful back-lit display

Next Page

Set-Up

NOTE: If this is the initial setup, please remove TWO (both) insulator strips located under the pill tray. One insulator is for the clock and the other insulator for the alarm.

CLOCK - Setting the Clock:

1. Slide the left side switch to T-SET.
2. Press the HOUR button to set the correct Hour. Please observe the PM indicator.
3. Press or hold the MINUTE button to set the correct Minutes.
4. If necessary press the 12H/24H button after setting the time to switch between 12 hour (AM/PM) format and 24 hour format.
5. Slide the switch back to TIME.

ALARM TIME(S) - Setting the Alarm Time(s):

For each alarm slide the left switch (Alarm#) to 1, 2, 3 or 4.

1. Slide the left switch to the alarm you want to set.
2. Press ALARM ON/OFF button. A “pill” symbol appear with the alarm #.
3. Press or hold the HOUR button to set the correct Hour. Please observe the PM indicator.
4. Press or hold the MINUTE button to set the correct Minutes.
5. Slide the switch back to TIME. Repeat for each alarm you want to set.

ALARM MODE - Setting the Alarm Mode:

Select how you want to be notified when the alarm goes off. The setting applies to all alarms.

1. Slide the right switch to set desired alarm mode. VIB (Vibrating Alarm), BEEP (Beeping Alarm) or VIB+BEEP (both).

Reminder Time – How to use

ALARM Beep/Vibrate - What to do when the ALARM goes OFF:

1. ALARM will beep/vibrate (or both) for 30 seconds. The alarm will repeat twice - 4 minutes and after 8 minutes after initial alarm. While the alarm is beeping/vibrating you will see:
 - a. A flashing “Pill Symbol”
 - b. The Alarm # (1, 2, 3 or 4) corresponding with the proper pill box on the back
 - c. The scheduled alarm time
2. TAKE your MEDICATION(S).
3. After you have taken your medication(s), press ROUND button (ALARM ON/OFF).
4. The alarm will go off again tomorrow at the same time.

Other

BATTERY – How to Replace the Batteries:

NOTE: If this is the initial setup, please remove TWO (both) insulator strips located under the pill tray. One insulator is for the clock and the other insulator for the alarm.

4Pocket uses four (4) replaceable batteries:

- One (1) AAA battery is used for the alarm (sound and vibration).
- Three (3) AG13 batteries (button cells) are used for the clock.
- Replace the AAA battery vibration/sound is getting weak.
- Replace the AG13 batteries when the display fades.
- Normal battery life is 12 months.

To open the battery compartment:

1. Slide the switch labeled OPEN in the direction of the arrow, and SLIGHTLY OPEN the top casing. *Be careful, if you use too much force to open you may risk breaking the hinges.*
2. Lift the back of the casing to separate it from the lower casing. The batteries are now visible.
3. Insert one (1) new AAA battery, take note of the (+/-) ends of the battery.
4. Insert three (3) button cell style batteries (type: GP76, LR44, SR44, L1154 or equivalent), with the positive (+) side up.
5. Close the casing.

RESET – How to reset:

1. Take a paper clip (not a sharp instrument) and insert it in the small hole above the word RESET. This located on back of the pill box, above the word “Medication” between pill box 2 and pill box 4. DO NOT use the speaker hole on the right side.
2. Press down until 12:00 is displayed.
3. Set the Time and Alarm(s).

BACKLIGHT – Turning on the backlight:

To turn on the backlight, press any button

Disclaimer: This Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. RxPerfect, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. RxPerfect, LLC's liability shall in no circumstances exceed the amounts paid by the user for the timer.

END