
R_xPERFECT

Perfect Dose. Perfect Health.

SlideBox User Guide



BENEFITS

Welcome to your brand new SlideBox prescription and vitamin reminder. SlideBox has been designed for those on a moderate to high complex dosing schedule but is flexible enough to be used with a wide range of regimens. The up to 37 daily auditory alarms are set either on the hour or half-hour with the attached 7 compartment pill box allowing patients the ability to fill their med device weekly or daily depending on the need. If a dose is not taken the “Missed Pill” message will appear giving the patient a visual indicator to further their medication adherence. Live a better life by having SlideBox help you in always remembering to take your prescriptions, supplements and vitamins as prescribed.

FEATURES

- Up to 37 daily auditory alarms
- ‘Missed Pill’ indicator
- Attached 7 compartment pill box
- Medical ID cover card

Next Page

Set-Up

NOTE: If this is the initial setup, please remove the clear battery insulator strip located above the clear 7 day pill box to activate.

CLOCK - Setting the Clock:

1. Press & hold the SET button for 3 seconds until the hour start to blink.
2. Press the ADJ button to set the correct Hour. Please observe the PM indicator.
3. Press the SET button until the Minutes start to blink.
4. Press the ADJ button to set the correct Minutes.
5. Press the SET button to finish setting the Clock.
6. If necessary press the ADJ button after setting the time to switch between 12 hour (AM/PM) format and 24 hour format.

ALARM TIME(S) - Setting the Alarm Time(s):

Note: The top row is used to set on-the-hour alarms only.

The bottom row is used to set half-hour alarms only.

To Set:

1. Open the clear cover.
2. For dosing times needing to be SET for on-the-hour, slide ONLY the corresponding Hour (top row) switch up.
3. For dosing times needing to be SET for on-the-half-hour, slide ONLY the corresponding 30 Minute (bottom) switch up.
4. Close the clear cover. Alarms will automatically repeat each day for the alarms set.

PILL BOX – Filling the 7 day pill box:

1. To Unlock: Push the white tab located on the left end of the SlideBox down until you hear a click.
2. To Open: Slide the white Reminder (top section) out towards the Right.
3. Place the clear 7 day pill box on the counter.
4. Fill the clear 7 day pill box according to your medication schedule taking note of the days of the week (M,T,W,T,F,S,S) written on the front edge .
5. To Close: Place the white Reminder (top section) on top of the Right end of the clear 7 day pill box and Slide the Reminder (top section) towards the Left.
6. To Lock: Push the white tab located on the left end of the SlideBox Up until you hear a click.

Reminder Time – How to use

ALARM Beep - What to do when the ALARM goes OFF:

1. Alarm will beep for 4 minutes.

2. TAKE your MEDICATION(S).

3. After you have taken your medication(s), press the STOP button.

Note: If the STOP button is not pressed within 4 minutes the “MISSED PILL?” indicator will be flashing on the display.

4. The alarm will go off again tomorrow at the same time.

‘Missed Pill’ Indicator - What to do when the ‘Missed Pill’ indicator is flashing:

1. TAKE your MEDICATION(S).

2. After you have taken your medication(s), press the STOP button.

Other

BATTERY – How to Replace the Batteries:

NOTE: If this is the initial setup, please remove the clear battery insulator strip located above the clear 7 day pill box to activate.

Replace battery when you hear a single beep alarm once every 15 minutes indicating low battery voltage. The low battery voltage indicator icon will also appear on the LCD display.

Normal battery life is 12 months.

1. Turn the device over and unscrew the two screws holding the battery compartment door.

2. Remove battery compartment and old batteries.

3. Insert two (2) new AAA batteries, take note of the (+/-) ends of the batteries.

4. Put the battery compartment door back on and tighten screws.

5. Press the reset button (gray button under clear cover) and set the time.

RESET – How to reset:

1. Open clear cover.

2. Press the reset button (gray button).

3. Set the time.

Disclaimer: This Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. RxPerfect, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. RxPerfect, LLC’s liability shall in no circumstances exceed the amounts paid by the user for the timer.

END