
R_xPERFECT

Perfect Dose. Perfect Health.

SlideMag User Guide



BENEFITS

Welcome to your brand new SlideMag prescription and vitamin reminder. SlideMag has been designed for those on a simple medication schedule wanting a stay-at-home solution. The up to 37 daily auditory alarms are set either on the hour or half-hour with a strong magnet adhered to the back of the device. This allows patients to secure their reminder to their refrigerator and as a result, have their reminder in a central location within their home. If a dose is not taken the “Missed Pill” message will appear giving the patient a visual indicator to further their medication adherence. Helps you Live Long and Healthy by making sure you never forget to take your prescriptions, supplements or vitamins again. *Due to SlideMag not containing an attached pill management system we recommend this device only for individuals living in a structured at-home environment.*

FEATURES

- Up to 37 daily auditory alarms
- ‘Missed Pill’ indicator
- Magnetic backing
- Medical ID cover card

Next Page

Set-Up

NOTE: If this is the initial setup, please remove the clear battery insulator strip located under the SlideMag.

CLOCK - Setting the Clock:

1. Press & hold the SET button for 3 seconds until the hour start to blink.
2. Press the ADJ button to set the correct Hour. Please observe the PM indicator.
3. Press the SET button until the Minutes start to blink.
4. Press the ADJ button to set the correct Minutes.
5. Press the SET button to finish setting the Clock.
6. If necessary press the ADJ button after setting the time to switch between 12 hour (AM/PM) format and 24 hour format.

ALARM TIME(S) - Setting the Alarm Time(s):

Note: The top row is used to set on-the-hour alarms only.

The bottom row is used to set half-hour alarms only.

To Set:

1. Open the clear cover.
2. For dosing times needing to be SET for on-the-hour, slide ONLY the corresponding Hour (top row) switch up.
3. For dosing times needing to be SET for on-the-half-hour, slide ONLY the corresponding 30 Minute (bottom) switch up.
4. Close the clear cover. Alarms will automatically repeat each day for the alarms set.

Reminder Time – How to use

ALARM Beep - What to do when the ALARM goes OFF:

1. Alarm will beep for 4 minutes.
2. TAKE your MEDICATION(S).
3. After you have taken your medication(s), press the STOP button.

Note: If the STOP button is not pressed within 4 minutes the “MISSED PILL?” indicator will be flashing on the display.

4. The alarm will go off again tomorrow at the same time.

‘Missed Pill’ Indicator - What to do when the ‘Missed Pill’ indicator is flashing:

1. TAKE your MEDICATION(S).
2. After you have taken your medication(s), press the STOP button.

Other

BATTERY – How to Replace the Batteries:

NOTE: If this is the initial setup, please remove the clear battery insulator strip located under the SlideMag.

Replace battery when you hear a single beep alarm once every 15 minutes indicating low battery voltage. The low battery voltage indicator icon will also appear on the LCD display. Normal battery life is 12 months.

- 1.** Turn the device over and unscrew the two screws holding the battery compartment door.
- 2.** Remove battery compartment and old batteries.
- 3.** Insert two (2) new AAA batteries, take note of the (+/-) ends of the batteries.
- 4.** Put the battery compartment door back on and tighten screws.
- 5.** Press the reset button (gray button under clear cover) and set the time.

RESET – How to reset:

- 1.** Open clear cover.
- 2.** Press the reset button (gray button).
- 3.** Set the time.

Disclaimer: This Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. RxPerfect, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. RxPerfect, LLC's liability shall in no circumstances exceed the amounts paid by the user for the timer.

END