

---

---

# R<sub>x</sub>PERFECT

Perfect Dose. Perfect Health.

---

---

## TimeCap & Bottle User Guide



### **BENEFITS**

Welcome to your brand new TimeCap & Bottle prescription and vitamin reminder. TimeCap and the included pill container have been designed as a per-prescription & per-vitamin reminder offering up to 24 on-the-hour auditory and visual alarms. The display indicates in large digits the day of the week along with the time of when you last took your dose. This reminder is ideal for those who just forget but is also great for live-in family members who need to be aware of the exact day and time their loved one last medicated. Live healthier by avoiding prescriptions, supplement, and vitamin errors. *For safety; doctors, nurses and pharmacists all strongly recommend using one and only one TimeCap & bottle per-prescription and per-vitamin. We see patients purchase multiple TimeCaps with the included pill container every day to insure they live longer, healthier lives by always being on a perfect dosing schedule.*

### **FEATURES**

- Up to 24 daily auditory and visual alarms
- Time indicating when the last dose was taken
- Included bottle to go along with the TimeCap

Next Page

## Set-Up

---

### **CLOCK - Setting the Clock:**

1. Remove the cap from the bottle.
2. Press the MODE button once and the display will start to blink.
3. Use the DOWN or UP buttons to change the time.
4. When desired time is reached, press the MODE button once and the weekday indicator will start to blink.
5. Use the DOWN or UP buttons to change the day of the week.
6. When desired weekday is reached press the MODE button once to store the settings and return to the time display.

### **ALARM TIME(S) - Setting the Alarm Time(s):**

1. Remove the cap from the bottle.
2. Press and hold the DOWN and UP buttons at the same time until the display blinks 12:-- PM.
3. Use the DOWN or UP buttons to reach the hour you wish an alarm to go off.
4. To set the alarm, press the CHECK button and an alarm symbol will appear in the upper left corner.
5. Repeat steps 3 & 4 to set additional alarms.
6. To remove an alarm, go to the selected alarm hour and press the CHECK button and the alarm symbol will disappear.
7. Press and hold the DOWN and UP buttons at the same time to store the settings and return to the time display.

## Reminder Time – How to use

---

*NOTE: The display on your TimeCap will always record the time that the cap was last opened to remind you of when you last took your medication, even if no alarms are set.*

### **ALARM Beep - What to do when the ALARM goes OFF:**

1. Remove the cap.
2. TAKE your MEDICATION.
3. Put the cap back on tightly until you hear a single 'beep'
4. The time and day of week of when you last took your medication (and put the cap back on) will automatically be recorded and displayed on the TimeCap.
5. If an alarm is overdue due to failure to take your medication the display will continue blinking; showing the last time a dose was taken. This is to remind you of when you last took your medication and the display will keep blinking until the TimeCap is opened and the dose is taken.

### **CHECK Current Time – How to check the current time set on the TimeCap:**

- 1.** Press and hold the CHECK button, and the display will show the current time and day of week.
- 2.** Let go of the CHECK button for the display to return to the last time you opened the bottle.

### **CHECK Alarm Time(s) – How to check the alarms time(s) set on the TimeCap:**

- 1.** DO NOT remove the cap from the bottle.
- 2.** Press and hold the DOWN and UP buttons at the same time until the display shows 12:-- PM.
- 3.** Use the DOWN or UP buttons to advance by one hour. The alarm symbol will appear on the display at each hour an alarm is set.
- 4.** Wait 5 seconds for the display to return to the last time you opened the bottle.

## **Other**

---

### **BATTERY – How to Replace the Batteries:**

Replace battery when you hear a single beep alarm once every 15 minutes indicating low battery voltage. The low battery voltage indicator icon will also appear on the LCD display. Normal battery life is 12 months.

- 1.** Turn the device over and unscrew the two screws holding the battery compartment door.
- 2.** Remove battery compartment and old batteries.
- 3.** Insert two (2) new AAA batteries, take note of the (+/-) ends of the batteries.
- 4.** Put the battery compartment door back on and tighten screws.
- 5.** Press the reset button (gray button under clear cover) and set the time.

Disclaimer: This Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. RxPerfect, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. RxPerfect, LLC's liability shall in no circumstances exceed the amounts paid by the user for the timer.

END