

CLAREBLEND MINI MICROCURRENT WARRANTY & WARNINGS

Upon receipt of your device, please go to <https://www.clareblendmini.com/warranty-submission> and register your purchase in order to ACTIVATE warranty. You MUST register in order to receive warranty repair or replacement. Device has a 6 month limited warranty. Warranty is void when device is transferred to another person, has been misused, subjected to accident, has improper maintenance or repair, or has undergone an unauthorized modification. If the Clareblend website malfunctions, email staff@clareblend.com. Christine Byer Esthetics does not own or maintain that website. **FAQs like how to know if the mini is working can be found at: <https://www.clareblendmini.com/faq>**

The MINI facial procedure is intended for facial stimulation and should not be used if the following indications are present: Pregnancy, seizures, cancer or tumors, cardiac pacemaker, or any electronic implanted defibrillator, simulators, or electronic devices as possible interference with these devices may occur. Do not use if prone to phlebitis, thrombosis, or on varicose veins. If you have any questions about use with a medical condition, please consult the applicable healthcare professional. We cannot and do not provide medical advice.

Use MINI only on normal skin. Do not use on skin that is swollen, infected, or has cancerous lesions. MINI is designed for cosmetic use only and for individuals in good health. As with any medical condition or after facial surgery of any kind, please consult your physician before use. MINI is for adult use only and should be kept out of the reach of children.

Wait a minimum of 10 days before using the MINI after the following procedures: Botox, injections and/or fillers, micro-needling (use after skin has fully healed).

To clean, wipe with a damp cloth. DO NOT SUBMERGE IN WATER.

BEFORE CONTACTING US FOR REPAIR OR REPLACEMENT, PLEASE TROUBLESHOOT AS FOLLOWS:

If no tone is heard when touching both electrodes to the skin, apply more gel; Ensure plug is secure in the adapter or power source; Ensure both electrodes are touching the skin. Movements must be SLOW and you must wait for the single beep before starting movements, and then wait for double beep (8-10 seconds) at end of movement. Rushing this will cause Mini to beep erratically. Unplug and start again if you mess up the beeping system.

When using the MINI, current should be sub-sensory, though some may experience a slight tingling sensation. Protective circuitry protects the user and internal circuits from surges of alternating current or power source failure.

No adverse effects have been reported from MINI use. Mild redness from the stimulation of circulation may occur. If irritation or burning is experienced at any time, discontinue use of the device.

Clareblend conductive gel is part of the skin enhancement process. Read all gel ingredients prior to application. Do not use conductive gel if any ingredient(s) may trigger a known allergic reaction. For a patch test, apply a small amount of gel on the skin for 15 minutes or longer. If allergic reaction occurs or develops after 24 hours, discontinue use and contact your skin care professional for alternative product suggestions. Apply gel as needed to maintain smooth glide and conductivity.

Any harm, to any area of the body or face, resulting from the misuse or application, because of improper voltage settings, contamination of conductive gel, lack of cleaning electrodes, or other improper applications, are not the responsibility of Clareblend, Christine Byer Esthetics, or any affiliates.

When procedure is complete, unplug device from power source, clean with damp cloth, and store. Apply sunscreen if going outside.

CLAREBLEND MINI MICROCURRENT INSTRUCTIONS

Cleanse the entire face and neck area. Apply a thin layer of conductive gel to one side of the neck & jaw line. Touch electrodes to the skin at the base of the neck & wait for a single beep. SLOWLY glide upward toward jaw line. Maintain electrode contact with the skin during movement & HOLD at the end of the movement until a double beep is heard. Move to the next section & repeat 3 times as shown. Avoid thyroid and the center area of the neck as indicated as a rectangle on the diagram.

Position MINI near the center of the neck and touch the electrodes to the skin. Wait for a single beep, then move the device SLOWLY toward the side of the neck. At end of movement, HOLD until a double beep is heard. Repeat movements as shown, working each section of the neck 3 times. Avoid thyroid and the center area of the neck as indicated in rectangle on the diagram.

Apply a thin layer of conductive gel to the jawline and cheek area. Position MINI at the jaw line and touch both electrodes to the skin. SLOWLY glide the MINI up and under the cheek bone. HOLD at end of movement until a double beep is heard. Repeat steps 3 times. Move to the opposite side of face and repeat jaw line steps 3 times.

Touch electrodes to the skin near the jaw line and SLOWLY glide across the skin toward the side of the face. Bottom movement follows the jaw line with one electrode above & one electrode below the jaw line. Maintain contact with the skin and HOLD at the end of the movement until a double beep is heard. Move to the next section as shown by the diagram by arrows and repeat until each movement has been completed 3 times.

Step 1



Step 2



Step 3



Step 4



Carefully apply conductive gel around the eye, temple and brow area. Touch electrodes to the skin near the nose waiting for a single beep before moving MINI toward the ear, continue movement gliding upward treating crows feet area. At the end of the movement HOLD until a double beep is heard. Repeat 2 more times, and repeat steps on opposite side of face.

Touch electrodes to the skin at the base of the brow bone. Wait for a single beep, then lift at the brow and HOLD until a double beep is heard. Repeat three times in each section as arrows show on the diagram. All eye movements should be on the orbital bone (never on the eyelid). Carefully reapply conductive gel as needed.

Apply conductive gel to one half of the forehead. Touch electrodes to the skin between the brows and wait for a single beep. Begin SLOW movement between the brows, moving up and over brow toward the side of the forehead. HOLD at the end of the movement until a double beep is heard. Movements for each half of forehead should be completed 3 times.

Starting at the top of the brow wait for single beep. Begin movement at the top of the brow and glide up toward the hairline. HOLD at end of the movement until a double beep is heard. Complete each movement 3 times. Repeat steps on opposite side repeating 3 times. When finished with all movements wipe face with damp cloth if needed. Unplug MINI, clean electrodes with damp cloth & store in safe place.

Step 5



Step 6



Step 7



Step 8

