

IMPORTANT SAFETY PRECAUTIONS & USAGE GUIDE

DO NOT USE THIS DEVICE IF ...

- You have recently tanned or sun-bathed.
- You have any type of skin irritation, rash, infection, burn, inflammation, active acne breakouts, rosacea, eczema, etc.
- You are recovering from a cosmetic procedure.
- You are using skin-thinning medications and/or topicals such as, but not limited to, Accutane, Tretinoins, and/or steroids. This could result in negative reactions in the skin.

READ ALL INSTRUCTIONS BEFORE USING

When using electrical appliances, basic safety precautions should be followed, including the following:

- Do not charge this device with any other charger.
- Only use the charger that was supplied.
- Do not place or store the device near water, such as a sink or tub.
- This device should be kept out of the reach of children.
- Keep the AC adapter and device dry.
- To clean after use, simply wipe with a damp (not wet), soapy cloth.
- DO NOT PUT THIS DEVICE UNDER RUNNING WATER OR IMMERSE IN WATER.
- NEVER OPERATE THE DEVICE WHILE IT IS PLUGGED IN.
- Never operate this device if it has a damaged cord or plug.
- Do not wrap the cord around the appliance.
- Remove contact lenses before using.
- If you are experience any negative reactions, discontinue use, and if you have serious skin problems, contact your physician or dermatologist before using this device.

Contents: Eye Rejuvenator Device | User Manual | Charging Cord

MAINTAIN CONTACT WITH SILVER CONTACT PLATE THAT SPANS 3/4 OF THE SIDE OF DEVICE TO MAINTAIN CURRENT.



BLUE LIGHT INDICATOR

USB CHARGING CORD



TITANIUM PROBE WARMS TO BODY TEMPERATURE



Christine Byer Esthetics | Website: <https://ChristineByerEsthetics.com>

Device has a 2-year warranty. If you have any issues with your device, please contact Support@ChristineByer.com to initiate the replacement process. You must provide your order number to verify your purchase from CBE.

IMPORTANT SAFETY PRECAUTIONS & USAGE GUIDE

HOW TO USE:

1. Charge your Eye Rejuvenator Microcurrent device for 3-4 hours PRIOR TO first use. Securely insert the USB end of the Charging Cord into a USB port OR you may use a cell phone USB adapter to plug into wall socket.
 - You may buy a power bank/charger compatible with a 5.0V device. The amount of mAh indicates how often you will have to recharge it. Using 10,000 mAh would be sufficient for this device. Over 40,000 mAh and you lose portability in the size of the power bank.
2. The red light will flash to indicate the battery is charging. When fully charged, the red light will stay on.
3. Once charged, unplug the cable. Do not keep the charging cable connected to a power source when not charging the device or while operating the device.
4. NEVER OPERATE THE DEVICE WHILE IT IS PLUGGED IN.
5. Once the battery begins to run low again, the red light will flash repeatedly to alert you it is time for a re-charge.
6. Turn on by holding down the on/off button for until it beeps and release to turn on. Once the BLUE indicator stops flashing, the device is ready to use. NOTE: The tip will heat and feel very warm.
7. Make sure skin is freshly cleansed, then apply the eye serum or cream (water-based, do not use oil-based – oil prevents conduction) of your choice around the eye area— eyelid, corner of the eye, and under the eye.
8. Hold the device, making sure your fingers are making direct contact with the silver contact plate (this maintains the microcurrent connection) and slowly and gently glide the device from the outer corner to inner corner of your eye, pausing for 3-5 seconds on acupoints (see diagram below). When the device beeps (at about 1 minute 30 seconds, switch to the other eye. It will shut off in about 1 minute 30 seconds. Total 3 minutes for the entire treatment. The device will turn off by itself after treatment. You can use this device daily.
9. Follow up by proceeding with your regular facial routine, such as additional facial serums, moisturizer, and daily SPF.



Christine Byer Esthetics | Website: <https://ChristineByerEsthetics.com>

Device has a 2-year warranty. If you have any issues with your device, please contact Support@ChristineByer.com to initiate the replacement process. You must provide your order number to verify your purchase from CBE.