



Z-SMART

Foot Massager Plus

USER MANUAL




Features:

1. Full Airbag double feature design with hinge cut-over feature on top of the device to allow full foot, calf and thigh length massage.
2. Airbag and roller massage offer multiple pressure varieties, with soft to strong intensities and stereoscopic targeting. Massage techniques are innovative and flexible allowing for a wide range of targeting options.
3. Three levels of intensities, five varieties of massage programs and 17 different massage combinations provide comprehensive massage options.
4. Multi-angle massage is possible with the 90° adjustable hinge system at the top of the unit and 20° adjustability in the base of the massager.
5. Jade Heating is included for comfort.
6. Full range of massage points from the foot to the calves to the top of the thigh.
7. Built in safety features such as auto deflate in the case that power is cut to the device the massage airbags will not hold.
8. Auto-Off function to extend the longevity of your device. Massage routines last 20 minutes at a time.
9. Easily removable and washable kneading cloth cover so your massager is fresh every time you use it.



Features Continued:



Power On/Off button

Power



Heat



Intensity

Massage strength from soft, medium and strong.



Foot



Leg

Heating On/Off button

Leg (Calf) Massage On/Off button

Foot Massage On/Off button




Revive




Reflexology

Reflexology Massage On/Off button

Revive Massage On/Off button



Hinge System



Lower Frame Adjustment



Easy Transport

Unit Diagram and Recommended Care

Avoid excessive temperatures especially if kept in a car. Avoid high temps during the summer and cold freezing temps during the winter.

Always take care to store in a clean and dry place free from dust, dirt, moisture and any corrosive chemicals.

To Power on plug, the power plug into the power port on the side of the massager. The power indicator light will flash and then the foot massager will enter standby mode.




• Power Plug Diagram

Notes:

Please switch power off and unplug power plug if not in use.

Do not stand during the massage before or after use to avoid damage to the machine.



Power On/Off switch

Usage Information



Power

Power on/off switch with indicator, doubles as reset switch.



Heat

Jade Heating function by Zarifa USA, lights up blue while in use.



Intensity

Three strengths included once massage routine begins. Green indicates soft massage, blue indicates medium massage and red indicates strong massage.



Foot

Lights blue when in use and provides a comprehensive bottom of the foot massage.



Leg

Lights blue when in use and provides a comprehensive calf, ankle to thigh



Revive

Lights blue when in use and provides a full Foot and Leg massage routine.



Reflexology

Lights blue when in use and provides a bottom of the foot and calf massage

Cleaning Instructions and Product Specifications:

- Unplug before beginning the cleaning process.
- Use a dry or minimally damp cloth to clean away dust, dirt, and lint from the massager.
- Do not use chemical cleaners on the product.
- Do not dry in direct sunlight for excessive periods of time.
- Do not place near open flame or heater.

Product Specification:

Name: Z-Smart Foot Massager

Model: Z801c19

Rated Voltage: 110 – 120V~ 60Hz 220 – 240V~ 50Hz/60Hz

Rated Power: 100W

Rated Working Time: 20 Minutes

*Massager is not waterproof image is for creative purposes only



Troubleshooting Tips

Troubleshooting Tips:

In rare cases that something goes wrong with your massager, the following are commonly encountered issues and how to solve them. Please review all procedures before calling our office for further instruction.

The Massager Is Not Working:

Problem: Not Plugged

Solution: Plug the product in

Problem: The Massager has Auto Shut off feature after a 20-minute cycle.

Solution: Press the On/Off button

Problem: The product has been overused past the 45 min recommended use time

Solution: Let product cool for 1-hour, especially if hot to the touch

Problem: There is too much weight placed on the massage heads (Such as standing on the massage units)

Solution: Remove heavy load from the program

The product is Hot to the touch (more than warm)

Problem: Continuous Over-use of the product has led to its overheating.

Solution: Unplug the device for 1 hour before use once more.

Problem: The product is being used in extreme temperatures

Solution: Place the product in a cooler area for 1 hour before use once more.

Soft Massage Sound Coming from Product:

Problem: This sound is normal and is quieter than any other massager on the market. It is caused by the massage motors.

Solution: The sound is normal and not defective.

*In the case of more severe issues such as defects, excessive noise and other related issues call our offices at Zarifa USA.

Prohibited Uses:

Do not use the massager in a highly humid environment such as a bathroom or sauna.

Do not use the massager immediately after a sharp change in temperature.

Do not use this massager in an industrial environment or a place with highly corrosive gasses and/or dust build up.

People Who Shouldn't Use Product:

Pregnant Women

People with inflammatory or rare skin conditions.

Patients with osteoporosis

Patients with heart disorders and/or Patients with a pacemaker.

Patients with a high fever

Minor's under the age of 18.

People who do not feel well or have low energy.

Those who are wet.

People under the influence of drugs or alcohol.

Consult with your doctor if you have recently undergone surgery.

People Sensitive to Heat.



www.zarifausa.com



Bonneville Business Center
Zarifa USA
461 West Century Drive
Suite D3-A
Murray, UT
84123

