





# ZARIFA USA

## Z-Fold

### Folding Shiatsu Massage Chair



Zarifa Usa  
(385) 645-0255  
465 W Century Dr #3  
Murray, UT 84123

User Manual  
M25

## Safety Instructions



**Please read the following instructions carefully before putting your massage device into use. Please retain these operating instructions for further use!**

- The massage device conforms to recognized technical principles and the latest safety regulations.
- Do not get wet, do not use pins, never remove chair cover.
- This item is NOT A TOY. Close supervision is necessary when this appliance is used by , or near children or disabled persons.
- This appliance should never be left unattended when plugged in.
- Never use any power sources other than the home adapter provided as original equipment with this machine.
- Any possible repairs may only be carried out by authorized expert staff, improper use and unauthorized repairs are not permitted for safety reasons and lead to loss of warranty.
- Never touch the power plug with wet hands.
- Avoid contact of the device with water, high temperatures and direct sunlight.
- Do not use any damaged cables, plugs or loose sockets.
- If plugs or cords are damaged, they must be replaced by the manufacturer, a service representative or by qualified personnel.
- In case of malfunction, disconnect immediately from the main power.
- Do not use if you have skin disorders, open wounds, or swollen or inflamed areas.
- Misuse or incorrect use excludes any liability for damage.
- Do not use while sleeping.
- In order to avoid excessively stimulating the muscles and nerves, the recommended massage time should not be exceeded 15 minutes at a time.
- Do not use this product in place of medical treatment.

### **Note!**

- Each massage – even a hand massage – must be refrained from during pregnancy or if one or more of the following complaints are present in the massage area: recent injuries, thrombotic diseases, all kinds of inflammations and swellings, and cancer. It is recommended to consult a doctor before having a massage for the treatment of ailments and illnesses.
- If you are dependent on electrical aids e.g. pacemakers, please be sure to consult your doctor for medical advice before taking a massage.
- The packaging material supplied should not be used as a toy.

Failure to observe the above instructions may constitute misuse of the product and can cause serious injury or burns.

## Instructions

Thank you for purchasing Folding Shiatsu Massage Chair.

With normal care and proper treatment, it will provide years of reliable service.

**Please read all instructions carefully before using this product.**

**Please retain this manual for future reference.**


## FEATURES

- Electrically height-adjustable Shiatsu neck massage
- Kneading Shiatsu back massage
- Customized Area back Massage
- Optional heat function on neck & back massage
- Seat vibration massage (3 intensity levels for options)
- Premium PU leather and breathable mesh cloth
- Adjustable backrest height
- Innovative no-assembly design sets up easily and conveniently folds away for storage

## CONTENTS

---

1. Folding shiatsu massage chair
2. 12V home adapter




## TECHNICAL DATA

---

<b>Measurements:</b>	18.5 x 29.5 inches (Back) 18.5 x17.3 inches (Seat)
<b>Weight:</b>	16 lbs (Net)
<b>Voltage:</b>	Input: AC 100-240V ~50/60Hz Output:12VDC 5.0A
<b>Nominal Power:</b>	max. 60Watt
<b>Default Automatic runtime:</b>	15 minutes


## Setup and Operation

1. Place the massage chair on the flat floor.
2. Open the folding chair (Figure 1).



3. Turn the handles on both sides of chair to unlock point (Figure 2)
4. Adjust the backrest to the desired position.
5. Rotate the handles back to lock position and make sure the backrest is stable to lean on.

Fig.2



**Unlock**

**Lock**

6. Connect the adapter cable to the corresponding cable on the massage chair
7. Plug the home adapter into an electric outlet between 100-240V.
8. Turn on the device using the controller ( instructions on page 6 ).
9. Folded away for storage when not in use for a long time (Figure 3).



## Setup and Operation

### Note !



1. Detachable flap for back lets you choose a softer or more intense back massage.
2. Detachable cover for neck lets you choose a softer or more intense neck massage.

### If the Shiatsu is too strong

- ① Place a blanket or towel between you and the massager to achieve a far softer massage intensity if you still feel shiatsu too strong after installing the flap onto the chair

### If the shiatsu can't reach your neck,shoulder,lower back area

- ② If you want to massage your neck, shoulder,lower back, but the shiatsu nodes can not reach the area in need, you may adjust the backrest height, or place a seat cushion or a pillow on the seat.




- The device will automatically switch off after 15 minutes when timer runs out
- Remove the adapter from the outlet when not in use.
- Do not setup or use the massager in a bathroom or similiary wet / damp areas.

### ! Attention

**You should allow your body some rest periods. To avoid straining your muscles too much, we recommend to not exceed a continuous massage of 15 Minutes per time and not exceed 3 times use one day.**

## Controller Instructions



### Note!

1. Press the power button, then only press the "HEAT" button once to activate the heat function on both the neck and back area, you may adjust the heat nodes around the neck up or down by pressing the arrow buttons on the neck area; Pressing the arrow buttons on the Back area would adjust the position of the heat nodes on the Back (Up & Down) and simultaneously activate the Shiatsu Function on the back.
2. After the whole massage experience is completed, press the Power button again to turn off the massager. The back massager rollers would be automatically reset and return to the lowest position for storage.

## Care and Cleaning Instructions

1. Do not wash the machine or immerse the product or the adapter in water.
2. Always switch off before cleaning and remove the adapter from the socket with clean, dry hands.
3. The massager should be wiped with a slightly moist cloth and dried thoroughly. Do not wipe the product with alcohol or harsh chemicals such as thinner or bleach.
4. Keep away from heat, fire or exposure to direct sunlight.
5. Do not apply excessive pressure or impact onto massager/massage heads. Application of a direct force in excess of 110 lbs. may damage your massager and invalidate your warranty. Refrain from over usage of the device.
6. Avoid contact with sharp edges or pointed objects which might cut or puncture the cover of massager.
7. Folded away for storage when not in use.
8. Store the massager in a cool, dry place.

## TROUBLESHOOTING GUIDE

Fault	Solutions
Massage heads rotating at reduced speed.	The load on the massage heads is excessive. Reduce load and try again.
Massage heads only rotate in the upper or lower region.	The massage has been activated for the upper or lower back region. Press the 'Full Back' button for the massage to cover both regions.
Device switches off suddenly	Device may have been in use for 15mins. Auto-shut-off function activated. Continue using by pressing the 'power' button after resting the unit for 15mins.
	Over-heat protection function activated. Device is switched off for your safety. Please rest device for 15-30mins before restarting.
Device fails to start	Check if device has been used for more than 15mins. If so, let device rest for 15mins before continuing. Device may require longer time to cool in warmer environments.
	Check if adapter is plugged into socket.
	Check if adapter output is connected to device power inlet wire.
Jerking motion of massage heads / Vibrations from device	Excessive pressure may have been applied onto the device. Reduce the load and try again.
Infra-red lights are on but no motion from massage heads	Only activate the heat function, press other buttons to activate the shiatsu function.