



ZARIFA USA

HOME SUITE

Shiatsu Massage Pillow

USER MANUAL



SPEED



HEAT



DIRECTION



POWER

USER MANUAL

PRODUCT USE BY ENVIRONMENT

- Do Not use this product in a damp or dusty environment to avoid electric shock
- Do Not Use this Product in Extreme Temperatures


CONSULT A DOCTOR BEFORE USE IF YOU ARE

- Receiving Medical Treatment
- A Cancer Patient
- Pregnant Woman
- A Person with a Heart Condition
- Osteoporosis Patients
- A Person with A Pacemaker




ACCESS IBI LITY:

- Do not pull off adapter with wet hands to avoid electric shock hazard
- When you pull off adapter, grasp the whole adapter instead of just the wire
- The timing function can automatically stop the massage cushion, please use the Power switch after each use
- Adapter shall not be packed by the product or other heavy objects
- If the adapter is damaged, immediately stop using the product, turn off the power Switch, pull the adapter out from the outlet and contact a qualified service Personnel
- When it is not in use for a long time, turn off the power switch







POWER ADAPTER



CAR ADAPTER



Keys

| | |
|--|--------------------------------|
|  POWER | Turn on and off the machine |
|  DIRECTION | Change the direction |
|  HEAT | Turn on the heating on and off |
|  SPEED | Change the massage speed |

- 1.) Insert the home power adapter into the port on the unit, then insert the power plug into a power outlet with 110-240V (50-60Hz). This will put the unit into standby mode.
- 2.) When using in your car, connect the dedicated cigarette lighter adapter to the 12V power source (car cigarette lighter port). This will put your unit into standby mode.
- 3.) Press the power button and the massage and heat function will start. If you press it again, all functions will stop.
- 4.) Press the heat button for the heating function to be turned on/off
- 5.) Press the direction button and the product will switch the direction of the rotation by itself while in use.
- 6.) Press the speed button and the product will change the massage speed. There is 3 different speeds you can use.
- 7.) The product has a 15 minute timer when in use. After 15 minutes of use, it will shut off.
- 8.) After massaging, please unplug the power and store your unit.

Massage Areas



Neck



Waist




Abdomen



Thigh



Arm



Calf

The full body Shiatsu massage cushion is able to perform a more in-depth massage of your neck, back, shoulders, legs, or waist. The heating treatment can relax your body, eliminate fatigue, relieve pressure, and enhance blood circulation. The compact/stylish design will not only meet household demand, but also while you are driving.

- 1.) After using the product, please turn off the power supply.
- 2.) Prevent the product from coming into contact with water or other corrosive chemical agents.
- 3.) The product should not come into contact with sharp/hard objects. It should not be pressed on hard, or stepped on.
- 4.) During use, if there is any abnormal situation, please cut off the power.
- 5.) Do not modify the power adapter or repair it.



7 98304 43032 5



461 West Century Drive
Suite D3-A
Murray, UT
84123

zarifausa.com

