

## Description

Our Light-Therapy Golden Facial Treatment device is a five-coloured light treatment that results in transformative radiant skin. The healing combination of the LED calms and rejuvenates the complexion while visibly diminishing hyper-pigmentation.

Each light performs a targeted treatment (150 diodes):

Red light (630nm) has profound anti-ageing benefits, it stimulates collagen and elastin production and strengthens skin. This light setting also reduces inflammation, heals blemishes, repairs sun damage and improves scarring.

Blue light (430nm) balances and clarifies the skin, treating acne. With repetitive use, this light improves blemished skin.

Green light (520 nm) calms skin and reduces pigmentation, broken capillaries and sunspots.

Yellow light (570-590nm) reduces redness and boosts circulation. It is ideal for sensitive skin as it boosts lymphatic flow and soothes the skin.

White light (450-460nm) penetrates deeply to promote wound healing and skin repair.

## THIS BOX INCLUDES:

1 x Gold face mask with adjustable strap

1 x Power cable

1 x Power control and remote

There are 3 plug options with this device, adjusted to the voltage in each country. Voltage for the EU plug corresponds to 220V, USA corresponds to 120V, UK corresponds to 230V.

## TECHNICAL SPECIFICATIONS:

Power: 15W

Max energy: 45J/Cm<sup>2</sup>

Output intensity: 35 mw/cm<sup>2</sup>

Frequency: 50Hz for UK & EU, 60 Hz for USA

Please note the Light-Therapy Golden Treatment Device is non-returnable unless faulty. If for whatever reason you are not happy with your purchase, you have 14 days to return your device. You will be eligible for a full refund if the item is new and unused,

with all the original accessories, packaging and manufacturing seals intact. Please note that we are unable to refund any delivery or postage costs incurred unless the device returned is faulty. The Light-Therapy Golden Facial Treatment Device has a 6-month warranty which covers the device should you experience any faults. MZ Skin reserves the right to repair or replace the device during this time.

## How to Use

Apply mask to cleansed, dry skin. Use 2-3 times per week, starting with 10 minutes of LED mask time and slowly building up to 30 minutes. Discover optimum results after around a month of regular use.

Each light performs a targeted treatment, so you can select the light you would like to use for a session. You can either use a single light to target a specific skin concern or you can use a combination of lights throughout the week. The device comes with an adjustable headstrap or it can simply rest it on your face if you're laying flat.

## WARNING:

Individuals suffering from eye-related disorders e.g. retinitis pigmentosa, ocular albinism, other congenital retinal disorders), repeated exposure may cause varying degrees of retinal damage that could be irreversible and could accelerate peripheral vision impairment or loss. Please consult your doctor before using.

Please do not exceed recommended usage time as prolonged exposure could cause skin or eye damage.

Individuals who are prone to migraines may want to approach this product with caution as bright lights can potentially trigger a migraine. At your own discretion, those that suffer from migraines may try this product for a few minutes at a time and gradually build up the time.

Do not use this LED mask if you are pregnant, have epilepsy, cancer, haemorrhage or if you are taking medications including, but not limited to tetracycline or steroids. In individuals suffering from eye-related disorders e.g. retinitis pigmentosa, ocular albinism, other congenital retinal disorders), repeated blue light exposure may cause varying degrees of damage that could be irreversible and could accelerate peripheral vision impairment or loss. Please consult your doctor before using.

There is no UV exposure with this treatment. LED does not increase sensitivity to sunlight but SPF is always recommended.

## Q. CAN THIS BE USED DURING PREGNANCY?

A. We would advise against the use of LED Mask if you're pregnant. While there are no known harmful effects of this, there also isn't any concrete research that suggests that the LED Device is safe to use during pregnancy, so as a brand we always tend to err on the safe side.

#### Q. HOW OFTEN CAN I USE THE LED MASK?

A. Use 2-3 times per week, starting with 10 minutes of LED mask time and slowly building up to 30 minutes. Discover optimum results after about a month of regular at-home treatments.

#### Q. IS THIS FDA REGISTERED?

A. Yes

#### Q. CAN I KEEP MY EYES OPEN DURING THIS TREATMENT?

A. We recommend covering your eyes with a cotton pad as repeated blue light exposure can cause varying degree of damage.

#### Q. HOW LONG IS THE GUARANTEE?

A. 6 months (Note from Emma: this is the manufacturer's warranty - we offer 2 years)

#### Q. ARE YOUR PRODUCTS CRUELTY FREE?

A. MZ Skin does not test on animals or commission any third parties to do so on our behalf. We do not use ingredients or suppliers who test on animals. Any animal-derived ingredients used in MZ Skin formulas, such as ovine placenta, are natural by-products and no sheep are harmed in the process.