

gandio™



## Treatment - Simple, Fast and Safe

A daily use of 5 minutes for users with short hair and 8 minutes for those with longer hair will bring a significant improvement within several months. It is recommended to use 1-2 times daily, however, for faster and better results use it twice a day for the first two months.

Once major potential has been achieved making hair appear thicker and fuller you may reduce the treatment to 3-4 times per week.

1. Lift the Applicator from the Base.
2. Lightly press the Applicator's button. After a quick self-test, the Applicator will emit a startup tune and start vibrating.
3. For effective treatment, roll the Applicator back and forth over the scalp area requiring treatment.
4. Profile II - 5 minutes treatment - recommended for short hair
5. Profile III (8 Minutes ) is intended for people with long or dense hair. For such people, it is recommended to part the hair and roll back and forth across the parting line for half a minute, then move on to another parting line.
6. If you find the intensity and sensitivity too strong, you can use profile I to let the scalp get use to the treatment.
7. The Base and Applicator should be cleaned with a damp cloth and then dried. The disks and the area between them can be cleaned with a brush. For sanitary purposes, roll the disks on an alcohol pad, and dry before use.