



Perfectio X™



1. Make sure the Perfectio X is fully charged before use (Green light indicator).
2. After you turn your Perfectio X on, it will start heating up slowly (no need to wait). Place the glass probe on your skin, then the LED light will automatically turn on.

***For pain management treatment:**

3. it is recommended to use the Perfectio X on clean skin for 10 minutes per treated area along with a small amount of the provided Lenio Red Algae Pain Relief Gel directly on the glass probe.

4. Treat each area of pain once or twice daily for 10 minutes each time as necessary to eliminate the pain. Red algae is well-known as a miracle ingredient full of incredible therapeutic benefits for the skin. It conditions, hydrates, softens and detoxifies the skin, replenishing essential vitamins and minerals, as well as reducing pain and inflammation.

5. Clean the Perfectio X after every use with a dry or slightly damp cloth, a cotton ball, or a paper towel.

***For skin rejuvenation treatment:**

3. it is recommended to first use the Perfectio X on clean skin for up to 3 minutes per treated area, then apply your preferred facial cream or mask and use the device for additional 90 seconds.

4. For best results, use the device 3 times a week in the first month, then use it as often as needed – usually once or twice a month (results differ by age).

5. Clean the Perfectio X after every use with a dry or slightly damp cloth, a cotton ball, or a paper towel.