



sapphire X™



Treatment - Simple, Fast and Safe

For best results, use it 3 times a week for the first month. After the first month, use as needed, normally once a week.

1. Make sure the Sapphire X is fully charged before use (Green light indicator).
2. After you turn your Sapphire X on, it will start heating up slowly (no need to wait).

Place the glass probe on your skin, then the LED light will automatically turn on.

1. Apply light pressure while slowly moving the Sapphire in small circular motions.
2. Treat each area for 3-5 minutes - face, forehead and around the eyes.
3. Do not use the device over the eye lids.
4. It is recommended to apply a small amount of moisturizing cream on the treated skin areas after each session and to massage it gently in circular motions.
5. Clean Sapphire X after every use with with a dry or slightly damp cloth, a cotton ball, or a paper towel.
6. Stop using if you feel the treatment head is too hot.