

Neck & Facial Beauty Device

HOW TO USE

USING THE MY SMOOTH SKIN NECK & FACIAL BEAUTY DEVICE IS A SIMPLE PROCESS AND ONLY MINIMALLY TAKES UP TO 5 MINUTES OF YOUR SKINCARE ROUTINE. FOR THE BEST RESULTS WE RECOMMEND BEING CONSISTENT FOR AT LEAST 4 DAYS A WEEK. BELOW ARE STEP-BY-STEP INSTRUCTIONS AND GUIDANCE ON HOW TO GET THE MOST OUT OF YOUR BEAUTY DEVICE.

BEFORE YOU START

- PLEASE ENSURE THAT YOUR DEVICE IS FULLY CHARGED PRIOR TO USING IT, THIS DEVICE TYPICALLY TAKES APPROXIMATELY **2 HOURS TO ACHIEVE A FULL CHARGE.**
- ONCE YOU ARE READY TO BEGIN, CHOOSE THE LIGHT SETTING THAT BEST SUITS YOUR **CURRENT SKIN CONCERN NEEDS.**
- ONLY USE YOUR DEVICE ON CLEAN SURFACES, INCLUDING BOTH YOUR **HANDS AND FACE. CLEANSE, THEN MASSAGE.**

READY TO USE!

* Now that your hands and face are clean, begin by taking oil, cream, serum of your choice and apply section by section of the area you want to target.

*Pro Tip! Reapply product as you go to ensure a smoothing gliding motion over each area!

* Be sure to remove the white cap off of the device, and wipe the unit to ensure a clean surface. Once completed, press and hold the power button to turn it on. This same button allows you to switch modes by pressing it to reach the desired setting. You will see the light change so you know which mode you are in!

* Now that you have your mode selected, you can begin gliding the device in an upward motion with repeated gentle strokes for about 6-10 times per area. The device will automatically shut off after 5 minutes of being in use.

* Once you are done having the most amazing and relaxing facial massage at home, you can clean the thermal head with a soft cloth and put the cap back on to protect the the device. It's really that easy!