

# Pro IPL Hair Removal Device

# HOW TO USE

*OUR IPL DEVICE IS USER FRIENDLY, AND HAS A SIMPLE TO UNDERSTAND MANUAL ATTACHED INSIDE THE PACKAGE. THIS GUIDE IS INTENDED TO HELP YOU MORE WITH HOW TO CARRY OUT YOUR HAIR REMOVAL TREATMENTS TO ENSURE A POSITIVE OUTCOME. AN IPL TREATMENT IS TYPICALLY 5 - 15 MINUTES. WE RECOMMEND DOING 2 TREATMENTS PER WEEK, FOR 12 WEEKS. TO MAINTAIN YOUR RESULTS THEREAFTER, DO A TREATMENT AT LEAST ONCE A MONTH. CONTINUE USE OF THE IPL, AS LONG AS NECESSARY.*

## BEFORE YOU START

- PREPARE YOUR SKIN FOR THE TREATMENT BY SHAVING ANY HAIR ON THE SURFACE OF YOUR SKIN. SHAVING WITH A RAZORS, OR MY SMOOTH SKIN DREAM HAIR ERASERS ARE THE ONLY RECOMMENDED METHODS OF HAIR REMOVAL.
- CLEAN TREATMENT AREA FROM ANY OILS, CREAM OR DIRT. EXFOLIATE IF POSSIBLE. THIS WILL MAXIMIZE HOW MUCH HEAT IS ABSORBED BY HAIR FOLLICLES.
- PLEASE ENSURE THAT YOUR DEVICE IS PLUGGED INTO A POWER OUTLET. A COMMON MISCONCEPTION IS THAT IPL UNITS ARE RECHARGABLE.
- SELECT YOUR PREFERRED HEAT INTENSITY AND GET READY TO ZAP.

## READY TO USE!

Place the IPL device flat against your skin. Press the zap button to emit the light pulse and move the device across the treatment area, ensuring overlapping sections slightly to cover the entire area. Go over the same area 4 - 6 times.

Split the area you're treating into smaller sections. This helps to ensure comprehensive coverage without missing any spots. Move the device systematically from one section to another, overlapping slightly for consistent exposure.

Certain areas might be more sensitive than others. Be mindful of these spots and adjust your approach accordingly. You might need to slow down or reduce intensity in sensitive areas to minimize discomfort.

---

## TIPS TO ENHANCE RESULTS

**Stay Consistent:** Adhere to the recommended treatment schedule without skipping sessions for optimal results.

**Adjust Intensity Carefully:** Start with lower intensity levels and gradually increase as tolerated for better comfort and safety.

**Patience is Key:** Results may not be immediate. It takes multiple sessions to notice a significant reduction in hair growth.

**Regular Maintenance:** Even after achieving desired results, occasional maintenance sessions might be needed to keep the hair growth in check.