# CLOUD NINE<sup>®</sup> BEAUTY

**REDEFINE**

NEXT-GENERATION SKINCARE

DETOX | REGENERATE | HEAL

CLOUD NINE® | BEAUTY

## WHAT'S IN THE BOX

CLOUD NINE REDEFINE  
USB CHARGER  
INSTRUCTION MANUAL  
DRAWSTRING BAG

BE FIERCELY YOU

CLOUD NINE® | BEAUTY

## CLOUD NINE REDEFINE

REDEFINE your routine with next-generation skincare. Inspired by Gua Sha - a massage technique used in Chinese medicine for thousands of years. Combining a unique shape designed to fit the contours of the face and body with advanced heat, light, and detoxing skincare technology to provide a nourishing, purifying, anti-ageing treatment from home.

CLOUD NINE Redefine lifts, firms and detoxes the skin from within. Purify and regenerate skin cells with lymphatic drainage, sculpt and tighten with electro muscle stimulation and target lines, blemishes and nourish the skin with Red, Blue, and Purple light therapy.

Own your glow and combine each setting to create your unique skincare treatment.

BE FIERCELY YOU

## USING YOUR REDEFINE

---

### ON/OFF

Press and hold the ON/OFF button to turn off

**GUA SHA:** Designed to mould to the contours of the face.

**DETOX** - Eliminate toxins in the skin cells

**BOOST** - Encourages collagen production and improves product absorption

**REGENERATE** - Penetrates deep into the skin to revive skin cells

**STIMULATE** - EMS lifts and tightens the skin

**PURIFY** - Unclogs the lymphatic and circulatory system

**REDUCE** - Red Light Therapy helps reduce fine lines and wrinkles

**HEAL** - Blue Light Therapy heals and soothes blemishes

## TREATMENT SETTINGS

### 1. Red Light, Vibrate & Heating

**2. Blue Light, Vibrate & Cooling** - with automatic shut off after 3 mins. We recommend leaving 10 mins between cold treatments for optimum results.

### 3. EMS & Purple Light

Apply your favourite skin care products such as serums, moisturisers and creams to use with the EMS setting. Ensure you are touching the conducting strip to activate the micro-current treatment setting. You will feel a tingling sensation and the intensity of the EMS felt will depend on skin type, thickness, and if products are being used. This is normal, safe and means the treatment is working.

### **LOW BATTERY**

When the battery is very low or out of power, the red light will flash fast.

### **CHARGE**

Charge using the USB cable provided. When your product is charging a red light will flash slowly. When fully charged the light will be solid blue.

### **HOW TO CLEAN YOUR REDEFINE**

---

- 1) Clean the treatment head with a cotton swab or soft cloth after use
- 2) Please do not use any of the following liquids to clean the device: alkaline lotions, diluents, essential oils, nail varnish remover, and other solvents
- 3) Do not immerse the main device into water

### **WARRANTY**

---

This CLOUD NINE Redefine device comes with a 12-month warranty. In the unlikely event that your product develops a fault, our warranty covers the replacement of your product for free within the guarantee period.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

## SAFETY PRECAUTIONS

### PLEASE READ CAREFULLY

- Before using the device, please read the manual carefully and keep it for future reference
- When you use the device for the first time, test it on your forearm before you use it on your face. Check your skin 1 hour after testing the treatment. If there is no irritation or reaction, you can safely use the device on your face
- Do not use if you have an allergy to metal
- Do not use the device on wounded or scarred skin
- Please store the device in a cool, dry place
- Do not allow the device to come into contact with water. Any water, liquid or foreign body entering the device will invalidate the warranty
- Do not use the device whilst charging
- Keep out of reach of children

**WARNING:** Do not use this appliance near bathtubs, showers, basins or other vessels containing water

This marking indicates that the product should not be disposed of with other household wastes. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources.

**CLOUD NINE® | BEAUTY**

## **OUR GIFT TO YOU**

Use this exclusive code to get 10% off your next  
CLOUD NINE styling tool or beauty device.

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CLOUD NINE International Ltd,
CLOUD NINE ANZ Pty Ltd,