CLOUD NINE[®]

BEAUTY

REJUVENATE

CLEANSE | BOOST | CLEAR | STIMULATE | REDUCE | HEAL

REJUVENATE YOUR SKIN FROM WITHIN

CLOUD NINE® | BEAUTY

WHAT'S IN THE BOX

Rejuvenate Beauty Device

Charging base

Power adaptor

Cotton cleansing pads

Travel pouch

Cotton cleansing pads ring

User Manual

BE FIERCELY YOU

CLOUD NINE® | BEAUTY

CLOUD NINE REJUVENATE

Rejuvenate your skin from within. CLOUD NINE Rejuvenate combines light and heat therapies with advanced facial cleansing technology.

Own your glow with advanced absorption technology, micro-current muscle massage, Red and Blue Light treatments, and Hot and Cold therapies. Use each function individually or combined depending on your skin type for a personalised at-home facial.

CLOUD NINE Rejuvenate is a home-use beauty tool for an effective, simple and safe treatment to deep cleanse, tighten and lift the skin and treat mild to moderate fine lines and wrinkles.

Rejuvenating technology triggers your skin's healing functions, resulting in the production of collagen and elastin fibers for a long-lasting youthful, radiant complexion.

BE FIERCELY YOU

REJUVENATE TREATMENT SETTINGS

Select mode to switch through each treatment setting in sequence:

CLEANSE - Select Clean to remove oil and impurities for a radiant texture*

BOOST - Select Import to stimulate collagen and improve product absorption*

CLEAR - Select Eye Care for gentle ion technology to reduce dark circles and puffiness around the eye area*

STIMULATE - Select EMS for Electro-Stimulating Massage to lift and tighten the skin

REDUCE - Red light therapy reduces fine lines and wrinkles

HEAL - Blue light therapy heals and soothes blemishes

Use your favourite skincare products with each setting to create your own signature treatment.

*Make sure you hold both sides of the ion electrode of the CLOUD NINE Rejuvenate Device when using "CLEAN", "IMPORT" and "EYE CARE" MODE for the current to activate.

USING YOUR REJUVENATE

1. CHARGE

Take the CLOUD NINE Rejuvenate from the box and remove the plastic adhesive films from the handle. Connect the power adaptor to the charging base then place the handle on it to charge. Continuously charge your Cloud Nine Rejuvenate for 3 hours before first use.

*Keep magnetic charging terminals clean and dry at all times. Product build up may result in the Rejuvenate not charging correctly. It is recommended to store on the charging unit at all times.

2. FIRST USE

Once fully charged, turn on by holding the MODE button for 3 seconds. Use MODE to switch on and off and to switch through the different skincare settings. When using CLEANSE mode use the ionic ring attachment to secure cleansing cotton pads. Remove when not in CLEANSE mode.

3. LEVEL

In the modes of "CLEAN", "IMPORT", "EYE CARE", "EMS UP" and "RF LED", you can use the LEVEL button to set your desired strength. There are a total of 5 levels to use.

Level 1	Level 2	Level 3	Level 4	Level 5
				

4. AUTO SHUTDOWN AND SWITCH OFF

After 1 minute of inactivity, the device will automatically turn off. You will be notified by a beeping noise.

SAFETY PRECAUTIONS

To prevent damage to persons and property, the following describes the precautions.

The following symbols are used to classify and describe the hazard levels and property damage improper use causes.

DANGER Denotes severe and likely risk of death or serious injury.

WARNING May denote risk of death or serious injury observed.

The following symbols classify and describe the type of instructions to be observed

Alerts users to a specific operating procedure that must not be performed.

Alerts users to a specific operating procedure that must be performed.

WARNING: Do not use this appliance near bathtubs showers, basins or other vessels containing water

This marking indicates that the product should not be disposed of with other household wastes. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources.

PREVENTING ACCIDENTS

Do not use this appliance near electronic medical devices, such as the following:

Failure to observe this precaution may cause incorrect operation of electronic medical devices, resulting in accidents or causing patients to feel unwell.

- Pacemakers and other electronic medical devices implanted in the body
 - Heart-lung machines and other life-supporting electronic medical devices
 - Wearable electronic medical devices, such as electrocardiographs and intravenous drip equipment
 - Keep away from flammable aerosol products being used or where oxygen is being administered
 - Keep the cord away from water and heat
 - DO NOT operate this device if the cord or the plug has been damaged or submerged in water
 - DO NOT put into the dishwasher or submerge in water or liquid
- 
- 

CLOUD NINE® | BEAUTY

SKIN TESTING FOR LIGHT SENSITIVITY

Certain medicines or products can cause light-sensitive skin, please consult a doctor prior to use. It is recommended to test the CLOUD NINE Rejuvenate on your arm or leg.

“RF LED” and “COOL” modes refer to Red or Blue Light.

TO PERFORM A SKIN TEST:

1. Turn on the CLOUD NINE Device.
2. Test the device on a small area of your arm or leg for 5 minutes.
3. Wait 24 hours. You should NOT experience any adverse reaction during the next 24 hours. If you have no reaction, you may continue with the treatment.

If you experience any pain or skin redness in the treatment area, please stop using the device and if necessary, consult your doctor.

BE FIERCELY YOU

To prevent skin and health problems the following persons should not use this product:

Failure to observe these precautions may in rare circumstances result in unwanted side effects.

- Persons undergoing medical treatment
- Persons with heart diseases
- Persons with facial neuropathy
- Persons with diabetes
- Persons with blood abnormalities such as hemorrhagic disease and hemophilia
- Persons with acute illnesses
- Persons with febrile diseases
- Persons with infectious diseases
- Persons with high blood pressure
- Persons undergoing dental treatment
- Persons with implanted metal, silicone, plastic, etc. on their faces due to fractures/plastic surgery
- Persons with capillary dilation due to prolonged use of steroid hormone drugs or liver dysfunction
- Do not use if you have epilepsy
- Do not use if using cancer fighting drugs or drugs increasing photo sensitivity
- Infants, growing children and persons who cannot express their intentions by themselves
- Persons with anterior or ophthalmic sequelae due to traffic accidents, etc.
- Persons who have undergone eye surgery (If you have undergone LASIK surgery, please consult with your family doctor)
- Persons with allergies
- Persons with sensitive skin
- Persons with skin diseases or atopic dermatitis
- Persons during menstruation/pregnancy or within one month after giving birth
- Persons who are inebriated
- Persons taking medicine
- Persons who cannot sense temperatures

Failure to observe this precaution may cause burns, electric shocks and injuries.

Do not use this appliance in the following body areas.

Failure to observe this precaution may cause problems with the skin and body, and impaired vision.

- Parts with scratches, eczema, swelling, etc.
- Parts with abnormalities such as dermatitis and excessive sunburn
- Parts with inflammation due to acne and pimples
- Parts where hair removal treatment or cosmetic surgery was performed today
- Parts with itching, hot flashes or pathologic spots due to physical stimulus
- Eyeballs, lips and Adam's apple
- On the eyelid
- Do not use over the thyroid gland (base of the neck)

If you feel pain or abnormalities on your skin or other parts of your body, immediately stop using this product.

Failure to observe these recommendations in rare occasions may result in unwanted side effects.

CARE AND MAINTENANCE

- Always unplug the CLOUD NINE Device from the power supply after each charge.
- Clean the device with a soft dry cloth only. Store it in a dry, safe place.
- The CLOUD NINE device will be fully charged after 3 hours and have a battery life of 45-60 minutes.

TECHNICAL FAULTS

Device won't turn on	<ul style="list-style-type: none">• Press [MODE] button for 3S• Charge the battery
Unable to charge	<ul style="list-style-type: none">• Remove transparent plastic adhesive films before charging.• Please keep both terminals clean and dry.• Firmly plug the AC adapter to the power supply• Ensure the device is secured firmly to the charging base
Not maintaining battery level	<ul style="list-style-type: none">• Make sure the device is fully charged prior to use.

WARRANTY

This CLOUD NINE Rejuvenate Device comes with a 12-month warranty. In the unlikely event that your product develops a fault, our warranty covers the replacement of your product for free within the guarantee period.

Customer Care team or visit
www.cloudninehair.com/pages/returns

BE FIERCELY YOU

CLOUD NINE® | BEAUTY

OUR GIFT TO YOU

Use this exclusive code to get 10% off your next
CLOUD NINE styling tool or beauty device.

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CLOUD NINE International Ltd,
CLOUD NINE ANZ Pty Ltd,