

Regrow^{MD}

MEDICAL DEVICE FOR HAIR GROWTH



LASER272

HAIR GROWTH LASER LIGHT DEVICE

..... **BIOLIGHT COMFORT**
DESIGN

USER MANUAL

TABLE OF CONTENTS

What you should know before you start

Introduction	1
Who is the best candidate? (Classifications male/female).....	2
Who is the best candidate? (Fitzpatrick Skin Types).....	3
Warnings and Cautions	4
Risks/Symbols/Testing for Sensitivity.....	5
How to Monitor Device Activity	6
Device Failure Mode	6
Expected Life of your Device	6
Summary of Clinical Study in Men.....	7
Summary of Clinical Study in Women	8

RegrowMD

Parts.....	9
Getting Started.....	10
Charging.....	10
Proper Use.....	11
Treatment Instructions	11
What is Included.....	12
Troubleshooting	13

Other Useful Information

Safety Label.....	14
Laser Parameters	14
Use Conditions	14

Power Supply Specifications.....	14
Maintenance/Cleaning/Disposal Instructions.....	15
Disposal Instructions.....	15
Electrical Safety.....	16
Electromagnetic Immunity Guidance	17
Electromagnetic Distance Recommendations.....	20
Warranty	21

Corporate Office – Customer Service – Warranty Repairs

Lexington Intl., LLC
 1040 Holland Drive
 Boca Raton, Florida 33487 USA
 Direct line +1-561-417-0200
 Email: info@regrowmd.com
www.regrowmd.com

United States 1-866-527-3726

United Kingdom 0808-134-9868

Australia 1-800-937-358

©Copyright 2020 Lexington Intl., LLC. All rights reserved.

Introduction of RegrowMD


Your RegrowMD Laser 272 (RegrowMD) is a home-use laser device that utilizes Clinically Proven Technology to treat male and female pattern hair loss, also known as hereditary hair loss in men and women. RegrowMD is an over-the-counter (OTC) device that can be sold directly to the public.

Indications

The RegrowMD Laser 272 is indicated to promote hair growth in females with Androgenetic alopecia who have Ludwig (Savin) Classifications I through II, and males who have Norwood Hamilton Classifications of IIa to V and for both with Fitzpatrick Skin Types I to IV.

How RegrowMD works

RegrowMD works by providing laser energy to stimulate hair follicles.



Treatment with RegrowMD

RegrowMD should be used 3 times per week, spaced out every other day. For example, Monday, Wednesday, Friday. The RegrowMD should not be used more than the recommended treatment time per day. Clinical studies demonstrate that most users see benefits in as little as 16 weeks.

RegrowMD Benefits

A clinical study using this technology was conducted in men and women over a six month time period. The study showed that 90% of the men and women had an increase in the number of hairs.

See pages 7-8 for detailed results of the study.


It is important that you read and understand the content of this user manual before using your RegrowMD.

WHO IS THE BEST CANDIDATE FOR THE REGROWMD?


How do you know if the RegrowMD is right for you?

RegrowMD can be used by men and women with thinning hair (including color treated hair) or pattern baldness caused by a hereditary condition. Doctors use a system known as the Norwood Hamilton Classification (men) and the Ludwig-Savin Classification (women) to describe the degree of hair loss. Below are the pictures of the scales. The shaded areas show the type of hair loss that can be treated with RegrowMD.

Norwood Hamilton
Classification - Males



Ludwig-Savin
Classification - Females



RegrowMD has not been clinically tested for hair loss outside of the shaded areas.

WHO IS THE BEST CANDIDATE FOR THE REGROWMD?

RegrowMD technology was studied in very white to brown skin tones (Types 1,2,3,4) as described in the Fitzpatrick scale below. This does not mean it will not work on the darkest skin tones. Due to limitations of digital photography and hair counting software we were unable to test the darkest skin tones.

Fitzpatrick Skin Type Classification Guide: A way to describe your skin type based on your complexion and what happens to your skin in sunlight

Type 1	Clinically Tested to be effective for these skin types	Very white or freckled skin, always burns with sun exposure (very fair; often in people with red or blond hair and blue eyes)
Type 2		White skin, usually burns with sun exposure (fair; often in people with red or blond hair and blue, green, or hazel eyes)
Type 3		White or olive skin tone, sometimes burns with sun exposure (fair; seen in people with any hair or eye color)
Type 4		Brown skin, rarely burns with sun exposure (common in people of Mediterranean descent)
Type 5	Not Clinically Tested for these skin types	Dark brown skin, very rarely burns with sun exposure (common in people of Middle-Eastern descent)
Type 6		Black skin, rarely burns with sun exposure

Adapted from Fitzpatrick's Dermatology in General Medicine. McGraw-Hill Professional; 5th edition.

WARNINGS AND CAUTIONS

WARNINGS

- ⚠ Do not stare directly at the red laser light or its reflection in a mirror because it can temporarily irritate your eyes. If you stare at the light for too long, it could harm your eyes. Also, never use a magnifying glass when using RegrowMD because it can cause temporary irritation to your eyes.
- ⚠ Keep the RegrowMD out of the reach of children and pets.
- ⚠ Do not drop your RegrowMD in water because you may get an electrical shock. If your RegrowMD gets damaged by water, contact customer service toll-free at the numbers provided below.
- ⚠ If the cord is damaged on your RegrowMD, do not use the device because you could get an electrical shock. If the cord gets damaged, contact customer service toll-free at the numbers provided below.
- ⚠ Do not use your RegrowMD near any heated surfaces because this could cause the device not to work properly and may cause an electrical shock. Contact customer service toll-free at the numbers provided below.
- ⚠ Laser light - avoid direct eye exposure.
- ⚠ This is a harmless low-level ESD (Electrostatic Discharge) sensitive device. Please do not touch the charging contacts or USB-port. (see page 9). Touching the contacts or USB-port is not of the essential performance and will not cause a basic safety issue. If touching the contacts or USB-port causes the device to malfunction, reset your device by powering off and back on.
- ⚠ No modification of this equipment is allowed.
- ⚠ This equipment has no user-serviceable parts. Do not attempt to modify or repair this equipment.
- ⚠ Do not attempt to modify or service while the equipment is in use.
- ⚠ If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use. Contact RegrowMD customer service toll free at the numbers provided below.

CAUTIONS

- ⚠ If you are taking any medications or products that make you sensitive to light, you should test your RegrowMD to make sure you do not have a reaction. Some of the things that can cause increased sensitivity to your RegrowMD could be medicines for colds, allergies, pain medicine and drugs to treat infections. To test your RegrowMD for sensitivity, turn the device on, hold the device on one spot, shining the light on your forearm. After two minutes, turn the device off and set aside. Look at the spot where you held the laser light over the next five minutes. If you see any reaction such as redness or rash, do not use your RegrowMD and contact our customer service toll-free at the numbers provided below for information on how to return your product and get a refund.
- ⚠ Only use your RegrowMD as directed. Do not try to change the way your RegrowMD operates. Altering controls, adjustments, performance or procedures other than those specified herein may result in hazardous radiation (laser) exposure. Contact customer service toll-free at the numbers provided below.
- ⚠ Never use acetone or other solvents on any part of your RegrowMD.

United States 1-866-527-3726















United Kingdom 0808-134-9868

Australia 1-800-937-358

Risks

RegrowMD technology has been clinically tested on both men and women. In the clinical study for men, 5% of the men said they had itching or tingling from using this technology. None of the women in the clinical study reported any of these side effects. If you have itching, tingling or any other side effect that bothers you, stop using your device and contact customer service toll-free at the numbers provided below.

Symbols

-  European Conformity
-  WEEE
-  Warning
-  Caution
-  Laser
-  Laser warning
-  Laser Aperture ↑ Laser aperture indication
-  Laser warning label
-  Power On/Off
-  Refer to instruction manual
-  RF interference
-  Type BF Applied Part
-  Serial Number
-  IP22 Degree of moisture resistance

Testing RegrowMD for Sensitivity

Before beginning treatment, be sure to test your RegrowMD for sensitivity. You only need to do this once to test your skin for any reaction to light therapy. Make sure your device is fully charged (see page 10). Turn on your device by pressing the power button on the battery pack. Place your forearm with your palm facing upward under the device. Allow your arm to stay in the light for 2 minutes. Turn your device off and set it aside.

If you don't experience any reaction, proceed with the treatment instructions on page 11. If you see any reaction such as redness or rash, do not use your device. Contact customer service toll-free at the numbers provided below or your place of purchase.

United States 1-866-527-3726

United Kingdom 0808-134-9868

Australia 1-800-937-358

How to Monitor Device Activity

Your RegrowMD is carefully inspected and tested during production. When working properly, a red laser light shines from inside the laser device.

DO NOT look directly at the laser light. Shine the laser light onto the palm of your hand to see if the red laser light appears.

If the laser light does not appear, contact customer service toll-free at the numbers provided below.

Device Failure Mode

Your device is designed to last for 5 years of regular use. The laser light does not fade or diminish over time. At the end of the laser's life or if it fails for any other reason, the red laser light will no longer shine out of the device.

Expected Life of Your Device

RegrowMD is a medical device manufactured under strict quality controls and validated through GMP (Good Manufacturing Practices) and quality standards as outlined by ISO 13485. RegrowMD contains high quality laser diodes.

Based on that rating and the durability of the components being used, the expected life of your RegrowMD is 5 years of regular use. Please take time to read this entire user manual before operating your device.

United States 1-866-527-3726

United Kingdom 0808-134-9868

Australia 1-800-937-358

SUMMARY OF CLINICAL STUDY IN MEN

Clinical studies on RegrowMD technology were conducted at multiple locations in the United States. Participants in the study were males diagnosed with Androgenetic Alopecia, which includes hereditary hair loss, and had light to brown skin tones. Subjects were divided into two groups, active group and a control group. The control device looked and sounded like the active device but did not include a laser light. The Clinical studies adhered to all GCP (Good Clinical Practice) guidelines, were approved by an IRB (Institutional Review Board) and listed on www.clinicaltrials.gov.

During the 6 month clinical study, subjects treated their hair once per day, 3 times per week on non-consecutive days. The number of thick, normal healthy hairs in the target zone were counted at Week 1, Week 16 and Week 26.

Results at the 26 week visit showed that over 90% of the men saw an increase in hair counts.*

No subjects in the study experienced any serious adverse events. In fact, the number and types of adverse events were similar in the active and control groups.


Hair Count Change

In the clinical study, after 26 weeks of treatment with the active device versus the control device, the following increases in hair count were found.

*Based on a minimum of 5 new hairs per square centimeter/32 new hairs per square inch being observed at the follow up visit



Mean Terminal Hair Count Changes
From Baseline in MALES
26 Weeks, Observed Hair Count



SUMMARY OF CLINICAL STUDY IN WOMEN

A clinical study on RegrowMD technology was conducted at multiple locations in the United States. All participants in the study were females, diagnosed with Androgenetic Alopecia, which includes hereditary hair loss, and had light to medium skin tones. Subjects were divided into two groups. One group used the active device and the other group used a control device. The control device looked and sounded like the active device, but did not include a laser light. The Clinical studies adhered to all GCP (Good Clinical Practice) guidelines, were approved by an IRB (Institutional Review Board) and listed on www.clinicaltrials.gov.

During the 6 month study, subjects treated their hair once per day, 3 times per week, on non-consecutive days. The number of thick, normal, healthy hairs in the target area were counted at Week 1, Week 16 and Week 26.

Results at 26 weeks showed that over 90% of the women saw increases in hair counts.*


No subjects in the study experienced any serious adverse events. In fact, the number and types of adverse events were similar in the active and control groups.

Hair Count Change

In the clinical study, after 26 weeks of treatment with the active and a control device, the following increases in hair count were found.

*Based on a minimum of 5 new hairs per square centimeter/32 new hairs per square inch being observed at the follow up visit

Mean Terminal Hair Count Changes
From Baseline in FEMALES
26 Weeks, Observed Hair Count



PARTS

- A LASER DEVICE
- B BIOLIGHT COMFORT INSERT
- C BATTERY PACK
- D INDICATOR LIGHTS
- E ON/OFF BUTTON
- F USB CHARGING CABLE
- G POWER SUPPLY
- H USB DEVICE CABLE
- I SPORTS CAP
- J SAFETY LABEL (on the left side of the laser device)
- K SERIAL NUMBER (on the left side of the laser device)
- L BELT CLIP FOR BATTERY PACK



Your RegrowMD comes with a rechargeable battery pack. Read this user manual completely before you start to use your device. The instructions in this user manual are intended for individuals with a minimum of an eighth grade education. RegrowMD is safe for home use. The patient is the intended operator of this equipment.



Use your RegrowMD for 30 minutes per treatment, three days per week. The treatments should be spaced every other day, for example Monday, Wednesday, Friday.

Most users see benefits in as little as 16 weeks (based on clinical studies).

When you receive your device, take all the parts out of the case. Make sure your device is clean and dry before beginning the treatment. You must charge the battery pack for 6-8 hours before initial use.


Charging the Battery Pack

Charge your RegrowMD by plugging the USB charging cable into the socket on the battery pack labeled 'DC-IN' (Figure A).

Plug the other end into an electrical socket.

The indicator lights on the battery pack will flash during the charging process. Four steady blue lights indicate your device is fully charged.



When your RegrowMD is turned on, the indicator lights will illuminate. As the battery power diminishes, the lights on the battery pack will disappear in sequence to indicate the status of the battery. When this happens, charge your laser device as previously described.



Proper Use

Before beginning treatment, test your RegrowMD for sensitivity (see page 5 of the user manual). Your RegrowMD should be used for 30 minutes, three times per week. The treatments should be spaced out every other day. For example Monday, Wednesday, Friday.

Before starting your treatment, make sure the battery is charged per the instructions on Page 10. Remove the USB charging cable from the battery pack. Connect the USB device cable to the twist lock connector on your RegrowMD device (Figure B). The cord connector will only fit when inserted the right way. Align the indicator arrows. Do not force the connector together, once aligned, the connectors will slide together easily. Once inserted, twist the retaining ring until the two cords are locked together. Connect the other end of the USB device cable to the battery pack socket labeled 'DC-OUT' (Figure C). You are now ready to start your treatment.



Your RegrowMD can be used by itself or inserted into the sports cap (provided).


Treatment Instructions

To start your treatment, place your RegrowMD on your head in a comfortable position (Figure D) and press the ON/OFF button on the battery pack (Figure E).

The blue indicator lights on the battery pack show that the device is in active mode and the red laser lights are on (Figure E).

When you press the ON/OFF button, the device will beep twice. During usage, your device will beep once after 10 minutes and 20 minutes to indicate the elapsed treatment time. After 30 minutes, the laser device will beep twice and turn off.

Disconnect the battery pack after the treatment is completed.



WHAT IS INCLUDED

					
Laser 272	Battery Pack and Belt Clip for Battery Pack	USB Device Cable	Power Supply	USB Charging Cable	User Manual

Notice: Entire RegrowMD assembly is considered as a Type BF applied part.

TROUBLESHOOTING

ANY ATTEMPT TO SELF-SERVICE YOUR DEVICE, UNLESS INSTRUCTED BY REGROWMD CUSTOMER SERVICE, WILL VOID YOUR WARRANTY.

RegrowMD has advanced electronics to ensure that the device is working like it is supposed to and that the proper amount of laser energy is being delivered. If for some reason you are experiencing difficulties with your device, please try one of the tips below:

Your laser device will not turn on

- Make sure your RegrowMD is fully charged
- Make sure the laser device and the battery pack are properly connected

Your laser device is turning off during use

- Make sure your RegrowMD is fully charged
- Make sure the laser device and the battery pack are properly connected

Your laser device is not lighting up, but the battery pack is fully charged

- Contact customer service to determine if we need to send you a replacement unit

Your laser device does not show the blue indicators

- Contact customer service to request a replacement battery pack

If you have trouble understanding any of these instructions, contact customer service toll-free at the numbers provided below.

Customer Service – Warranty Repairs

Lexington Intl., LLC
1040 Holland Drive
Boca Raton, Florida 33487 USA

United States 1-866-527-3726


United Kingdom 0808-134-9868

Australia 1-800-937-358

SAFETY LABEL | LASER PARAMETERS | USE CONDITIONS

Wavelength:	655nm +/- 10nm visible red light
Laser Power for classification:	<5mW CW
Divergence:	57 mrad
Number of lasers:	272
This Laser Product is designated as Class 3R during all procedures of operation.	
Note: RegrowMD laser devices have no user-serviceable parts.	
Safety label location:	on the left side of the laser device
MTBF: mean time between failures:	10,000 hours at 25°C ambient temperature
Operating temperature/altitude/atmospheric pressure/humidity:	5°C to 35°C/0 to 2,000 m/80kPa to 106kPa/ 0% to 90%, relative humidity
Storage & transportation temperature/humidity/atmospheric pressure:	-20°C to 60° C / 0% to 90%, relative humidity/ 50kPa to 106kPa
Ingress of water or particulate matter:	IP22
Battery charging cycles:	300 cycles
ROHS: restriction of hazardous substances:	2011/65/EU
Device size and weight:	22 x 18 x 9 cm/440 grams

Warning label

Serial number
(on the left side of
the laser device)

Power Supply Specifications

Output = 5V DC-1.0A

The power supply for the RegrowMD can be used in the United States and internationally. This power supply will adapt to all power supplies from 100 volt/60 cycle (US) to 240 volt/50 cycle. If outlets in your country have different plugs than the power supply that comes with the RegrowMD, use an international plug adapter (not supplied with the device).

Maintenance

Your RegrowMD is designed to last for 5 years. It does not require any maintenance other than periodic cleaning.

Cleaning your RegrowMD

- Keep your RegrowMD clean.
- Turn off the RegrowMD by pressing the ON/OFF button until the light turns off.
- Unplug the USB device cable from the laser device and battery pack.
- Clean your device with a soft, slightly moistened cloth before and after each treatment.

Disposal Instructions


















This symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

Electrical Safety

Basic safety precautions should be observed when using this, or any other electrical device, especially when children are present.

 **WARNING:** Keep away from water. Contact with water may result in an electrical shock.

To reduce the risk of hazard by electric shock:

-  Do not use while bathing.
-  Do not place or store your RegrowMD where it can fall or be pulled into a tub or sink.
-  Do not expose your RegrowMD to running water, or drop into water or any other liquid.
-  If your RegrowMD falls into water, unplug it immediately from the power outlet.
If using your RegrowMD near water, the use of an electrical outlet protected by a Ground Fault Circuit Interrupter is recommended.
-  Do not dispose of your RegrowMD in a fire or store in high temperatures greater than 140°F (60°C).
-  Do not turn on or use your RegrowMD while it is charging.
-  Power supply cannot be separated during use. Do not position the equipment in a manner that will make it difficult to plug in or remove the power cord.
-  If you have any problems with this device, such as setting up, maintaining or using, please contact customer service.
-  Please report to customer service if any unexpected operation or events occur.
-  After drop/shock that may cause changes in the performance, please contact customer service.
-  Be careful, do not put the power cord around your neck, it may cause strangulation.
-  Keep the device out of the reach of children/pets to avoid inhalation or swallowing of small parts.
-  Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
-  The device should be used only with the accessories recommended for use by the manufacturer.
-  When not in use, store the device with the adapter in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight.

Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

The RegrowMD Laser 272 is intended for use in the electromagnetic environment specified below. The customer or the user of the RegrowMD Laser 272 should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The RegrowMD Laser 272 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The RegrowMD Laser 272 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complies	

ELECTROMAGNETIC IMMUNITY GUIDANCE

Guidance and manufacturer's declaration-electromagnetic immunity

The RegrowMD Laser 272 is intended for use in the electromagnetic environment specified below. The customer or the user of the RegrowMD Laser 272 should assure that it is used in such an environment.


Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	± 2kV for power supply lines ± 1kV for input/output lines	± 2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV differential mode ± 2kV common mode	± 1kV differential mode ± 2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle g) At 0°, 45°, 90°, 135°,180°, 225°, 270°, and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycle	0% UT; 0,5 cycle g) At 0°, 45°, 90°, 135°,180°, 225°, 270°, and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the RegrowMD Laser 272 requires continued operation during power mains interruptions, it is recommended that the RegrowMD Laser 272 be powered from an uninterruptible power supply or a battery.
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

ELECTROMAGNETIC IMMUNITY GUIDANCE

Guidance and manufacturer's declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

The RegrowMD Laser 272 is intended for use in the electromagnetic environment specified below. The customer or the user of the RegrowMD Laser 272 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 KHz to 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 10 V/m 80 MHz to 2.7 GHz	3 V 6 V 10 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the RegrowMD Laser 272 including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{12}{V_2} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 2,7 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,765 MHz to 6,795 MHz, 13,553 MHz to 13,567MHz; 26,957 MHz to 27,283 MHz; and 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,0 MHz to 21,4 MHz.
- b Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephone and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the RegrowMD Laser 272 is used exceeds the applicable RF compliance level above, the RegrowMD Laser 272 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the RegrowMD Laser 272.
- c Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

ELECTROMAGNETIC DISTANCE RECOMMENDATIONS

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

The RegrowMD Laser 272 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the RegrowMD Laser 272 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the RegrowMD Laser 272 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz outside ISM and amateur radio bands $d = [\frac{3.5}{V_1}] \sqrt{P}$	150 kHz to 80 MHz outside ISM and amateur radio bands $d = [\frac{12}{V_2}] \sqrt{P}$	80 MHz to 800 MHz $d = [\frac{3.5}{E_1}] \sqrt{P}$	800 MHz to 2.7 GHz $d = [\frac{7}{E_1}] \sqrt{P}$
0.01	0.12	0.20	0.035	0.07
0.1	0.38	0.63	0.11	0.22
1	1.2	2.0	0.35	0.70
10	3.8	6.32	1.10	2.21
100	12	20.00	35	70

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

WARRANTY

The Regrow Laser Device has a 2 year warranty against defects in materials or workmanship from Lexington International, LLC. Items that are damaged or abused are not covered under the warranty and will be charged at their regular repair rates. Any attempt to self-service your device, unless instructed by LEXINGTON INTERNATIONAL CUSTOMER SERVICE, will void your warranty. Please refer to the user manual for complete warranty terms.

To validate your warranty, register online at
www.regrow.com/warranty

If you do not have a computer you can register by phone at the following numbers:

Direct line +1 561-417-0200

United States and Canada 1-866-527-3726

United Kingdom 0808-134-9868

Australia 1-800-937-358

Monday - Friday 9AM - 5PM ET

WARRANTY

ANY ATTEMPT TO SELF-SERVICE YOUR DEVICE, UNLESS INSTRUCTED BY REGROWMD CUSTOMER SERVICE, WILL VOID YOUR WARRANTY.

THE REGROWMD IS WARRANTED BY LEXINGTON INTERNATIONAL, LLC AGAINST DEFECTS IN MATERIALS OR WORKMANSHIP FOR 24 MONTHS FROM DATE OF PURCHASE. IF THIS PRODUCT IS DETERMINED TO BE DEFECTIVE, LEXINGTON WILL REPAIR OR REPLACE IT, AT ITS OPTION, AT NO CHARGE DURING THE WARRANTY PERIOD.

AFTER THE WARRANTY PERIOD, THE USER MUST PAY FOR ALL PARTS AND LABOR CHARGES. IF THE PRODUCT IS REPAIRED OR REPLACED, EITHER WITHIN OR OUTSIDE THE WARRANTY PERIOD, LEXINGTON RESERVES THE RIGHT TO PROVIDE EQUIVALENT PARTS IN DIFFERENT COLORS TO YOUR ORIGINAL PRODUCT.

THIS WARRANTY DOES NOT COVER DAMAGES CAUSED BY ACTS OF GOD, MISUSE, NEGLIGENCE, ACCIDENT, MODIFICATION OF, OR TO ANY PART OF THE PRODUCT. THIS WARRANTY DOES NOT COVER DAMAGE DUE TO IMPROPER OPERATION OR MAINTENANCE, CONNECTION TO IMPROPER VOLTAGE SUPPLY, OR ATTEMPTED REPAIR BY ANYONE OTHER THAN LEXINGTON. THIS WARRANTY IS VOID IF PRODUCT IS TAMPERED WITH.

TO OBTAIN WARRANTY SERVICE, YOU MUST CALL REGROWMD CUSTOMER SERVICE TO REQUEST AN AUTHORIZATION NUMBER. YOU MUST DELIVER THE PRODUCT FREIGHT PREPAID, IN EITHER ITS ORIGINAL PACKAGING OR PACKAGING AFFORDING AN EQUAL DEGREE OF PROTECTION, TOGETHER WITH YOUR NAME AND ADDRESS, TO:

Lexington International, LLC, 1040 Holland Drive, Boca Raton, FL 33487, USA

LIMITATION OF LIABILITY; INDEMNITY: EXCEPT AS EXPRESSLY SET FORTH ABOVE, LEXINGTON INTERNATIONAL, LLC SHALL NOT BE LIABLE FOR ANY SPECIAL, INDIRECT, OR CONSEQUENTIAL DAMAGES THAT ARE DIRECTLY OR INDIRECTLY RELATED TO THE USE OF, OR THE INABILITY TO USE, THE REGROWMD, EVEN IF LEXINGTON INTERNATIONAL, LLC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES WHETHER IN AN ACTION OF CONTRACT, NEGLIGENCE OR OTHER TORTUOUS ACTION. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

GOVERNING LAW: THE USAGE OF THE REGROWMD AND THE INTERPRETATION OF THE TERMS OF THIS AGREEMENT SHALL BE GOVERNED BY AND CONSTRUED IN ACCORDANCE WITH THE LAWS OF THE STATE OF FLORIDA APPLICABLE TO AGREEMENTS NEGOTIATED, EXECUTED AND PERFORMED WHOLLY AND SOLELY WITHIN THE STATE OF FLORIDA. THE PARTIES AGREE THAT THE STATE OF FLORIDA'S CONFLICTS OF LAWS RULES SHALL NOT APPLY TO THIS AGREEMENT OR ANY TRANSACTION PERFORMED BY OR UNDER THIS AGREEMENT.

JURISDICTION AND VENUE: THE PARTIES IRREVOCABLY SUBMIT AND CONSENT TO THE SOLE AND EXCLUSIVE JURISDICTION AND VENUE OF THE FLORIDA STATE COURTS IN AND FOR PALM BEACH COUNTY, FLORIDA, U.S.A. AND THE FEDERAL COURTS IN AND FOR THE SOUTHERN DISTRICT OF FLORIDA, U.S.A. THE PARTIES AGREE NOT TO RAISE THE DEFENSE OF FORUM NON CONVENIENS.

AUTHORIZED REPRESENTATIVE IN EUROPE (REGULATORY AFFAIRS ONLY): EMERGO EUROPE, MOLENSTRAAT 15, 2513 BH THE HAGUE, THE NETHERLANDS. TEL: (31) (0) 70 345-8570 FAX: (31) (0) 70-345-7299

To validate your warranty please register online at www.regrowm.com/warranty

Corporate Office – Customer Service – Warranty Repairs

Lexington Intl., LLC
1040 Holland Drive
Boca Raton, Florida 33487 USA

Direct line +1 561-417-0200

United States and Canada 1-866-527-3726

United Kingdom 0808-134-9868

Australia 1-800-937-358

Monday - Friday 9 AM - 5 PM ET

Authorized Representative in Europe (Regulatory Affairs Only):



Tel: +31 (0) 70 345-8570 Fax: +31 (0) 70-345-7299

Australian Sponsor: Emergo Australia
Level 20, Tower II Darling Park, 201 Sussex Street
Sydney, NSW 2000 Australia

Regrow^{MD}

MEDICAL DEVICE FOR HAIR GROWTH

