Masderm

RF Face Device Plus
DISPOSITIVO RADIOFRECUENCIA
FACIAL PLUS

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2. FUNCTIONALITY OF THE DEVICE

Our RF FACE DEVICE PLUS device (facial radiofrequency device plus) has the main objective of improving the general appearance of facial skin as it will help you reduce sagging caused by the passing of the years, wrinkles, spots and also , will improve the luminosity, uniformity and softness of the skin, also helping to delay the signs of aging.

In addition, it combines several functions, it is a 7-in-1 device. It has specific functions to treat the eye area, wrinkles, sagging, jaw, among other technologies that will help your skin regain a younger appearance. With RF FACE DEVICE PLUS you do not need another type of specific device since it covers all your needs.

2. PARTS OF THE DEVICE

3. PACKAGING CONTENTS

When you open the packaging box you will have found the following:

1. The RF FACE DEVICE PLUS device.
2. This instruction manual.
3. A USB cable for device charging with USB A connector on one end and USB C on the other.

4. HOW TO PLUG IN THE DEVICE

Before starting to use our RF FACE DEVICE PLUS you need to charge it for approximately 3 hours to ensure a full charge.

Plug the Type C cable into the device and the USB A portion into a mains adapter.

IMPORTANT: When you are going to use the device, you must unplug it from the power. You should never use the device while it is charging.

5. HOW TO START IT AND SHUT IT OFF

1. To turn on the device, press and hold the on/off button for a few seconds.
2. The device turns on automatically using the CLEAN function at the lowest level the first time, but it has memory, meaning that if the level is changed, the next time you turn it on it will continue at that same level.
3. Each performance lasts 4 minutes. All functions warn with a beep halfway through the time. At the end of the function, the device moves on to the next function automatically, you will see that the next function starts flashing and you just have to press the on/off button to activate it. If you do not press it, in 1 minute the device will turn off automatically.
4. If you want to change the function manually, press the same on/off button.
5. To change the power level, press the levels button. Our device has 3 levels but you should always start with the lowest level. Except for the COOL function which has no levels.
6. All functions except the COOL function vibrate but you can stop the vibration by pressing and holding the level button. Likewise, you can turn the vibration back on.

6. FUNCTIONS OF THE RF FACE DEVICE PLUS

Our RF FACE DEVICE PLUS device uses different technologies such as radio frequency, electrostimulation, long frequency vibration, pulse vibration, positive and negative ion and LED light.

CLEAN: Cleans residual dirt from the pores, improves acne reduction and improves the skin's own luminosity thanks to cellular oxygenation.

IMPORT: Improves the appearance of spots and wrinkles by promoting and accelerating the absorption of nutrients from the skin, improving the effectiveness of our products and the overall appearance of skin care.

EMS: Promotes microcirculation and improves skin firmness thanks to muscle toning.
RF: Promotes the formation of collagen and elastin among others. Firms and improves fine lines of expression.

EYE: Helps reduce and blur wrinkles around the eyes, bags, dark circles, edema and also relaxes the area.

JAW: Improve the jawline to relax and firm the contour.

COOL: Soothes and softens the skin, reduces enlarged pores and reduces inflammation.

The CLEAN, IMPORT, RF and EYE functions use RF technology.

The EMS, JAW and COOL functions do not use RF technology.

IMPORTANT: The device must always be in continuous movement and with enough product to slide the head correctly. If necessary, apply more product.
Do not stop the device at any time in the same area.

7. APPLICATION TIPS

Our device has a total of 7 functions, you can use the device for a maximum of 20 total minutes per session (1 session per day maximum) divided into: 2 functions with RF + 2 functions without RF + COOL function.

The COOL function is the only one that can be used every day.

The rest of the functions can be used every other night, that is, about 3 times a week.

Example:

8 minutes of RF (4 +4) and 8 minutes others functions (4 JAW + 4 EMS) + 4 COOL = 20 MIN

IMPORTANT:
(4min per Function / Day)

	L	M	X	J	V	S	D
WEEK 1	RF + EYE + EMS + JAW + COOL	COOL	RF + EYE + EMS + JAW + COOL	COOL	RF + EYE + EMS + JAW + COOL	COOL	COOL
WEEK 2	CLEAN + IMPORT + EMS + JAW + COOL	COOL	CLEAN + IMPORT + EMS + JAW + COOL	COOL	CLEAN + IMPORT + EMS + JAW + COOL	COOL	COOL
WEEK 3	RF + EYE + EMS + JAW + COOL	COOL	RF + EYE + EMS + JAW + COOL	COOL	RF + EYE + EMS + JAW + COOL	COOL	COOL
WEEK 4	CLEAN + IMPORT + EMS + JAW + COOL	COOL	CLEAN + IMPORT + EMS + JAW + COOL	COOL	CLEAN + IMPORT + EMS + JAW + COOL	COOL	COOL

This is a practical example, but the device is very versatile. Each person can use it according to their needs as long as they respect the times and uses of each function)

8. HOW AND WHEN TO USE OUR DEVICE

We recommend using our device 2 or 3 times a week, alternately and always at night, except for the COOL function which you can use every day. Please respect these instructions for use as they are recommended for this type of device and validated by skin care experts.

Step 1: Press and hold the power button for a few seconds to turn on the device. Choose the function you prefer or let them happen automatically.

Below we indicate the steps you must follow to apply our RF FACE DEVICE PLUS correctly:

CLEANING function

If you have makeup, we recommend that you apply a cotton pad to the head and hold it with the accessory piece included with the device. In this case:

Place a cotton ball on the head of the device.

Apply the RF Cleanmas makeup remover gel.

Always turn the device on the CLEAN function at the lowest level and pass it over your face in circular movements.

When finished, remove the excess product with a cotton ball or clean towel. No need to clarify.

If you don't have makeup, you can do without cotton and use RF Cleanmas directly on your face.

This function uses Radio Frequency, High Frequency Vibration, Positive Ion and Red LED technology.

IMPORT Function

Apply enough product to cover your entire face and neck. Depending on your needs, we recommend different products (such as our firming facial cream RF Facial Cream) that we detail on the website so that you have updated information.

Slide the device slowly and continuously over the entire face and neck, avoiding the Adam's apple. Ascending and circular movements to stimulate circulation and improve the penetration and effectiveness of our products.

This function uses Radio Frequency, Pulsed Vibration, Negative Ion and red LED technology.

EMS Function

Apply enough product to cover your entire face and neck. Depending on your needs, we recommend different products (such as our firming facial cream RF Facial Cream)

that we detail on the website so that you have updated information.

Slide the device slowly and continuously over the entire face and neck, avoiding the Adam's apple. Ascending movements with special emphasis on the jaw towards the cheekbones and upward movements in front.

This function uses variable frequency EMS technology, long frequency vibration and red LED.

RF Function

Apply enough product to cover your entire face and neck. Depending on your needs, we recommend different products (such as our firming facial cream RF Facial Cream) that we detail on the website so that you have updated information.

Slide the device slowly and continuously over the entire face and neck, avoiding the Adam's apple. Ascending and circular movements with special emphasis on cheekbones, eye contour and zig zag on the forehead.

This function uses Radio Frequency, EMS, Pulsed Vibration and red LED technology.

EYE Function

Apply enough product to cover the entire eye contour area, in our case we recommend our RF Eyemas eye contour.

Slide the device slowly and continuously with semicircular movements from one side to the other, covering the entire dark circle to drain bags and improve the appearance of the dark circles.

This function uses Radio Frequency technology, EMS variable frequency, long frequency vibration and red LED.

JAW Function

Apply enough product to cover the neck and jaw. Depending on your needs, we recommend different products (such as our firming facial cream RF Facial Cream) that we detail on the website so that you have updated information.

Slide the device slowly and continuously across the neck, avoiding the Adam's apple. Ascending movements and reviewing the jaw until reaching the cheekbone.

This function uses EMS technology, long frequency vibration and red LED.

COOL Function

On skin moistened by the products used previously, slide the device slowly and continuously over the entire face and neck. Ascending and circular movements to stimulate circulation and improve the penetration and effectiveness of our products.

If you want to use only the COOL function, we recommend doing it with clean skin. If you want to decongest bags and dark circles you can apply our RF Eyemas eye contour, and if what you are looking for is to improve the general appearance, you can use our Triphasic Serum.

This function uses low temperature technology at 16 degrees.

This function does not vibrate or have levels.

Step 2: At the end of the treatment, turn off the device or, once you have finished using the COOL function, the device turns off automatically after the second beep (4 minutes of use). You can clean the head part with a damp towel but never immerse the device in water.

Step 3: If you have leftover product on your face, take advantage of it by massaging until completely absorbed.

Step 4: You can now continue with your usual routine. We recommend using photoprotection the next morning.

Avoid very irritating products on the same day of treatment with RF FACE DEVICE PLUS, with the exception of the COOL function as it calms and decongests the skin.

Visual Application Example::

1. Place a cotton ball on the head. Apply RF Cleanmas to the skin. Slide the device to remove impurities. No need to rinse with water

2. Apply the RF Facial Cream. Gently slide the RF FACE DEVICE PLUS. Changes function every 4 minutes. Ends the cold function.

9. CONTRAINDICATIONS

1. Pregnant, lactating women and children.
2. Those who carry electronic medical devices.
3. People with allergic constitution and skin diseases (such as ulcers, ringworm, eczema, rashes, burns, scalds, sunburn, abnormal wound healing, etc.).
4. People with diabetes, lupus erythematosus, porphyria or congestive heart disease, hemorrhagic disease, and a history of immunosuppression (including HIV infection or AIDS).
5. Have a pacemaker or internal defibrillator installed, suffer from arrhythmia or any other known serious heart disease.
6. Wearing metal implants or perforated ornaments and presenting any form of suspicious injury.
7. Those receiving radiotherapy or chemotherapy.
8. Any disease caused by thermal stimulation, such as recurrent herpes simplex virus.
9. Do not use it in case of wound, plastic part, mouth, upper eyelid, eyeball, laryngeal node and lymph node.

Please, if you have any questions, before using our RF FACE DEVICE PLUS or if you have several Masderm devices and want to know how to combine them, contact us at hello@masderm.com and our specialists will be happy to assist you.

10. PRECAUTIONS AND MAINTENANCE

1. Before using the beauty device for the first time, please make sure to fully charge the device.
2. You cannot use the device while it is charging.
3. After using it, make sure you clean the heads correctly with a slightly damp cotton or towel for correct use again.
4. Do not store the device within the reach of children and avoid placing it in humid places and/or exposed to high temperatures.
5. Do not wash it with water directly, nor use corrosive products such as alkaline

detergents, abrasive cleaners, alcohol, thinners, gasoline, nail polish remover, corrosive liquids, etc. to clean it.

6. Please keep the RF FACE DEVICE PLUS device in its original box with all its accessories, since in the event of an incident, we cannot guarantee that we will be responsible for it if all the equipment is not complete.

7. Do not repair, disassemble, or modify the device without authorization. Do not use the beauty device for other purposes.

8. In case of device failure or abnormal reaction, please contact us at hello@masderm.com

11. TECHNICAL INFORMATION OF THE DEVICE

Name: RF FACE DEVICE PLUS

Model: JC-F2610

Size: 52x52x175mm

Net weight: 165g

Input: 5V=1A

Power: 8W

Battery: 1050mAh

Battery time: 2 hours

Levels: 3 gears

Execution standards: GB4706.1-2005

GB4706.15-2008 GB/T26572-2011

Produced: China

*Masderm is not responsible for possible collateral damage in the use or maintenance of the device if the instructions in these instructions have not been correctly followed.

These instructions are subject to change for better customer understanding.

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