



# MANUAL

THE CLEOPATRA MASK

# USAGE GUIDE

**STEP 1:** Charge your mask fully by using the provided USB cable. Once fully charged as indicated by a green circular light, you may unplug the mask. (Do not attempt to use while plugged in)

**STEP 2:** Insert provided eye protection in the eye slots.

**STEP 3:** Cleanse your face and remove any makeup.

**STEP 4:** Flip the "ON" switch located at the top of the mask.

**STEP 5:** Power on the mask by tapping the middle of the forehead with your palm.

**STEP 6:** Select desired color by tapping above the right eye (if looking at the mask)

**STEP 7:** Select desired power level by tapping above the left eye (if looking at the mask)

**STEP 8:** Set adjustable strap to fit your face. (Easiest to do when off your face and to test from there)

**STEP 9:** Enjoy for 15 minutes/day

**STEP 10:** Power off device by tapping the middle of the forehead again. Then flip the OFF switch at the top of the mask.

# SETTINGS

Note: It is always easiest to change settings with the palm of your hand instead of your finger tips.

Your mask comes with 7 different colors and 4 different power levels. Additionally you may also select the 8th color option which cycles through all colors.

- 1** USB/On-Off Switch
- 2** Power On
- 3** Switch Color
- 4** Adjust Power Level



Studies have shown these colors may assist with:

**Red** - Anti-aging, fine lines, wrinkles, collagen and elastin production.

**Green** - Anti-inflammation, hyper-pigmentation, dark spots, and oily skin.

**Blue** - Anti-acne and psoriasis.

**Purple** - Combination of red and blue light.

**Yellow** - Balancing of skin texture, redness, and rosacea.

**Cyan** - Reduces swollen capillaries and inflammation.

**White** - Anti-inflammation, boosts skin metabolism and tightening/firming.

Citations:

Wunsch A, Matuschka K. A controlled trial to determine the efficacy of red and near-infrared light treatment in patient satisfaction, reduction of fine lines, wrinkles, skin roughness, and intradermal collagen density increase. *Photomed Laser Surg.* 2014;32(2):93-100. doi:10.1089/pho.2013.3616

Gold MH, Andriessen A, Biron J, Andriessen H. Clinical Efficacy of Self-applied Blue Light Therapy for Mild-to-Moderate Facial Acne. *J Clin Aesthet Dermatol.* 2009;2(3):44-50.

Lubart R, Lipovski A, Nitzan Y, Friedmann H. A possible mechanism for the bactericidal effect of visible light. *Laser Ther.* 2011;20(1):17-22. doi:10.5978/islsm.20.17

Avci P, Gupta A, Sadasivam M, et al. Low-level laser (light) therapy (LLLT) in skin: stimulating, healing, restoring. *Semin Cutan Med Surg.* 2013;32(1):41-52.

Lee MW. Combination 532-nm and 1064-nm lasers for noninvasive skin rejuvenation and toning [published correction appears in *Arch Dermatol.* 2004 May;140(5):625]. *Arch Dermatol.* 2003;139(10):1265-1276. doi:10.1001/archderm.139.10.1265

Rohringer S, Holnthoner W, Chaudary S, et al. The impact of wavelengths of LED light-therapy on endothelial cells. *Sci Rep.* 2017;7(1):10700. Published 2017 Sep 6. doi:10.1038/s41598-017-11061-y

Kleinpenning MM, Otero ME, van Erp PE, Gerritsen MJ, van de Kerkhof PC. Efficacy of blue light vs. red light in the treatment of psoriasis: a double-blind, randomized comparative study. *J Eur Acad Dermatol Venereol.* 2012;26(2):219-225. doi:10.1111/j.1468-3083.2011.04039.x

# SAFETY

## CONTRAINDICATIONS

- DO NOT** use if you are pregnant, planning on becoming pregnant or breastfeeding.
- DO NOT** use if you suffer from lupus erythematosus, photosensitive eczema or albinism.
- DO NOT** use if you suffer from any photosensitive disorder (sensitization to light).
- DO NOT** use if you are taking any medication that can cause photosensitivity.
- DO NOT** use if you suffer from light induced headaches.
- DO NOT** use if you suffer from any genetic conditions of the eye.

## OPERATIONS

1. Use the mask only after thoroughly reading directions in this manual.
2. Keep eye protection on.
3. If any redness, itching, discomfort, etc. occurs during or after use, discontinue use immediately.
4. Do not attempt to use while the usb cable is plugged in.
5. Do not use more than 15min/day

## STORAGE & MAINTENANCE

1. Keep away from direct light, high temperatures, and excessive moisture.
2. To clean your mask you can use a makeup removal wipe or similar product that is meant for facial cleansing.
3. Store the mask in the provided soft shell case when not in use.
4. Keep away from children.

### LEGAL MEDICAL DISCLAIMER

No products sold through this website, cleopatramask.com, Cleopatra LLC, or any social media, including Facebook, Instagram, Snapchat, Twitter, Pinterest, Google, Youtube, TikTok, and LinkedIn, advertisements, testimonials, or reviews associated with cleopatramask.com or Cleopatra Mask (all collectively, "Cleopatra Mask"), are intended for use in the diagnosis, cure, mitigation, treatment, or prevention of disease or any other condition, or to affect the structure, function, or condition of the body. All products sold through Cleopatra Mask do not claim to impart any health benefits or cure disease. Always consult with your physician or other qualified healthcare professional before using any products sold through Cleopatra Mask. The contents of Cleopatra Mask such as text, graphics, images, information, and other material, are for informational purposes only ("Content"). All Content is not intended to be (1) a substitute for or (2) used as professional medical advice, diagnosis, or treatment. Always seek the advice of your physician or other qualified health provider with any questions you may have regarding any condition. Never disregard professional medical advice or delay in seeking it because of something you have read through Cleopatra Mask You must not rely on the information on cleopatramask.com or any Cleopatra Mask Content as an alternative to medical advice from your doctor or other professional healthcare provider. If you have any specific questions about any medical matter, you should consult your doctor or other professional healthcare provider. If you think you may be suffering from any medical condition, you should seek immediate medical attention. You should never delay seeking medical advice, disregard medical advice. If you think you may have a medical emergency, call your doctor or 911 immediately. Cleopatra Mask does not recommend or endorse any specific tests, physicians, procedures, opinions, or other information that may be mentioned on the Site. Reliance on any information provided by Cleopatra Mask is solely at your own risk. Nothing in this disclaimer will:(a) limit or exclude any liability for death or personal injury resulting from negligence;(b) limit or exclude any liability for fraud or fraudulent misrepresentation; (c) limit any liabilities in any way that is not permitted under applicable law;