

LIGHT THERAPY FACIAL MASK

BY

BEAUTY BY LIGHT



OVERVIEW

LED light therapy has been widely recognised as a safe, effective treatment for many skin conditions. Each LED light has a unique wavelength that penetrates the skin at a different level, providing cells with usable energy. This results in various reactions in the cells including improved cellular metabolism, reduction of inflammation and enhanced production of collagen.

As the light waves only enhance your skin's own natural abilities, the treatment is completely safe. Clinically proven and FDA approved, the technology has previously been confined to doctor's offices and specialist treatment salons. BEAUTY BY LIGHT brings this technology into an at-home mask, with the same benefits delivered by in-salon treatments.

What does it do?

Depending on which wavelength is used, light therapy can increase collagen production & cell renewal, diminish fine lines & wrinkles, reduce pigmentation & even skin tone, shrink large pores & smooth skin texture, fight acne causing bacteria and minimize scarring.

Is there proof that this technology works?

The efficacy of light therapy has long been acknowledged in the medical community, frequently used by dermatologists and beauty practitioners alike.

Initially developed by NASA, this technology has been scrutinised in clinical studies, scientific and medical research and independent testing, and has been approved for use by the FDA.

Previously available only through a doctor or specialised salons, BEAUTY BY LIGHT is now making this technology available for home use.

Are there any adverse effects from this treatment/is it safe?

LED devices are classified by the FDA as possessing a non-significant risk; they are completely safe and possess no possibility of harm from over-exposure.

The LEDs emit no UV light, so there is no risk of sunburn or skin cancer. This treatment has no adverse reactions and no down-time.

We recommend speaking to your doctor first if you are pregnant, have epilepsy or photoallergy or are taking any of the following medications; steroids injections, cortisone injections, medications causing light sensitivity (such as Tetracycline).

Is the treatment abrasive or painful?

Light therapy is a gentle, painless, non-invasive and non-abrasive treatment. You'll have no down time after each treatment.

How often should I use the mask?

We recommend using the red mode for 15 minutes followed by the blue mode for an additional 15 minutes daily. This combination will give the most visible results. After 30 days you can try combining the red mode with another colour as per your requirements.

How soon can I expect to see results?

Light therapy stimulates your own cells natural healing and collagen production, so the most visible results are not immediate. The skin will look different after the first treatment; smoother and more hydrated with a natural radiant glow, and with each treatment you will see improvement.

Acne scars and scars from other types of injuries can start to improve in as little as two weeks and should gradually fade over time with continued use.

Should I continue my skin care routine as normal?

Absolutely! The LED mask will help to strengthen and rejuvenate the skin at a cellular level so will amplify the results achieved in other facial treatments. We recommend using the mask as an addition to your existing skincare regime.

Is LED for everyone?

Light therapy will have the greatest impact on skin cells with sluggish metabolism and slow cellular activity, showing visible signs of aging. It can be used on all skin types, colours and conditions.

What is the difference between the BEAUTY BY LIGHT facial mask and in-salon treatments.

Not much! All LED technology will differ slightly in the specific wavelengths and strength of light emitted, but, the effects will mostly be the same. The light emitted by in-salon devices is the same light emitted by the BEAUTY BY LIGHT mask. BEAUTY BY LIGHT just gives you the freedom to complete your treatments on a more regular basis, in the comfort of your own home.

Can I use red light therapy after Botox injections?

We recommend that you wait two weeks after Botox injections before beginning to use the mask. If you have any concerns, speak with your doctor or dermatologist.