MAGNET-PERFECT™ THERAPEUTIC NECK MASSAGER

USER MANUAL

Magnetinnov
MEDICAL INNOVATION

TABLE DES MATIÈRES

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Magnetinnov
MEDICAL INNOVATION

TABLE OF CONTENTS

Operational Mode	15
Battery installation, storage, and maintenance	17
FAQs	18
Troubleshooting and solutions	19
Parameters and technical specifications	20
Use	22
Usage precautions	23
Warning	26

Operational Mode

FR **01 Power supply operation**

EN Open the battery compartment lid and place the batteries in the compartment

02 Start-up

1. The power switch must be “on” (**black button**) beforehand. Wipe your neck with a moist towel as well as the conducting sheet and the areas near the neck. The parts of your neck that are in contact with the device must be clean before each use.

2. Use the second "on" button (**red button**) of the handle to turn the massager on. Light 1 will turn on. Then, you can choose the intensity using the "high" or "low" buttons at the left side of the handle. The following predefined settings will turn on when the product is first turned on:

- A. Automatic mode on (mode 1)
- B. The lowest power level
- C. Session duration : 15 minutes

3. Mode button use

The switch button is accessible by pressing the “mode” button on the righthand side. There are six modes to choose from. The default setting is mode 1 (automatic mode) and the mode light will turn on to indicate which mode has been selected. Each mode’s intensity will be at the lowest speed by default. The six modes are:

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Mode 3 : mixing, the mode light 3 is lit;

Mode 4 : pressure massage, mode lights 1 and 3 turn on simultaneously;

Mode 5 : kneading, mode light 1, 2 turned on simultaneously ;

Mode 6 : beating, mode light 1, 2, 3 turned on simultaneously.

Automatic massage mode :

we recommend using the automatic massage mode on loop for best results

4. "High" button use (intensity increase)

Intensification of key intensity, each time you press the left "high" button, the intensity increases by one speed for a total of 16 speeds. This physical therapy involves the user choosing an appropriate intensity according to how they feel. We recommend using a low intensity the first time of use in order to get used to the sensation.

5. "Low " button use (intensity decrease)

The low button allows the user to decrease the intensity. Each time the low button is pressed, the speed decreases the applied speed intensity for a total of 16 speeds.

Suggested use: 10 to 20 minutes per day.

Battery installation

FR
EN

1. The product uses "AAA LR03" alkaline batteries that can be used for up to three months (15 minutes per day).
2. If you aren't using the product for a long period of time, remove the battery. Otherwise the battery leak could cause the device to malfunction.
3. Do not install the batteries in the wrong direction. Check the "+" sign. Use the battery disposal treatment method according to your city's applicable environmental protection laws.
4. Replace the batteries while the device is turned off.

Storage and maintenance

1. Keep out of reach of children.
2. Do not store in direct sunlight, high temperatures, and humid environments.
3. Keep the device in a dry, well-ventilated place.

To avoid causing an accident or malfunction, do not dismantle, repair, or change the product.

Q : What causes the tingling on my neck when I use the electrical impulse feature?

A : Tingling on your neck can cause skin dryness and can happen when the device isn't entirely placed in contact with your neck. When this type of discomfort from the electric impulse takes place, we recommend that you wipe your neck with a moist towel. The electrical force should be weak during the first application and then will progressively adapt (the pulsing electrical massage is used to massage the neck meridian, not the cervical vertebrae.)

Q : Why is my skin irritated?

A : The use of electrical impulses can cause dryness on the area of the skin in contact with the neck massager which can lead to irritation. We recommend that you wipe your neck with a moist towel while the device is in contact with your skin. This irritation generally disappears after two days.

Q : When do I have to replace the batteries?

A : The product uses "AAA LR03" alkaline batteries that can be used for up to three months (15 minutes per day).

Q : Why does my neck hurt during use?

A : The product is ergonomically designed and perfectly adapted to the human body, especially to the neck. Wear the product correctly during use and try to maintain good posture. Use the massager every day for 10-20 minutes to relax your neck.

Troubleshooting and solutions

FR

EN

Error	Check the following	Solution
The product doesn't turn on	Are the two main machine switches in the "on" position?	Check the number of switches that are in the "on" position and change accordingly
	Was the battery installed correctly?	The battery must be correctly installed according to polarity.
	Is there a low battery ?	Replace the battery if it's low.
	Are the battery and spring contact defective?	Use the appropriate tools to clean the battery spring.
Showing normal signals but they are imperceptible or non-reactive	Does the treatment instrument electrode come into direct contact with the cervical vertebrae?	The therapeutic instrument electrode should be directly exposed to the cervical vertebrae skin.
	Is there any hair, clothing, or other foreign objects caught between the electrodes and your skin?	Remove the foreign objects.
	Are the treatment instrument electrode and your skin firm?	Put the electrodes on your skin and adjust the treatment intensity and mode accordingly. Dampen your skin with a moist towel.
	Is it too long ?	The appropriate amount of time is 10-20 minutes.
Tingling	Is it too intense ?	Adjust the level accordingly.
	Is your skin sensitive or allergic to the electrodes?	Check if you have a history of skin allergies. If the sensitivity is slight, you just need to decrease the treatment duration. Severe allergies require stopping use since it requires desensitizing treatment first.
	Cold weather, dry skin?	Dampen your skin and the electrode surfaces with a damp towel.
La stimulation est très faible	Is the electrode clinging onto your skin?	Detach the electrode from your skin and adjust the treatment modes and intensity accordingly.
	Does the electrode have dirt or oil on it?	Clean the electrode surface.
	Is the battery low ?	Replace the battery if it's low.

FR

EN

1. Product specifications

Power supply : 3VDC30mA (2 batteries)

Equipment safety category: BF type applied part

Packaging box size: 220X200X60mm

2. Normal conditions of use

Temperature : +5 C ~ +40 C .

Relative temperature: less than 80%

Atmospheric pressure : 86kPa ~ 106kPa ;

3. Electrical properties

Impulse frequency : 0hz-1000hz.

Impulse duration: 20 μ s ~ 400 μ s.

Maximum amplitude of a single pulse output: >7 μ C

Maximum energy of a single pulse output: \leq 300 m J
Maximum output amplitude: < 50 V (50 m A)

Peak output voltage during opening measurement: \leq 500 V
Rated load impedance: 500 Ω , +10% tolerance

Surface magnetic force: no less than 0,01 T ;

Recovery time : 15 minutes, + 10% tolerance.

Utilisation

FR

EN

Remove the packaging and assemble the pulsing patches

Plug in the patches with the cable provided

Insert batteries

Turn the device on

Adjust the device to your liking

Turn on the patch switch

Change modes with the "mode" button

Use the "high" button to increase intensity.

Use the "low" button to decrease intensity.

Precautions

FR
EN

Danger

- Never use the product with the following electronic medical devices:
 1. Pacemakers and other implantable electronic medical devices;
 2. Artificial hearts, lungs, and other electronic medical devices that maintain vital functions
 3. Electrocardiographs and other electronic medical devices.

* These electronic medical devices can cause malfunctions and turn out to be fatal.

Warning

- If a high-frequency surgical device is applied at the same time as the massager, the massager's electrode could cause burns and damages to the device.
- If the product is used near a shortwave or microwave therapy (1M), its output may be unstable.
- The electrode should not be in direct contact with wounds or used on scars during treatment.
- To avoid the increasing the risk of cardiac fibrillation not place the electrode near your heart or chest while you use the product. Do not use it over your neck, head, mouth, genitals, or on rashes, and do not try to place it on the soles of your feet.

Warning

- Product use is recommended for only one user. If the product is being used by multiple people, clean it before use.
- Attach the electrode evenly and firmly on your skin while using the massager.
- If you don't feel well during or after use, stop using the product immediately and consult your doctor.
- To avoid accidents or physical illness, consult your doctor before using the product if any of the following cases apply to you:
 1. If you are receiving medical treatment and are dizzy.
 2. If you have heart or nerve problems.
 3. If your body temperature is over 100.4° F
 4. Infectious diseases.
 5. Rashes or skin sensitivity.
 6. DO NOT use while driving.
 7. Elderly people should use the product under the instruction of their doctor.
 8. DO NOT use the product after drinking or taking sleeping pills.
 9. NEVER try to inspect or open the product.
 10. DO NOT put heavy objects the product.
 11. DO NOT continue to use the product if it's damaged.
 14. DO NOT drop the product from a high place.
 12. DO NOT submerge or drop the device in water or any other liquid.
 13. NEVER use pins or other metal clips with this device.

Precautions

FR

Notes

EN

- If your skin or body feel strange, consult your doctor and follow their instructions.
- To avoid electric shocks, cut off the power supply before moving the device if you want to move the massager or put it on someone else during the massage.
- To avoid accidents, keep out of reach of children.
- Do not use while sleeping. Use while sleeping could cause bodily dysfunction, or the massager could move to other parts of the body and cause a loss of consciousness.
- Do not put the electrode in the same place for a long time. Doing so can cause skin inflammation.
- To avoid electric shocks, do not allow metal objects to come into contact with the electrode during the massage (necklaces, belts, etc.)
- Do not use the massager in the same place for an extended period of time (more than 30 minutes). Prolonged use can cause muscle fatigue and discomfort in the massaged area.
- Wear the product correctly by putting the device on the cervical spine in a symmetrical manner. Improper use of this product could cause accidents or physical illness.

Avertissement:

This symbol on a device or its packaging means: do not throw away electronic devices in unapproved municipal landfills: use separate collection facilities. Contact your local authority to obtain information on available collection systems in your area. If electronic devices are thrown away in landfills, dangerous substances can seep into underground water and end up in the food chain which is dangerous to your health and wellbeing.

Contact Us

If you have any questions, feel free to contact us, and we'll do the best we can to address your questions or concerns.

contact@magnet-innov.com

Magnet-Innov

[magnet_innov](https://www.instagram.com/magnet_innov)

magnetinnov
MEDICAL INNOVATION

*Thank you
for purchasing the
Magnet-perfect™
therapeutic neck massager*

Discover now our best products

Magnet-Relax 2.0
The Smart Eye
Massager

Magnet-Zen™
The Shiatsu
Massage Pillow

Magnet Relax™
The Massage
Mask

 Magnet-Innov

 magnet_innov

magnetinnov
MEDICAL INNOVATION
www.magnet-innov.com

 Scan me
for a gift