

**ULTRA BRACES  
ULTRA DOUBLE BRACES AND  
ULTRA BRACES FLAVOR FUSION**



**SHOCK  
DOCTOR**

DENTAL WARRANTY  
DISCLAIMER

**CERTIFICATE OF LIMITED WARRANTY AND DISCLAIMER**

**WARRANTY:** Shock Doctor, Inc. ("Shock Doctor") warrants the within mouthguard unit to be free from defects in material and workmanship which could cause damage or injury to the sound, natural teeth of the purchaser herein when properly used in a supervised athletic contest or coach-supervised training session. This warranty shall extend for a period of one (1) year from the date of purchase.

Shock Doctor's obligations under this warranty shall be limited to the lesser of:

- i. The actual dental expense incurred by the user; or
- ii. \$1563 per injured or replaced tooth, up to a maximum liability of \$50,000 per incident.

Shock Doctor's obligations hereunder shall be reduced by the amount of any insurance collected or payable for the damage or injury claimed. This certificate does not provide payment or indemnity for orthodontic treatments or preparation of a tooth to be used as an anchor for a bridge.

Shock Doctor shall have no responsibility or liability hereunder unless:

- i. The mouthguard was fitted by claimant in accordance with the printed instructions accompanying this warranty; and
- ii. A complete written notice of the injury or damage is received by Shock Doctor within thirty (30) days after the occurrence of the injury or damage; and
- iii. The notice of injury is substantiated by a statement of a licensed dentist providing necessary dental treatment within thirty (30) days after the date of injury or damage and a statement of a participating game official, if a supervised contest, or a supervising coach or school superintendent, if a practice session, that the mouthguard unit was in proper use at the time of injury; and
- iv. The within mouthguard unit is returned to Shock Doctor with the written notice of injury. The mouthguard unit shall become the property of Shock Doctor upon receipt of return.

This warranty is made in lieu of all other warranties, express or implied, including but not limited to the implied warranty for a particular purpose and any implied warranty arising out of a course of a dealing, custom, usage or trade. Except as is otherwise expressly provided herein, Shock Doctor disclaims all liability for loss, injury or damage based upon negligent design or defects in materials and workmanship for the within mouthguard.

*Keep this dental warranty, disclaimer and store receipt in your records. In case of injury, mail mouthguard, store receipt, the notice of injury and accompanying statements to Shock Doctor, Inc.*

Shock Doctor, Inc., 11488 Slater Avenue, Fountain Valley, CA 92708  
Phone: (800) 233-6956 | Fax: (888) 446-5999 | [www.shockdoctor.com](http://www.shockdoctor.com)

# ULTRA BRACES ULTRA DOUBLE BRACES AND ULTRA BRACES FLAVOR FUSION



**SHOCK  
DOCTOR**

FITTING INSTRUCTIONS

## WARNING!

### DO NOT CHEW YOUR MOUTHGUARD

In order for a mouthguard to protect, fit, and last it must be worn properly on your teeth and gums as indicated in the fitting instructions. If any mouthguard is chewed or not worn properly, its protection, fit and durability will be greatly diminished.

### TO FIT THIS MOUTHGUARD YOU WILL NEED:

- > Boiling water
- > Slotted spoon
- > Saucepan
- > Cold tap water
- > Mirror
- > Timer

### FITTING INSTRUCTIONS:

1. Boil 4 cups or more of water in a saucepan.
2. Bring water to full boil.
3. Remove from heat and let sit and still for 30 seconds.
4. Attach tether to Ultra Braces mouthguard and submerge in stilled hot boiled water.
5. Use tether to hold mouthguard under water for 20 seconds.
6. Remove from water. Shake off excess water. Cool for 2-3 seconds.
7. Place in mouth facing a mirror so mouthguard is centered.  
**\*NOTE: For Ultra Double Braces, place in mouth with blue translucent gel facing up.**
8. With mild pressure, bite down along full arch of mouthguard all the way to back molars making sure to push with fingers the inner brace channel of mouthguard over the top of brace brackets on teeth. (Double Braces lower channels sits naturally on top of lower brackets)
9. Keep biting down while sucking in air and remove mouthguard after 20-25 seconds.
10. Drop into glass of cold water with ice cubes and let sit for at least two minutes and remove and insert in mouth to ensure mouthguard channel is aligned over the brace brackets and fit overall is snug.

### For Re-fitting if first attempt was not satisfactory.

- Follow Steps 1-4 from fitting instructions.
- Step 5. Immerse in hot water for 15 seconds.
- Follow Steps 6-10 from fitting instructions.

### For Re-fitting after brace brackets are adjusted by orthodontist.

- Follow Steps 1-4 from fitting instructions.
- Step 5. Immerse in hot water for 20 seconds.
- Follow Steps 6-10 from fitting instructions.

NOTE: Mouthguards are best stored in mouthguard cases while not in use. Clean your mouthguard regularly by washing under hot water. If using soap, rinse the mouthguard thoroughly before wearing.