


PRODUCT OVERVIEW

Thank you for choosing the OneCompress Smart Hand Massager. For ease of use, please read and understand all safety instructions prior to using the product. The images in this manual are for illustration only.



QUICK START

- STEP 1:** Insert your hand inside the machine using the hand entrance.
- STEP 2:** Place your palm in the center of the massage pad.
- STEP 3:** Press the power button (👉), all key displays will light up.
- STEP 4:** Choose a Size Option: Large (👉) or Small (👉)
- STEP 5:** Press the Intensity mode button (📶), There are 3 different strength modes: Gentle mode, Soothing mode and Strong mode.
- STEP 6:** For heating function, press (🔴) and adjust the appropriate temperature.
- STEP 7:** Please press (👉) to turn off the machine.



TIPS

- If no other button is selected after pressing the power button - product will turn off automatically after 5 seconds.
- When the product is charging, the power button will flash until fully charged.
- Please do not use the machine while charging.

TROUBLESHOOTING

If the product does not work normally, please check whether the following conditions fix the issue. If the suggestions provided cannot solve the problem, please return the product to Onecompress for an exchange.

Issue / Malfunction	Possible Cause	General Solution
The massager cannot be turned on	The battery has ran out of charge or is exhausted	Please charge. Restart after charge is restored
The product can be turned on, but does not work	Power may be weak	Please charge
	Did not press Size Option selection key	Press Size Option mode button within 5 seconds of Power On
Massager cannot be recharged	Power adaptor connected incorrectly	Re-plug until LCD power display flickers
	Is data line and USB socket intact?	Check socket. If there is damage, please contact us

Attention:

For best results, limit use from 15-30 minutes per hand. Use charger provided or DC5V 1 ~ 4A charger. If the power cord is damaged, stop using it immediately. Contact service advisor for maintenance. Store in a dry ventilated place.

DO NOT USE: when hands are wet, while driving, under influence of alcohol/sleeping pills. Follow professional medical advice if: Recently had hand surgery, swelling and inflammation, congestion, severe arthritis, trauma and other hand diseases; people with heart, brain and lung diseases; blood diseases with bleeding tendency, osteoporosis and soft tissue diseases; stroke; infants and pregnant women.

Please review safety precautions here:
www.onecompress.com/pages/hand-massager-safety

Questions? Email us at: support@onecompress.com

