

# *Osaki* KneeMedic



## Company Information

*Osaki*

Contact Number: 1-888-848-2630

Website: [www.titanchair.com](http://www.titanchair.com)

# *Osaki* KneeMedic




Superconducting Acupuncture + Infrared Physiotherapy

+ Vibration Massage







Superconducting Acupuncture + Infrared Physiotherapy+ Vibration Massage  
**KneeMedic**

**In order not to subject with danger or material loss, please use after careful reading.**










-  Danger: Ignorance to this sign may result in serious body injury or death.
-  Warning: Ignorance to this sign may result in serious body injury or death.
-  Attention: Ignorance to this sign may result in body injury or material loss.

**Legends**






-  Restricted
-  No decomposition
-  No touch
-  For strict compliance

**Attentions for users**

**The following items must be complied:**

-  Please ensure power supply from customized socket.  
(Both 110V and 220V compatible. Can be conveniently used anywhere anytime)
-  Unplug from the socket in case of no use.
-  Never touch power plug with wet hand.
-  Never sleep during use of the product (possibly subject with burning injury).
-  Please consult doctor first in case of acute swelling or bleeding with the knee.
-  Consult doctor first in case of rheumatoid arthritis.
-  Patients suspected with swelling, bacterial infection or other special inflammation shall consult doctor first.
-  Upon more severe pain, stop use and contact the customer service center.
-  Don't directly touch skin, but wear on underpants or the given safety pad. (Temperature above 60°C may result in burning injury, for your attention)  
The shining part may result in burning upon touch to electric wire, cotton fabric, plastic or other flammable materials, for your special attention.

**The following items are for restriction:**

-  Don't expose to sharp items.
-  Don't throw or strike, otherwise the Unit may not be turned on properly.
-  Don't use within sauna room or bathhouse.
-  Don't wash within water. (Dry cloth may be used to sweep off, if in need of cleaning)
-  Don't self-decompose the product or change its shape (no post-sale service will be provided in case of self-renovation)  
Don't bend the power plug or put below heavy items.

Superconducting Acupuncture + Infrared Physiotherapy+ Vibration Massage  
**KneeMedic**

**Characteristics**

- New-concept medical treatment unit for muscle joint, penetrating thermal energy via warm radiation from the shining radiating light and infrared lighting into the deep zone of knee/shoulder/elbow, which can rapidly improve pain at the muscle.
- Powerful built-in infrared light (700~50000nm), effectively penetrating the internally generated thermal energy into the deep zone of knee/shoulder/elbow and thus improving the symptom.
- High-concentration silver with high thermal conductivity on the warm guide can rapidly pass the thermal energy into the knee/shoulder/elbow.
- According to ergonomic method, comfort wearing will improve the performance.
- The control furnished with the temperature memory device can adjust the temperature between 45°C and 65°C, or maintain at a temperature suit to the user.  
-----In case of long time of use, the set temperature may be exactly maintained.
- Strong vibration massage treatment and magnetic therapy.

Deliver thermal energy into the joint,  
through infrared lighting.

Strong vibration massage  
treatment and magnetic therapy.

Internal infrared light and warm guide  
that deliver thermal energies by  
penetration into the knee/shoulder/elbow  
joint at 700-50000nm wavelength




**Parts/ accessories**

(the Adaptor and Power Cable may be changed according to different market)



Adaptor



Safety pad



Shoulder wear tape

**Method of use**

1. After wearing the safety pad, wear up the temperature controller on the knee. Adjust the knee belt until feeling comfortable.  
-----For best performance, try to be close to the knee.
2. Connect the AC adaptor with the power plug, subsequently connect to the power port of the control (machine body).  
-----Control on display, upon power plug into socket.
3. Press once the most left button on the control (START/STOP), for commencement of heat radiation and infrared lighting.  
-----Each press of the button will cause switch between ON/OFF.
4. Adjust to body comfortable temperature by pressing the middle button (Up) and most right button (Down).  
-----Initial treatment is suggested to be under 55oC.
5. In the later stage the Up and Down buttons are still to be used to adjust the temperature between 45 C - 65 C .  
-----For temperature within 45 C -55 C , the unit increment is 5 C ; for temperature within 56 C - 65 C , the unit increment is 1 C , for better accuracy.
6. Upon completion of use, press once the most left button START/STOP to automatically stop.
7. There is a built-in auto timer within the control, which will automatically turn off power after use of 30 minutes. For any extension, press one more times of the START/STOP button.  
-----For the time being, the temperature memory device will save the temperature set at power off. No need re-set temperature.
8. Click VIBRATION button on the inner side of control panel to start/stop the vibration massage.
9. Same usage for shoulder and elbow treatment.

Superconducting Acupuncture + Infrared Physiotherapy+ Vibration Massage  
**KneeMedic**

Front elevation



Ergonomic, convenient to used, close touch, and easy to be operated



Ergonomically designed infrared light and warm guide, for convenient penetration

Superconducting Acupuncture + Infrared Physiotherapy+ Vibration Massage  
**KneeMedic**

**Suggestions for use:**

1. Wear for a while before movement or going out, which will save the sense of pain at the knee/shoulder/elbow and muscle and help various activities in daily life.
2. For expected performance, it's suggested to wear for 2-3 times everyday at warm temperature, 30 minutes each time.  
-----Warm temperature is the most proper temperature for the user. Don't set over-high temperature.
3. Acute temperature change will discount the performance. Try to progressively adjust the temperature.  
-----Don't make too high increment or decrease at temperature.
4. Upon relief of the sense of pain, it's suggested to conduct extension exercises by continuously unbending and bending the knee/shoulder/elbow joint.
5. Shall wear on slim underwear (pant) or the safety pad provided together with the product.  
-----Long time of direct exposure to skin may result in burning injury or skin sensibility.

**Target population:**

- With sense of pain during walking on stairs or moving the knee joint
- Extreme pain at the knee/shoulder/elbow joint during climate change or cloudy weather (at low pressure)
- Frequent sense of discomfort at the joint or lack of flexibility
- Knee/shoulder/elbow swelling or rigidity in the afternoon or night after body movement
- 

A new kind of physiotherapy unit totally different from the old simple knee protecting belt or normal hot-laid knee heater.

**Can effectively relieve pains at knee/shoulder/elbow muscle**