

• ABOUT US •

Driven by the passion for innovation, we at Dr Trust endeavour to provide our customers with the latest medical inventions with an objective to promote good health and wellness all around the world. All the medical devices and health monitors provided by Dr Trust are supported by accurate, latest and ground breaking technologies, innovated at our headquarter in NY, USA. All our products adhere to the most stringent CE and FDA guidelines and are strongly recommended by doctors and health practitioners. Our products are designed in the utmost exemplary ways to ensure that their accuracy and convenience are unrivalled. The ease of their use and operation makes them even more suitable for users of all age groups.

Dr Trust strives to enhance the quality of lifestyle by providing with the most trusted and innovative health care and wellness products. Being a renowned global leader in health care products, Dr Trust ensures that our technically efficient team works dynamically and tirelessly to provide the best of the medical devices to our clients. The products that we have to offer are suitably designed for use at homes, laboratories and hospitals.

Our ground breaking solutions allow you to monitor your health in the easiest ways possible. In today's era when all of our lives are too hassled to handle, it becomes a bit difficult to pay attention to our health. But it has now become easier with the coming of the monitoring devices which can be conveniently used at homes and even on the go.

We bring to you a variety of best self medical devices, trusted and used by Doctors, medical professionals and home users all over the world.

Dr Trust USA  
Grow Buddy Weighing Scale - 510



• QUICK START GUIDE •

Weighing and Measuring A Child (Up to 2 Years)

**Step 1**  
Clean and put together the baby tray on the platform of the weighing scale.

**Step 2**  
Place a soft towel/ blanket on the tray and press the ON/OFF/TARE key to turn on the scale.

**Step 3**  
When the scale shows 0, place your baby on the blanket gently.

**Step 4**  
Check the screen of the scale for weight measurement.

**Step 5**  
To measure the baby's height, use the scale given on the baby tray.



Adult Weighing

**Step 1**  
Detach the baby tray from the platform to use it as a normal scale.

**Step 2**  
Step on the scale and stay still until the screen locks your weight.

*Note: This is a baby cum adult weighing scale with max 150 kg weight capacity. It is suitable for weighing babies, children, grownups, and pets.*

• OVERVIEW •

Dr Trust Grow Buddy Scale – 510 can be used to monitor the body weight of all family members on a daily basis. It has high precision sensors so that users can get very accurate results every time when they use it. The maximum weight capacity of it is 150 kg and shuts off automatically when you do not use it. The weighing scale comes with a removable tray so that it can be used to monitor a baby's weight in an easy way. It has all the needed features with amazing looks and durable built.

• ON/OFF/ TARE BUTTON •

Press the ON/OFF button to start the scale switch it ON/OFF. It also switches ON automatically when you step on it or keep your baby on the weighing tray. The scale switches OFF automatically after you have weighed or step off the scale.

• LCD DISPLAY •

The power saving blue colour backlight turns ON when you switch it ON and it switches OFF within 1 minute.

• THE UNIT BUTTON •

Press the UNIT button to switch between kilograms (kg), pounds (lb), and ounce (oz) easily.

• THE HOLD BUTTON •

You can select the MECHANICAL HOLD function by pressing the "HOLD" button. Within 3 seconds, the word "HOLD" stops flashing and the weight is displayed for a further 20 seconds.

• WEIGHING TRAY •

The weighing tray is easy to remove from the machine.



• BATTERY COMPARTMENT •

Fit in 3 AAA batteries to power up the scale. The battery needs to be replaced immediately after 'LO' appears in the display.

• GETTING STARTED THE SCALE •

Installing Batteries

1. Open the battery cover on the backside of the scale.
2. Put three "AAA" batteries in the right polarity.
3. Place the battery compartment cover back and the scale powers ON automatically.

Start Weighing

- Press "ON/OFF" button to switch on the scale.
- As soon as "0.000 kg" appears on the display, the scale is ready to take measurement.
- Now place the baby on the plastic tray.
- The weight will be shown firstly, the weight and D-value will be shown alternately twice, and then the scale will be switched off automatically.

Using Mechanical Hold Button

- This baby scale had a mechanical HOLD function. You can select the MECHANICAL HOLD function by pressing the "HOLD" button.
- Within 3 seconds, the word "HOLD" stops flashing and the weight is displayed for a further 20 seconds.

- If you do not want to place the baby directly on the weighing surface but prefer to use a towel.
- If the towel (and any other items which you want to deduct again during weighing) weighs over 50 g, you can also use the TARE function. Once you have pressed the TARE button, 'TARE' is shown in the display.

Resetting the Scale

- To reset the scale to 0.000 g, press the TARE button after "HOLD" has stopped flashing.
- On delivery, the scale is set at "kg" units.

Unit Conversion

You can choose to display your weight in kilograms (kg), pounds (lb), and ounce (oz) easily. The Unit button allows you to switch units as per choice.

• CAUTION •

- The scale is mainly intended for kids' use. Although adults can also use it.
- Do not use any chemical cleaning agent.
- Always use the scale on a flat surface; do not use it in an upright position.
- Keep the unit out of reach of infants or children.



- Avoid using this device in wet conditions.
- Do not use it on slippery wet floors.
- Never put heavy objects on the scale.
- The scale can accept a maximum load of 150 kgs

• TECHNICAL SPECIFICATIONS •

Technology	Equipped with a High Precision "Strain-Gauge" Sensor
Display	Large LCD Display
Max Capacity	150KG/330LB/529IOZ
Units	Weighing Units Switch Between KG/LB/OZ
Power	3 AAA batteries

• CUSTOMER SUPPORT •

CONTACT ADDRESS

USA  
Nureca INC.USA  
276 5th Avenue, Suite 704-397,  
New York (NY) - 10001, USA

INDIA

Corporate Office (Mumbai)  
Nureca Limited  
Office No. 101, 1st Floor, Udyog Bhavan  
Sonawala Lane, Goregaon East,  
Mumbai, Maharashtra - 400063

Contact us

India: +91-7527013265 / +91-9356658436  
Website: [www.drtrust.in](http://www.drtrust.in)  
Corp Website: [www.nureca.com](http://www.nureca.com)  
Customer Care Executive Email: [customercare@nureca.com](mailto:customercare@nureca.com)  
Connect with us on social networks  
Facebook: @drtrust  
Instagram: @drtrustisn  
Youtube: NurecaUSA

COPYRIGHT©2022 NURECA LTD ALL RIGHTS RESERVED



Scan to View  
Product Demo Video  
[www.drtrustusa.com/510](http://www.drtrustusa.com/510)