

**Function**

- Portable one key to start design, easy to use
- Four high-frequency massage head, massage head can emit light and heat
- The positive and negative kneading humanoid massage
- Flexible and convenient to move, can achieve full touch of your body , the effect is more obvious
- This neck massage pillow has a built-in safety electronic system.
- Smart design equipped with rollers to offer a kneading massage for neck, shoulder, waist etc.
- Two different kneading massage and no-step speed regulation relaxing your body
- Infrared heating therapy promotes blood circulation.
- Relieves muscle pains & fatigue and improve sleep, etc.

**Cleaning and maintenance**

- Be sure to clean the machine after use and before storage
- Unplug the machine from the socket
- Store in a clean, dry space.
- Avoid dust from entering or accumulating on the massager.
- Use a clean, dry cotton cloth to dust or clean the machine
- Avoid any kind of liquid intrusion or contact with the machine
- Please do not use volatile oil, bleach or other corrosive solvents to wipe the machine.

For any other query Contact: [customercare@nureca.com](mailto:customercare@nureca.com)

**SHIATSU DEEP KNEADING  
MASSAGER PILLOW WITH HEAT**



**Instruction Manual**

- A. Warning and Attention**
- Women should not use massager during menstrual period and pregnancy
  - Patients with severe heart disease or who have undergone relative surgeries should refrain from using the massager
  - Patients of Hemorrhagic disease, sepsis, patients with fever should not use the product
  - Do not use if skin is damage or inflamed
  - This product cannot be used by physically distempere, sensorially or mentally disordered, or inexperienced persons (including children), unless they are under the supervision or instruction by adult(s) in charge of their security.
  - This machine has heating function, please exercise caution and pay attention during usage to adjust heat to a comfortable level for your usage
  - Please use massager carefully not adding too much body pressure on the massager which may cause motor damage to the product
  - After a massage please drink 500 ml of water or tea, so as to promote the body metabolism
  - Do not insert your fingers into rotation massage head bottom, in order to avoid the task
  - Do not to remove the internal components of the unit, it may cause malfunction and also void the warranty of the product
  - Please do not try to self repair the unit in case of issues, immediately contact Lifelong Customer Care. Please read the Warranty details carefully.
  - When massager is not in use or post usage, please unplug the massager from the power socket
- Instructions for use of Product**

- Do not use wet hands to operate the massager, especially when plugging or unplugging it from the electric socket.
  - The product have owner timer and goes off automatically after 15 minutes.
- B. Power Supply and Cord**
- Do not use beyond the rated voltage of the power supply, it may cause fire or damage to the product.
  - When pulling the power cord, one should hold the power cord and then unplugging the power cord. Do not pull the cord from the socket with a jerk.
  - If the power cord is damage or exposed, please bring it to immediate attention of professionals or the manufacturer for maintenance or replacement. Please refer the Warranty section for more details.
- C. Operating Instructions**
- Insert the plug into power socket.
  - Press the ON/OFF key to open the massager function with heat.
  - Press the ON/OFF key again, switch to reversing mode
  - Press the ON/OFF key for 3 seconds, massager is only kneading without heat.
  - Press the ON/OFF key one more times, turn off the massager.
- D. Packing - Products include: 1 Pcs massage pillow, 1 pcs Adapter, 1 Manual**
- E. Technology**
- Voltage: DC12V
  - Power: 24W