

Quick Start Guide

* When using shot mode and booster mode at the same time, use booster mode after using shot mode.

SHOT MODE



STEP 01

After washing your face, in a dry skin condition rub very lightly along the contour using the head ball.

*Shot mode does not work on wet surface.



STEP 02

Contact the head ball on the skin surface and rub lightly along the contour where you want to improve.

BOOSTER MODE



STEP 01

Apply skincare product on the contour where you want to improve. Turn on the device, and switch the mode to booster mode, then use it.



STEP 02

Gently massage with the head ball on the skin surface, and then apply your daily routine skincare products.

Firm eye contour effect and innovative dual volumizing technology
Use **MEDICUBE AGE-R LINE SHOT** as following.

SHOT MODE



After washing your face, in a dry condition rub very lightly along the contour using the head ball.

BOOSTER MODE



After the SHOT mode, apply skincare product, and massage with the head ball on the skin surface.


Touching Method



SHOT MODE

Lining

Contact the head ball on the skin surface, rub very lightly along the contour where you want to improve.



BOOSTER MODE

Circling

Contact the head ball on the skin surface, massage as if drawing a circle.

Read before use!

It is recommended to start from level 1 when using it for the first time. Please adjust the level appropriately while continuously using it.




Caution

Avoid using around the Adam's apple or on areas with cuts or open wounds of any kind.



Recommended starting level

Start slowly from the lowest setting and find the appropriate level for your skin



LINE SHOT Use Cycle

SHOT MODE : 1-2 times a week
BOOSTER MODE : Everyday