

## Quick Start Guide



### STEP 01

After washing your face gently, applying the skincare product on your face.



### STEP 02

Fully attached the head to the skin surface and rub the skincare product to be absorbed into the skin.



### STEP 03

Refer to the boosting method, please use the device until the skincare product on your face is fully absorbed.



### STEP 04

When the skincare product is fully absorbed, finish with intensive care.

Bring out the natural glow of the skin with AGE-R's deep absorption technology



After cleansing, apply skincare product thoroughly. Massage gently with the head of device for boosting effect **in closed contact** with the skin surface.

## Boosting Method



### Watering

Push the head of device against the skin surface and rolling it like a drawing the circle on the skin.



### Boosting

Push the head of device against the skin surface and sliding it on the skin.



### Focus

Push the head of device against the skin surface and pressing it one spot at a time on the skin.

## Read before use!

It is recommended to start from level 1 when using it for the first time. Please adjust the level appropriately while continuously using it.



### Caution

Please avoid using uvula (Adam's apple), wounds, and surgical area.



### Recommended starting level

Start slowly from the lowest setting and find the appropriate level while continuously using it.



### BOOSTER-H Use Cycle

You may use the Age-R Booster-H 1-3times a day. Please gradually increase the frequency of using.