

Quick Start Guide

Face Slimming → **SLIM**

Face Up → **UP**

Body Contouring → **BODY**



STEP 01

After cleansing and applying toner, apply the Booster Gel to the roller heads of the Age-R Device.
**Gel must be applied in order to feel the energy.*




STEP 03

Apply the device to the desired area and massage in an upward direction.
**Recommended order of use: Jawline -> Laugh Lines -> Cheeks -> Crow's Feet -> Temples -> Forehead*



STEP 02


Select desired mode and intensity level. **When using the product for the first time, start at the lowest setting.*
**When the device is in operation, it will emit a notification sound every 5 seconds.*



STEP 04

After 10 minutes of use, pat the skin to allow the gel to fully absorb.
**If the skin feels dry, apply a moisturizer to finish.*

Getting the Best Results from Your Age-R Device



Jowl & Jawline

Massage from the tips of the chin to the bottom of the ear, pulling the skin gently in an upward motion.



Smile Lines

Divide the area into 2 horizontal zones, and massage both areas in an upward motion.



Cheek Care

Divide the area into 3 vertical zones, and massage each area in an upward motion from the jawline to just below the eye.



Crow's Feet

Gently roll in an upward motion to relax the muscles.



Temples & Forehead

Divide the area into 3 vertical zones, and massage each area in an upward motion.



Shoulder & Trapezius Muscle

Massage the muscles by rolling back and forth in upward and downward motions.

Mode Settings

Each mode designed to replicate the feeling of an esthetician's acupressure massage techniques.



Slim Mode



Up Mode






Body Mode

Mode Description	Mode
<p>SLIM</p> <p>Slim</p> <p>Stimulates face and helps with slimming the face</p> <p><i>*This mode is designed to simulate the Tapotement massage technique.</i></p>	
<p>UP</p> <p>Up</p> <p>Help to stimulate face to improve the appearance of the face's contours</p> <p><i>*This mode is designed to simulate the Petrissage massage technique</i></p>	
<p>BODY</p> <p>Body</p> <p>Helps to relax and improve the appearance of the body's contours</p> <p><i>*This mode is designed to simulate the Friction massage technique.</i></p>	

Read Before Use

The Age-R Device works by moving and massaging your skin.



- When using the device for the first time, start slowly from the lowest setting.
- It is highly recommended to test the device before using on the face.

**Both roller heads must maintain contact with the skin simultaneously for the medium-frequency-energy to be emitted successfully.*