

WARNING

DO NOT CHEW ON YOUR BRACES MOUTH GUARD

Your guards protective properties will diminish significantly if it is chewed or not worn properly.

FITTING INSTRUCTIONS:

No fitting/molding required. The guard is ready to use out the box. Gently fit over braces while making sure the interior channel properly aligns with your brackets. If you feel any pressure points after fitting, carefully trim the guard with a pair of scissors to achieve a perfect fit. When removing, peel the guard off instead of pulling it out.

USING THE TETHER:

Insert the thinnest end of the tether into the front locking system of the guard. The locking system can be cut with scissors if not used.

CLEANING METHODS:

Use a toothbrush and toothpaste to thoroughly scrub and clean your braces mouth guard or use mild soap and wash mouth guard in cool water. If you have mouthwash on hand, rinse after every use and let dry before storing.

Make sure to read the following safety instructions

1. Made of FDA approved silicone.
2. Visually check your mouth guard after every use. Cracks, splits or separations could render your mouth guard useless.
3. Visit the doctor immediately if you experience any unusual irritations after using the mouth guard.
4. This mouth guard does not contain any latex or materials that are currently known to cause allergic reactions.

SHARE YOUR NXTRND EXPERIENCE

Show us what you've been up to with your Nxrnd gear. Snap a photo of your Nxrnd product and share it with the community in our Amazon review section!

Product review link:

<http://amazon.com/ryp>

CONTACT INFO

OUR TEAM OF PROFESSIONALS IS ALWAYS HAPPY TO HELP



Support@nxtrndusa.com

OUR SOCIALS



@NXTRND



@NXTRNDUSA



@NXTRNDUSA

WWW.NXTRNDUSA.COM

NXTRND