lumie[®] Brightspark[™]


USER GUIDE

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BRUKSANVISNING

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USER GUIDE

Brightspark

Introduction	3
Precautions	4
Safety	5
Getting started	6
When should I use it?	6
What time of day?	7
How long?	7
Problem-solving	8
Maintenance	9
Technical description	10
Contact us	11

Introduction

Using Brightspark will provide you with the bright light your body needs during the dark winter months. It can help to boost your energy levels, put you in a better mood and make you feel more awake¹. Light therapy is safe and is proven to help alleviate the winter blues and SAD (Seasonal Affective Disorder)².

Please read these instructions carefully to get the most out of your Brightspark.

¹ J Affect Disord. 2000 Jan-Mar;57(1-3):55-61.

² Am J Psychiatry. 2005 Apr;162(4):656-62.

	Serial number
	European Declaration of Conformity: Class IIa Medical Device
	Manufacturer: Lumie, 3 The Links, Bar Hill, Cambridge, CB23 8UD, United Kingdom.

	Do not place in household waste. Dispose of at a recycling centre.
	Read instructions for use
	Caution
	Keep dry

	Fuse
	AC (alternating current) power supply
	Push ON/Push OFF switch
	For indoor use only

	Member of the 'Green Dot' recycling scheme
	Recyclable

Precautions

Research shows that light therapy is very safe. However, please check with your doctor first if any of the following apply:

- You have a major eye problem or have had eye surgery;
- You have been advised to avoid bright light or are taking medication that may cause photo-sensitivity;
- You have been diagnosed with depression or other mood/sleep disorders. This is especially important if you are also taking medication for this, since light therapy can also affect your mood.

A few people experience side effects but these are generally mild (including headaches, eye strain, nausea and hyperactivity) and can usually be solved simply by switching off the light. Over the next few days, gradually build up the time you spend in front of it to find a treatment time that works best for you. Some people complain of difficulty sleeping and often this is down to using the light too close to bedtime.

If you have any problems while using Brightspark, or if you've been using the light for a week but haven't noticed any benefit, please call Lumie or speak to your doctor.

IMPORTANT: the guidelines in this leaflet are not medical directions for treating any condition. If your doctor has suggested using light therapy, follow their instructions and discuss any effects with them.

Safety

- This product must be connected to a supply mains with protective earth.

- Keep the unit away from water and damp.

- For indoor use only.

- Keep the cable away from heated surfaces. If the cable is damaged it must be replaced by Lumie, a Lumie distributor or suitably qualified person in order to avoid a hazard.
- Children using the light, or playing near it, should be supervised.
- Do not attempt to service any part of the unit.
- Not to be used by or on individuals who are unconscious, anaesthetised or incapacitated e.g. under the influence of alcohol.

Electromagnetic compatibility

- This product conforms to the electromagnetic emissions compatibility (EMC) standard for medical devices and is unlikely to cause interference affecting other electrical and electronic products in its vicinity. However, radio transmitting equipment, mobile phones and other electrical/electronic equipment can interfere with this product; if this happens, simply move the product away from the source of interference. In rare cases where product function is lost due to EMC effects, disconnect the product from power source for 5 seconds and then restart. If the problem persists, please contact support. For further detailed information on the electromagnetic environment in which the device can be operated refer to the website: www.lumie.com/pages/e-m-c.

Disposal - products and bulbs

- At the end of its life, do not throw the unit out with normal domestic waste. Please take it to an official recycling point to be disposed of responsibly.

Getting started

If the product is damaged please contact us before use.

- ① Position Brightspark at arm's length (about 50cm/20in) with the plastic screen angled towards your face. Plug into the power outlet and switch on.

Although you don't need to stare at it, the light must reach your eyes to have any effect. This means you can use Brightspark whilst you're working at the computer, watching TV, reading, etc. Dark glasses, tinted lenses or closing your eyes will all limit the effects of light therapy.

You should notice positive benefits after using Brightspark for three or four days in a row.

When should I use it?

It's best to tackle the winter blues and SAD as soon as possible. Look out for early signs - often in September or October - and if you start to feel lethargic, anxious or irritable, start using your Brightspark.

Most people can skip a day here and there, especially when the weather is bright. As spring comes around you will probably start to cut back your light therapy; your body will tell you if it's too soon and if you notice your symptoms returning simply use your Brightspark again for a few days.

Light therapy is proven to put you in a better mood, to boost alertness and to help you feel more awake. Use Brightspark wherever and whenever you need to restore concentration and give yourself a mental boost.

What time of day?

If you're using Brightspark to help with the winter blues or SAD, bright light in the morning has been shown to be most effective so try that first. For some people, early evening light exposure is beneficial. Whether you're using Brightspark for winter blues, SAD or as a general pick-me-up, once you've found a pattern that's right for you, try to use it around the same time every day.

I find it difficult to get going in the morning

Switch on as soon as you can. If you don't have much time in the morning, then you should aim for a top-up or full session in the afternoon.

I feel sleepy or fall asleep too early

Use Brightspark in the late afternoon/early evening. It's best to avoid bright light too close to bedtime - within about three hours - as you may find it difficult to get to sleep.

How long?

Everybody's different but, if you use Brightspark at about 50cm/20in as suggested, 40-90 minutes will probably be enough for effective SAD treatment. There's no harm in using it for longer and you can even make up this time by dipping in throughout the day.

If you're more comfortable with it further away, then less light will reach your eyes and you'll need to use it for longer.

If you find your eyes ache slightly, try moving it further away or switch off for a while.

Problem-solving

The light is less bright or doesn't come on at all

- One or both of the bulbs may have failed. Switch off and unplug the lightbox. Let it cool if necessary.
- Lay it on a flat surface and gently squeeze the plastic screen lengthways to remove it from the case.
- Check that the bulbs are firmly in the sockets. If you have to adjust them, plug in the lightbox and switch on afterwards to see if this has worked.

The bulbs are firmly in place but it hasn't fixed the problem

- You may need to replace the bulb(s). To remove a bulb, push the red button down hard and slide the bulb out of the socket.
- Gently slide the replacement bulb into the socket until it clicks into position.

I've replaced the bulbs but it's still not working

-  It is possible that the plug fuse might have blown. The correct fuse is a 3A fuse to BS1362.
- If your Brightspark keeps blowing fuses, there is a problem and you should contact us.

If this doesn't work or if you have any other problems with your Brightspark, please call Lumie.

Do NOT remove the reflector as there are no user-serviceable parts under the reflector.

Maintenance

Replacement bulbs

The bulbs will slowly lose their brightness so, to maintain maximum output and benefit, we recommend you replace them every two to three years. Replacements are available from Lumie.

Cleaning

Clean only with a dry cloth or duster.

Warranty

This unit is covered under warranty against any manufacturing defect for three years from the date of purchase. This is in addition to your statutory rights. If possible, use the original packaging (including internal pieces) to return your Brightspark for service. Make sure the plug is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please call Lumie for further advice.

Technical description

Specification subject to change without notice.

Lumie Brightspark

230Vac, 50Hz, 50W

Lamp: 2 x PLL/24W/865

Operating temperature: 5 to 35°C

Transport/storage temperature: -25 to 70°C

Class IIa Medical Device

Contact us

Lumie is Europe's leading specialist in light therapy. Our products are based on the latest scientific studies and we focus solely on light therapy and its benefits. If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on +44 (0)1954 780500 or email info@lumie.com.

Lumie
3 The Links, Trafalgar Way
Bar Hill
Cambridge
CB23 8UD
UK

Lumie is a registered trademark of:
Outside In (Cambridge) Limited, incorporated in England and Wales.
Registered Number: 2647359. VAT No: GB 880 9837 71.

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© Lumie 2015
Outside In (Cambridge) Limited, met rechtspersoonlijkheid in Engeland en Wales.
Outside In (Cambridge) Ltd