

Daffodil®

User Manual

Model No.: HPC300

EN Digital Thermometer

FR Thermomètre Digital

DE Digital Thermometer

IT Termometro Digitale

ES Termómetro Digital


Digital Daffodil Ltd
United Kingdom
Tel: +44(0) 845 293 9358
www.digitaldaffodil.com

User Manual

1. Features

- 1) Easy to read LCD screen
- 2) Suitable for Oral, Axillary and Rectal use
- 3) Power Saving Auto-Off function
- 4) Safe for children. No mercury, no glass

2. Parts & Functions



3. How to Use

- 1) Disinfect the probe with ethyl alcohol.
- 2) Press the on/off button.
- 3) Once Lo°C shows on the display screen, position the probe for measuring.
- 4) Hold in place until the beep sounds, then remove to read measurement.

Note: To recall the last measurement taken, press and hold the power button for 5 seconds whilst the thermometer is turned off.

4. Methods of Measurement

ORAL USE

Place the probe under the patient's tongue and instruct the patient to keep their mouth closed during measurement.
Normal Oral Temperature: 37°C

AXILLARY USE

Dry the patient's armpit with a towel. Place the probe against the patient's armpit and keep the patient's arm pressed firmly against their side.
Normal Axillary Temperature: 36.7°C

RECTAL USE

Lubricate the probe with a water soluble lubricant. Gently insert the probe approximately 1cm into the patient's rectum.
Normal Rectal Temperature: 37.6°C

5. How to Replace the Battery

When the Battery is running low, the  icon will show in the display window.

- 1) Remove the battery cover.
- 2) Use a toothpick or insulated, non metal probe to remove the battery.
- 3) Insert the new battery with the positive side facing up.
- 4) Replace the battery cover.

NOTE: Do not attempt to recharge the battery. Dispose of used batteries safely.

WARNING: Placing batteries in fire can cause them to explode.

6. Precaution

- 1) Do not drop or strike the thermometer.
- 2) Only the probe on this unit is waterproof, clean only by wiping with a dry cloth.
- 3) Do not allow the thermometer to come into contact with any thinners or chemical solvents.
- 4) Do not store the unit in direct sunlight or at high temperatures.
- 5) Do not bend the thermometer.
- 6) Do not disassemble the unit other than to replace the battery.