





QuickPore™ Anti-Aging Rejuvenation & Repair NanoPen Lux Device Kit
BD-DM-K702



Step 1
Charge NanoPen Lux before use




Step 2
Begin with freshly washed skin




Step 3
Choose Daily NanoChip or Booster Nanochip




Step 4
Install NanoChip to NanoPen Lux



Step 5
Open serum, Place rubber head over the bottle



Step 6
Apply Deminuage Blanq Serum



Step 7
Turn on the NanoPen Lux by pressing and holding down the indicator light



Step 8
Use NanoPen Lux with NanoChip



Step 9
Remove and Throw Nano Chip