



**QuickPore™ Anti-Aging Rejuvenation & Repair NanoPen Prime Device Kit**  
BD-DM-K706

**Step 1**  
Install battery to NanoPen Prime

**Step 2**  
Begin with freshly washed skin

**Step 3**  
Choose Daily NanoChip or Booster Nanochip

**Step 4**  
Install NanoChip to NanoPen Prime

**Step 5**  
Open serum, Place rubber head over the bottle

**Step 6**  
Apply Deminuage Blanq Serum

**Step 7**  
Turn on the NanoPen Prime

**Step 8**  
Use NanoPen Prime with NanoChip

**Step 9**  
Remove and Throw NanoChip