



QuickPore™ Daily Moisturizing and Lifting Essential NanoPen Prime Device Kit

BD-DM-K707

Step 1
Install battery to NanoPen Prime

Step 2
Begin with freshly washed skin

Step 3
Install Daily NanoChip to NanoPen Prime

Step 4
Open serum, Place rubber head over the bottle

Step 5
Apply Deminuage Blanq Serum

Step 6
Turn on the NanoPen Prime

Step 7
Use NanoPen Prime with Daily NanoChip

Step 8
Remove and Throw Daily Nano Chip