




QuickPore™ NanoPen Prime Acne Scar Removal and Acne Spot Treatment Kit


BD-DM-K710




Step 1
Install battery to NanoPen Prime




Step 2
Begin with freshly washed skin




Step 3
Install Booster NanoChip to NanoPen Prime




Step 4
Open serum, Place rubber head over the bottle




Step 5
Apply Deminuage Blanq Serum



Step 6
Turn on the NanoPen Prime



Step 7
Use NanoPen Prime with Booster NanoChip



Step 8
Remove and Throw Booster NanoChip

Boost your inner glow

