
CONTENTS: 1 Thermometer, 1 Battery, 1 Owner's Manual, 1 Storage Case

- WARNING:**
- Read instructions thoroughly before using digital thermometer.
 - Choking Hazard: Thermometer cap and battery may be fatal if swallowed. **DO NOT** allow children to use this device without parental supervision.
 - DO NOT** use thermometer in ear. Designed use is for oral, rectal, or armpit (axilla) readings only.
 - DO NOT** place thermometer battery near extreme heat as it may explode.
 - Remove battery from the device when not in operation for a long time.
 - The use of temperature readings for self-diagnosis is dangerous. Consult your doctor for the interpretation of results. Self-diagnosis may lead to the worsening of existing disease conditions.
 - DO NOT** attempt measurements when the thermometer is wet as inaccurate readings may result.
 - DO NOT** bite the thermometer. Doing so may lead to breakage and/or injury.
 - DO NOT** attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
 - After each use, disinfect the thermometer especially in case the device is used by more than one person.
 - DO NOT** force the thermometer into the rectum. Stop insertion and abort the measurement when pain is present. Failure to do so may lead to injury.
 - DO NOT** use thermometer orally after being used rectally.
 - For children who are two years old or younger, please **DO NOT** use the devices orally.
 - If the unit has been stored at temperatures over 104°F (50°C), leave it in 41°F – 104°F (5°C – 40°C) ambient temperature for about 15 minutes before using it.
 - Babies younger than 6 months old should see a doctor when they have a fever.
 - Use of a new disposable cover is recommended; **DO NOT** reuse disposable covers.
 - DO NOT** smoke or eat/drink anything hot or cold for 10 minutes before taking an oral temperature.
 - ALWAYS** stay with your child when a thermometer is in place.

- PLEASE READ CAREFULLY BEFORE USING**
- This digital thermometer provides a quick and highly accurate reading of an individual's body temperature. The digital thermometer is intended to measure the human body's temperature in regular mode orally, rectally or under the arm, and the device is reusable for clinical or home use on people of all ages. To better understand its functions and to provide years of dependable results, please read **ALL** instructions first.

- PRECAUTION:**
- The performance of the device may be degraded should one or more of the following occur:
- Operation outside the manufacturer's slated temperature and humidity range.
 - Storage outside the manufacturer's slated temperature and humidity range.
 - Mechanical shock (for example, being dropped).
 - Patient temperature is below ambient temperature.
 - Portable and mobile RF communications can affect the device. The device needs special pre-cautions regarding EMC.
 - Keep thermometer away from Magnetic Resonance equipment.

SPECIFICATIONS:

Type: Digital Thermometer (Not Predictive)
Measure Range: 90.0°F -109.9°F (32.0°C -42.9°C)
(C°/F chosen by manufacturer)

Accuracy: ±0.2°F(±0.1°C) during 95.9°F ~ 107.6°F (35.5°C ~ 42.0°C) at 64.4°F ~ 82.4°F (18°C ~ 28°C) ambient operating range ±0.4°F(±0.2°C) for other measuring and ambient operating range


Operating mode: Direct Mode
Display: Liquid crystal display, 3 1/2 digits
Memory: For storing the last measured value
Battery: One 1.5 V DC. button battery (size LR41 or SR41, UCC 392)

Battery life: Approx. 200 hours
Dimension: 13.9cm x 2.2cm x 1.3cm (L x W x H)
Weight: Approx. 10 grams including battery
Ambient operating range: Temperature: -4°F ~ 131°F(-20°C ~ 55°C)
Relative humidity: 15% ~ 95%RH
Atmospheric Pressure : 800hPa ~ 1060hPa
Ingress Protection Rating: IP22
Classification: Type BF

Storage and transportation condition: 15% ~ 95%RH
Relative humidity: 800hPa ~ 1060hPa
Atmospheric Pressure : 1P22
Classification: Type BF




- DIRECTIONS**
- Press the On/Off Button next to LCD display. A tone will sound as the screen shows 188.8% followed by last recorded temperature. After showing the self-test temperature, the thermometer is now in the testing mode. ***Note:** If the measured temperature is less than 90°F or 32°C, the LCD will display Lo, if the measured temperature is higher than 109.9°F or 42.9°C, the LCD will display Hi.

- Position thermometer in desired location (mouth, rectum, or armpit.)
 - Oral Use: Place thermometer under tongue as indicated by "↖" position shown in Figure 2. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. Normal temperature between 95.9°F - 99.5°F (35.5°C to 37.5°C).
 - Rectal Use: Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1 cm (less than 1/2 ") into rectum. Normal temperature between 97.9°F - 100.4°F (36.6°C to 38°C).
 - Armpit Use: Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side. From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required. Normal temperature between 94.5°F - 99.1°F (34.7°C - 37.3°C).

- The degree sign flashes throughout the testing process. When flashing stops an alarm will beep for approximately 10 seconds. The measured reading will appear on the LCD simultaneously. Keep the thermometer in place until it beeps. The measurement continues even after the alarm notification. The displayed temperature will not change when thermometer is removed from its testing position. ***Note:** Normally the alarms are "Bi-Bi-Bi-Bi-"; Alarm beeps more rapidly when temperature reaches 100°F (37.8°C) or higher, and the alarms are "Bi-Bi-Bi-----Bi-Bi-Bi----- Bi-Bi-Bi-".
- To prolong battery life, press the On/Off Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around 10 minutes.
- Clean and disinfect the thermometer and then store the thermometer in its protective case.

- BATTERY REPLACEMENT**
- Replace battery when "■" appears in the lower right corner of LCD display.
 - Pull battery cover off as shown in Figure 3.
 - Gently pull out plastic circuit board with battery chamber approximately 1 cm (slightly less than 1/2") (See Figure 4).
 - Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 1.5V DC button type LR41 or SR41, UCC392, or equivalent. Be sure battery is installed with "+" polarity facing up (See Figure 5).
 - Slide battery chamber back into place and attach cover.

- CLEANING AND DISINFECTION**
- Use of probe covers is recommended to help avoid spreading germs, and to keep the thermometer clean. The thermometer should be cleaned after use even when probe covers are used.
 - Clean the thermometer using isopropyl alcohol or soap and water, being careful not to wet the window or multi-function button.
 - Store the unit in the protective case while not in use.
 - DO NOT** store or keep in direct sunlight, or in high temperatures.
 - DO NOT** drop the thermometer, as this may break the unit.
 - To prevent damage to the thermometer please note and observe the following:
 - DO NOT** use benzene, paint thinner, gasoline or other strong solvents to clean the thermometer.
 - DO NOT** use ultrasonic washing to clean the thermometer.

- FCC INFORMATION**
- Caution:** Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user authority to operate the equipment.
- *Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:
- Reorient or relocate the receiving antenna.
 - Increase the distance between the equipment and the receiver.
 - Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.

CALIBRATION

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used according to the use instruction, periodic readjustment is not required. However, we recommend checking calibration every two years or whenever clinical accuracy of the thermometer is in question. Turn on the thermometer and insert into the water bath which the temperature is steady. ASTM laboratory accuracy requirements in the display range of 98.6 to 102.2 °F (37.0 to 39.0 °C) for electronic thermometers is ±0.2°F (±0.1 °C). The above recommendations do not supersede the legal requirements. The user must always comply with legal requirements for the control of the measurement, functionality, and accuracy of the device which are required by the scope of relevant laws, directives or ordinances where the device is used.

COMPASSHEALTH
Manufactured for / Fabricado por:
Compass Health Brands Corp.
Middleburg Heights, OH 44130
www.compasshealthbrands.com
Phone / Teléfono 800.376.7263
Made in China / Hecho en China
42-70033B_06 © 2018

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