Mood and Energy
Enhancing Light
User Guide

Thank you for purchasing the TheraLite. Please read this manual in its entirety. For more information on light therapy, please visit www.carex.com.

Bright Light Benefits

Daily exposure to bright light helps to maintain your circadian rhythm, boost energy and manage sleeping patterns. Light therapy is also the first line of recommended treatment for the Winter Blues.

For many, morning sun exposure is limited by work, school or other commitments, even in the summer. The TheraLite is a bright light therapy lamp designed to fill this gap.

Included Items

TheraLite lamp, AC power supply and User Guide.

Before using your TheraLite:

1. Check that unit is in good condition.
2. Unfasten the travel cover and set-up the stand.
3. Insert the power adapter into the TheraLite and plug into a standard power outlet to ensure that it is working properly.

Safety Precautions

1. As with any electrical device, do not use the TheraLite in or near water.
2. Do not overuse the TheraLite. Overuse may cause irritability, excessive energy, or difficulty falling asleep at bedtime.
3. It is not necessary to stare into the light. Read, work or talk on the phone during your light therapy session.
4. During the first few uses, your TheraLite may give off an odor; this is harmless, will not affect the use and will quickly disappear.
5. Always use in a well ventilated area away from other heat sources and ensure that air vents are not blocked when the product is in use.
6. Only use the supplied power adapter.

Disclaimer

TheraLite is an innovative light supply system and is not a listed medical device in the USA.

Warranty

Your lamp is maintenance-free and should provide years of service. The TheraLite has a one-year limited warranty from the date of purchase by the original purchaser. Please retain proof of purchase to ensure eligibility for warranty coverage.

Carex Health Brands, may choose to repair or replace a defective unit. Shipping costs are the responsibility of the product owner.

Misuse, abuse, accidental or incorrect operation, contrary to the instructions, will void a warranty claim. Only use the supplied power adapter with the TheraLite.

-
1. Open the cover by gently detaching the 2 retention clips.
 2. Swing the hinged cover so it flips under the lamp.
-

3. Select one of the 3 angle slots and insert the stand tab. Lowest is 20°, middle is 25° and the top angle of the lamp is 30° from vertical.
4. Loosen the back knob to adjust the height and orientation of the lamp and re-tighten at desired position. Be careful not to overtighten and damage the lamp.
5. When the lamp is no longer needed, it may closed and stored. Always ensure that the both the lamp and plug are cool before storage.

Using Your TheraLite for Therapy


For those with mild to moderate Winter Blues, sleep or circadian rhythm issues, regular daily use of the TheraLite can help provide relief. For most conditions, use shortly after waking or in the morning.

Place the TheraLite in your line of sight, where possible, to maximize the light reaching your eyes. You should be facing the lamp with eyes open during sessions, but you do not need to look directly at the light.

Over the first few days of use, assess the effect and, if necessary, lengthen or shorten the session time based on your response.

Symptoms of Overuse

Negative effects are rare with light therapy. However, if you experience increased irritability, excessive energy, and/or any consistent discomfort during or after use, decrease session time or move the light further away. Using the lamp late in the afternoon or in the evening is not advised, as this may disrupt your sleep. It is recommended that you consult your doctor before beginning any bright light regimen.

Specifications


UL Listed Power Supply

Power:	24V, 1A output, 100-240V power adapter
Size:	9.2" x 5.3" x 2.1" / 23.3 cm x 13.5 cm x 5.2 cm
Lens:	High-impact polycarbonate
Light Intensity:	10,000 LUX at 8.5"/22 cm
Light Source:	104 LED's, 3000 Kelvin, CRI - 80, UV-free
Angles:	20°, 25°, 30° from vertical

Light Therapy Basics

LUX or the light intensity at a specific distance from the light is the common term used in bright light therapy. The chart below details the average treatment times for the Winter Blues based on your eyes proximity to the lamp.

DISTANCE	SESSION TIME
10" / 25 cm	30 minutes
14" / 36 cm	60 minutes
20" / 50 cm	120 minutes

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