



## FAQs

**Q: Will the SpaceCHQR be painful?**

A: Pain should not be felt while using this product. The SpaceCHQR is a series of medical grade stainless steel gauges that have been highly polished on the edges to protect your gums . There maybe some soreness due to tooth movement but contact your doctor if you feel any sharp pain during insertion of the SpaceCHQR.

**Q: How do I know if the contact points between my teeth is tight?**

A: The SpaceCHQR should be able to be inserted and removed passively. **DO NOT FORCE THE GAUGE BETWEEN THE TEETH.** If you find the gauge does not fit between any of the teeth, move to a finer gauge.

**Q: What do I do if the Green, Amber and Red does not go in between the teeth?**

A: Ensure to brush and floss your teeth first. If you still cannot get the Red gauge in call your doctor.

**Q: What do I do if there is food or foreign object that prevents the gauge from going in the proximal contact points?**

A: Use floss to clear the contact points to remove. If you cannot correct this with floss, call your doctor.

**Q: What do I do if the gauge is bent and cannot be inserted?**

A: The stainless steel gauge can be distorted during use. If it cannot be inserted due to damage, call your doctor and replace the device.

**Q: What do I do if the gauge is broken or torn?**

A: Do not use the SpaceCHQR if the gauge is torn or broken as this may cause injury. Call your doctor and request for a replacement device.